

amsbio

magtivio

# MagSi-DNA Animal

FAST & COST-EFFECTIVE  
EXTRACTION OF  
GENOMIC DNA FROM  
ANIMAL SAMPLES



## DNA EXTRACTION FROM ANIMAL SAMPLES

---

MagSi-DNA Animal allows fast and cost-effective extraction of DNA from various samples like blood, semen, hairs, saliva/swabs or lysed tissue. This universal DNA purification kit is optimized to extract DNA from these sample materials with the highest purity and delivering DNA which is suitable for genotyping assays or other qPCR based analysis. The extraction chemistry is validated on different species, e.g. horse, swine, dog, cattle and can be customized to meet any specific requirements of yields, purity, working volumes. The kit can be used on KingFisher™ Flex or similar systems or on many other automated DNA purification instruments.

## GENERAL FEATURES

- Short protocols, complete processing at room temperature possible (after sample lysis)
- High DNA yield and purity
- Suitable for many genomic applications such as SNP genotyping, DNA sequencing, NGS
- Suitable for animal parental testing / breeding identity checks
- Preparation time for 96 samples: 20 min after lysis


## QUALITY

- Validated procedures for many sample materials (e.g., blood, semen, saliva, swabs (saliva/oral fluid), hair, lysed tissue)
- High molecular weight and pure DNA suitable for long-time storage
- Efficient removal of PCR inhibitors

**Table 1.** Results of typical DNA concentrations obtained from various sample materials. DNA purification was performed using a MagMax EXPRESS 96 magnetic bead processor

Sample	DNA concentration <sup>1</sup> , ng/μL (Qubit)				
	Dog	Cat	Horse	Cattle	Swine
Blood, 20 μL	5.5	3.5	2.6	9.3	5
Semen, 30 μL	74	nt	170	26	11
Swab/Saliva2, 100 μL	2.2	5.6	nt	nt	nt
Hair, 8 pcs	nt	nt	11.8	5.7	13.7
Tissue, 4 mm	nt	nt	nt	17.8	3.1

nt: not tested; 1 100 μL elution volume; 2 from DNA Genotek PG 100 tubes



**Figure 1.** Exemplary STR profile, DNA purified from 20 μL of dog blood

## FLEXIBILITY

- Suitable for processing various sample materials
- Suitable for small, medium and high-throughput automation
- Small elution volumes


## EASY TO AUTOMATE

- Minimal accessory requirements
- KingFisher™ / Biosprint 96 protocols available on request
- Consumables for KingFisher™ Flex 96 / MagMax EXPRESS 96 / Biosprint 96 available
- Compatible with liquid handling robots (e.g. Hamilton®, TECAN®, oKtopure™)
- Magnetic separators for microtubes and microplates available

Description	Pack Size	Catalogue No.
MagSi-DNA Animal	96 preps	MDKT00140096
MagSi-DNA Animal	10 x 96 preps	MDKT00140960

AMS BIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

V1

 **UK & Rest of the World**  
184 Park Drive, Milton Park  
Abingdon OX14 4SE, U.K.  
T: +44 (0) 1235 828 200  
F: +44 (0) 1235 820 482

 **North America**  
1035 Cambridge Street,  
Cambridge, MA 02141.  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218

 **Germany**  
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880

 **Switzerland**  
Via Lisano 3  
(CP.683)  
6900 Massagno  
T: +41 (0) 91 604 55 22  
F: +41 (0) 91 605 17 85