


CellO-IF: The Smarter Way to Label

All-in-one Reagent for Immunofluorescence Labelling of Organoids, Spheroids and Cells



Liver Cancer Spheroid labelled using CellO-IF technology - Tok et al (2022)

Simplify the steps of immunofluorescence labelling and maximize your results. The CellO-IF innovative technology allows you to stain organoids, spheroids or cells in just 8 hours, while still in hydrogel or extracellular matrix, preserving delicate structures and cellular integrity.

This all-in-one reagent removes the need for harvesting, clearing, transferring and centrifuging, preventing sample loss and ensuring consistent accurate results.

Are you having issues with unsuccessful labelling of precious samples?

Features

- High quality images in hours, not days.
- Reduced number of reagents and steps.
- Label organoids, spheroids, and cells while still in hydrogel or extracellular matrix.
- Reduces background noise for high-resolution images.
- Boosted antigenicity for clearer visualization of target.

Benefits

- Faster Results.
- Streamlined Workflows.
- Effortless labelling.
- Enhanced Clarity.
- Optimize detection.

Request a free sample: info@amsbio.com


PROTOCOL FOR USE

1. Fix in 4% paraformaldehyde (PFA) for 1 hr at 37°C.
2. Rinse with PBS X3, 10 min at 37°C.
3. Incubate for 1.5-2 hrs at 37°C with primary antibodies diluted in Cello-IF.
4. Wash with Cello-IF X3, 10 min at 37°C.
5. Incubate for 1.5-2 hrs with secondary antibodies diluted in Cello-IF.
6. Wash with PBS X3, 10 min at 37°C.
7. Mount & image under microscope.



	Traditional Workflow	Workflow with CELLO-IF
Reagents	>3	1
Steps	24	8
Transferring & Pipetting	>3	0
Time	~44 hours	8 hours

Available in two formats for labelling both 3D and 2D samples producing high quality images.



Left: Induced Pluripotent Stem Cell Derived Neurons (MAP2 - green, CREB - red, DAPI - blue).

Right: Brain Organoid (PAX2 - red, DAPI - blue).

"We really didn't think imaging the neuronal organoids would be as easy because they are so fragile and difficult to section with traditional methods. Using Cello™-IF simplified our workflow and it has great potential to make high-throughput imaging of organoids feasible. We got beautiful images using Cello-IF on both 2D and 3D cultures."

Didem Demirbas, PhD - Jenny Lai, MD-PhD Candidate, Boston Children's Hospital, Harvard Medical School



AMSBIO LLC
USA & Canada
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
+1 (800) 987-0985
F: +1 (617) 945-8218



AMSBIO Europe BV
EU
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142



AMSBIO Biotechnology (Europe) Ltd
UK & Rest of the World
184 Park Drive, Milton Park,
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482



AMSBIO Europe BV
Deutschland
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

AMSBIO™ / AMS Biotechnology (Europe) Ltd. AMSBIO™ is the global source of all the products listed here.
All trademarks mentioned here are the property of their respective owners.