

O-Glycosidase

O-Glycosidase (O-glycopeptide endo-D-galactosyl-N-acetyl- α -galactosaminohydrolase, EC 3.2.1.97) cleaves only unsubstituted Gal β (1-3)GalNAc α disaccharides attached to the serine or threonine residues of glycoproteins or glycopeptides. Substitutions such as sialic acid, galactose, fucose or N-acetylglucosamine must first be removed with the appropriate exoglycosidase prior to treatment with O-Glycosidase. At minimum, a neuraminidase such as α (2-3,6,8,9) neuraminidase is almost always required to remove sialic acid.

There is no activity on α -GalNAc linked either to protein or carbohydrate. Although limited by its strict specificity, O-Glycosidase, in conjunction with other exoglycosidases, is still the method of choice for removing O-linked sugars from glycoproteins. The protein remains intact, as does the disaccharide.

O-Glycosidase is isolated from a clone of *Streptococcus pneumonia* (formerly *Diplococcus pneumonia*).

O-Glycosidase is useful for:

- Determining O-glycosylation in proteins
- Studying the effects of O-glycosylation in binding and antigenic studies
- Removing O-linked sugars for X-Ray crystallography and protein sequencing

Specifications

Activity

≥16 U/mg, ≥17 U/mL

Storage

Store at 4°C. Do not freeze

Formulation

Enzyme is provided as a sterile solution in 50 mM Sodium Phosphate pH 7.5.

Stability

Stable at least 12 months when stored properly. Several days exposure to ambient temperature will not reduce activity.

Product Description

Molecular Weight

~ 180,000 Daltons

Purity

O-Glycosidase is tested for contaminating protease as follows; 10 μ g of denatured BSA is incubated for 24 hours at 37°C with 2 μ L of enzyme. SDS-PAGE analysis of the treated BSA shows no evidence of degradation.

The production host strain has been extensively tested and does not produce any detectable glycosidases.

Specificity:

Serine- or threonine-linked unsubstituted Gal β (1-3)GalNAc α .

pH Range:

Optimum: pH 5 Range: pH 5 - 7

The supplied buffer concentrate provides the optimal pH for enzyme activity with the standard substrate. If glycosidase treatment is performed at suboptimal pH because of glycoprotein solubility or activity requirements, expect some diminution in enzyme activity.

Assay

One unit of O-Glycosidase is defined as the amount of enzyme required to produce 1 μ mole of *p*-nitrophenol (pNP) in 1 minute at 37°C pH 5 from *p*-nitrophenyl-2-acetamido-2-deoxy-3-O-(β -D-galactopyranosyl)- α -D-galactopyranoside.

Reagents

- 5X PNGase buffer 5- 250 mM sodium phosphate pH 5

Note: O-Glycosidase cleaves methylumbelliferyl- α -D-N-acetylgalactosamide (but not pNP- α -D-N-acetylgalactosamide) at 0.1% the rate of *p*-nitrophenyl-2-acetamido-2-deoxy-3-O-(β -D-galactopyranosyl)- α -D-galactopyranoside.

Suggestions for Use

Procedure for Deglycosylation

1. Add up to 100 µg of glycoprotein to tube.
2. Add water to 13 µL and 4 µL 5X Reaction Buffer 5
3. Add 1 µL α(2-3,6,8,9) neuraminidase
4. Add 2 µL O-Glycosidase
5. Incubate at 37°C for 1 hour

Cleavage may be monitored by SDS-PAGE if the size differential between native and de-O-glycosylated protein is sufficient for detection.

References

1. Bhavanandan, V.P., J. Umemoto and E.A. Davidson. Characterization of an endo-α-N-acetyl-galactosaminidase from *Diplococcus pneumonia*. *Biochem Biophys* 70:738-745 (1976).
2. Fan, J. Q., K. Yamamoto, H. Kumagai and T. Tochikura. Induction and efficient purification of endo-α-N-acetyl-D-galactosaminidase from *Alcaligenes* sp. *Agric Biol Chem* 54:233-234 (1990).
3. Glasgow, L.R., J. C. Paulson and R. L. Hill. Systematic purification of five glycosidases from *Streptococcus pneumonia*. *J. Biol Chem* 252:8615-8623 (1977). Iwase, H. and K. Hotta. Release of O-linked glycoprotein glycans by endo-α-N-acetyl-D-galactosaminidase. *Methods Mol Biol* 14:151-159 (1993).
4. Umemoto, J., V. P. Bhavanandan and E. A. Davidson. Purification and properties of an endo-α-N-acetyl-D-galactosaminidase from *Diplococcus pneumonia*. *J. Biol Chem.* 252:8609-8614 (1977).

Order Information

Catalog No.	Product Description	Package Size	Temp. °C
100686-1	O-Glycosidase (<i>Streptococcus pneumonia</i> recombinant)	60 µL	+4

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85