

# PRODUCT INFORMATION

**Coomassie® Brilliant Blue G 250 Powder**

**Cat. No. 17524**

## Product description:

<b>General</b>	The "G" in the name refers to a blue dye with a greenish hue
<b>Application</b>	<ul style="list-style-type: none"> <li>• Protein assay with Bradford reagent<sup>1</sup></li> <li>• Staining of proteins in polyacrylamide and agarose gels</li> <li>• Determination of CMC of detergents<sup>2</sup></li> </ul>
<b>Storage</b>	Store powder dry and dark at +15 °C – +30 °C. If stored as indicated, shelf life of powder is ca. 3 – 5 years.
<b>Solubility</b>	Soluble in water (1 mg/ml, dark blue solution). For complete solution heating may be required.

## Staining protocol:

<b>Solutions</b>	Fixing solution	10 % glacial acetic acid, 40 % EtOH (10 ml glacial acetic acid/40 ml EtOH/50 ml H <sub>2</sub> O)
	Stock solution I	0.2 % Brilliant Blue G in 90 % EtOH (100 mg in 50 ml 90 % EtOH) Solution is stable at RT for several months.
	Stock solution II	20 % Acetic acid
	Destainer I	conforms to fixing solution
	Destainer II	10 % glacial acetic acid, 20 % EtOH
<b>Staining</b>	<ol style="list-style-type: none"> <li>1. Fix gel for 30 min. in 100 ml fixing solution.</li> <li>2. Prior to use add 50 ml stock solution I to 50 ml stock solution II.</li> <li>3. Stain gel in solution for 20 min.</li> </ol>	

<b>Destaining</b>	<ol style="list-style-type: none"> <li>1. For a first rapid destain, incubate gel for ca. 0.5 - 5 min. (dependent on gel thickness) in destainer I.</li> <li>2. Destain gel further as long as needed in destainer II (visual inspection).</li> <li>3. Rinse destained gel 2x for 5 min. in H<sub>2</sub>O dist.</li> </ol>
-------------------	---

<sup>1</sup>Bradford, M. M., Anal. Biochem., 72, 248 – 254 (1976), <sup>2</sup>Samsonoff, C., et al., J. Colloid and Interface Science, 109(2), 325 – 329 (1986).

\*Coomassie is a tradename of ICI, UK **Ver. 04/07**



**UK & Rest of World**  
184 Milton Park, Abingdon  
OX14 4SE, Oxon, UK  
Tel: +44 (0) 1235 828 200  
Fax: +44 (0) 1235 820 482

**Switzerland**  
Centro Nord-Sud 2E  
CH-6934 Bioggio-Lugano  
Tel: +41 (0) 91 604 55 22  
Fax: +41 (0) 91 605 17 85

**Deutschland**  
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
Tel: +49 (0) 69 779099  
Fax: +49 (0) 69 13376880

**United States**  
23591 El Toro Rd, Suite #167  
Lake Forest, CA 92630  
Tel: +1 800 987 0985  
Fax: +1 949 265 7703

**amsbio**  
info@amsbio.com