

## Arthrogen-CAIA® Arthritogenic Monoclonal Antibody

Catalog # 53040

*For Research Use Only - Not Human or Therapeutic Use*

DESCRIPTION:	Mouse anti-type II collagen 5-clone monoclonal antibody cocktail kit with LPS from <i>E. Coli</i> O111:B4
APPLICATION:	Use for collagen antibody induced arthritis (CAIA) in mice
QUANTITY:	Cocktail: 10 mg/ml x 4 ml LPS: 0.5 mg/ml x 3.5 ml
FORM:	A cocktail of 5 monoclonal antibodies dissolved in 0.05M phosphate buffered saline, pH 7.4 which recognize the conserved epitopes on various species of type II collagen: clone A2-10 (IgG2a), F10-21 (IgG2a), D8-6 (IgG2a), D1-2G (IgG2b), and D2-112 (IgG2b). Clones A2-10, D1-2G, and D2-112 recognize individual epitopes clustered within the 167 amino acid peptide fragment called LyC1 (124-290) of the CB11 fragment (124-402). Clones F10-21 and D8-6 recognize epitopes within the 83 amino acid peptide fragment called LyC2 (291-374) of the CB11 fragment.
SOURCE:	Mouse
CROSS-REACTIVITY:	Cross-reacts to most species of type II collagen including mouse, porcine, chick*, bovine, human, rat, monkey, rabbit, equine, and dog.  *D2-112 does not cross-react with chick type II collagen.
ANIMALS:	7-8 weeks or older, high responder mice: DBA/1, Balb/c, B10.RIII, C.B-17, scid/scid, or 129/Sv. Low responder mice: C57BL/6 or C57BL/6 background. Sensitivity to CAIA may vary depending on mouse vendor and facility. Specific pathogen free (SPF) housing conditions are strongly recommended over conventional housing conditions as bacterial contamination may reduce the immune response of the animals resulting in an attenuated arthritis. Therefore, running a pilot study using a small number of animals before conducting a large-scale study is strongly recommended. Please contact AMSBIO customer support ( <a href="mailto:info@amsbio.com">info@amsbio.com</a> ) for guidance regarding animals and housing conditions.
USAGE:	Administration of the monoclonal antibody cocktail by IV injection (i.e. tail vein) is recommended; however, intraperitoneal (IP) injection may also be used. Moreover, if LPS use is not desired, please note that the severity of arthritis tends to be lower.
PROTOCOLS:	

- A. Inducing arthritis with a combination of the monoclonal antibody cocktail and LPS in CAIA susceptible mice (DBA/1, Balb/c, B10.RIII, C.B-17, scid/scid, or 129/Sv)

Day 0: Inject 1.5 mg of 5-clone cocktail by IV or IP injection.

Day 3: Inject 25-50 µg of LPS by IP injection.

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

UK & Rest of the World   
184 Park Drive, Milton Park  
Abingdon OX14 4SE.  
T: +44 (0) 1235 828 200  
F: +44 (0) 1235 820 482

North America   
1035 Cambridge Street,  
Cambridge, MA 02141.  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218

Europe   
Berenkoog 41,  
1822 BH Alkmaar,  
Netherlands  
T: +31 (0) 72 8080244  
F: +31 (0) 72 8080142

Switzerland   
Via Lisano 3,  
(CP.683)  
CH-6900  
T: +41 (0) 91 604 55 22  
F: +41 (0) 91 605 17 85

Germany   
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880

NOTE: Moderate to severe arthritis is observed on day 3-4 and peaks around day 7-10.

- B. Inducing arthritis with a combination of the monoclonal antibody cocktail and LPS in CAIA low responder mice (C57BL/6 or C57BL/6 background)

Day 0: Inject 5 mg of 5-clone cocktail by IV or IP injection.

Day 3: Inject 25-50 µg of LPS by IP injection.

NOTE: Moderate to severe arthritis is observed on day 3-4 and peaks around day 7-10. An injection of LPS on day 10-14 can be used to re-stimulate inflammatory arthritis.

- C. Inducing arthritis with the monoclonal antibody cocktail without LPS in DBA/1 and Balb/c mice

Day 0: Inject 6-10 mg of 5-clone cocktail by IV injection.

NOTE: Arthritis can be observed on day 2-3. This protocol has not been confirmed in other CAIA susceptible mouse strains such as B10.RIII, C.B-17, scid/scid, or 129/Sv mice.

STORAGE:	-80°C
STABILITY:	2 years
NOTES:	N/A
REFERENCES:	K. Terato <i>et al. J. Immunol.</i> <b>148</b> : 2103-2108 (1992) K. Terato <i>et al. Autoimmunity</i> <b>22</b> : 137-147 (1995) S. Yoshino <i>et al. J. Immunol Methods</i> <b>343</b> : 49-55 (2009)

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

UK & Rest of the World   
184 Park Drive, Milton Park  
Abingdon OX14 4SE.  
T: +44 (0) 1235 828 200  
F: +44 (0) 1235 820 482

North America   
1035 Cambridge Street,  
Cambridge, MA 02141.  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218

Europe   
Berenkoog 41,  
1822 BH Alkmaar,  
Netherlands  
T: +31 (0)728080244  
F: +31 (0)728080142

Switzerland   
Via Lisano 3,  
(CP.683)  
CH-6900  
T: +41 (0) 91 604 55 22  
F: +41 (0) 91 605 17 85

Germany   
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880