

Notch Signaling Pathway
Notch CSL Reporter – HEK293 Cell line
Catalog #: 60652

Background

The Notch signaling pathway controls cell fate decisions in vertebrate and invertebrate tissues. Notch signaling is triggered through the binding of a transmembrane ligand to Notch transmembrane receptor (NOTCH1/ NOTCH2/NOTCH3/NOTCH4) on a neighboring cell. This results in proteolytic cleavage of the Notch receptor, releasing the constitutively active intracellular domain of NOTCH (NICD). NICD translocates to the nucleus and associates with transcription factors CSL (CBF1/RBPJk/Suppressor of Hairless/Lag-1) and coactivator Mastermind to turn on transcription of Notch-responsive genes.

Description

The Notch CSL Reporter – HEK293 cell line contains the firefly luciferase gene under the control of Notch-response elements (CSL responsive elements) alone with expression construct for Notch1ΔE (NOTCH1 that has a deletion of the entire extracellular domain) stably integrated into HEK293 cells. Inside the cells, the Notch1ΔE can be cleaved by γ-secretase. The active Notch1 NICD is released to nucleus and induces the constitutive expression of luciferase reporter. The cell line is validated for the inhibition of the expression of luciferase reporter using a known inhibitor of the Notch signaling pathway.

Application

- Monitor Notch signaling pathway activity.
- Screen for inhibitors of the Notch signaling pathway.

Format

Each vial contains ~1.5 X 10⁶ cells in 1 ml of 10% DMSO.

Storage

Immediately upon receipt, store cells in liquid nitrogen.

General culture conditions

Cells should be grown at 37°C with 5% CO₂ using MEM medium (Hyclone #SH30024.01) supplemented with 10% FBS (Invitrogen #26140-079), 1% non-essential amino acids (Hyclone #SH30238.01), 1 mM Na-pyruvate (Hyclone #SH30239.01), 1% Penicillin/Streptomycin (Hyclone SV30010.01), plus 400 µg/ml of Geneticin (invitrogen #11811031) and 100 µg/ml of

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
 184 Park Drive, Milton Park
 Abingdon OX14 4SE, UK
 T: +44 (0)1235 828 200
 F: +44 (0) 1235 820 482



North America
 1035 Cambridge Street,
 Cambridge, MA 02141
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218



Germany
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880



Switzerland
 Centro Nord-Sud 2E
 CH-6934 Bioggio-Lugano
 T: +41(0) 91 604 55 22
 F: +41(0) 91 605 17 85

Hygromycin B (Hyclone #SV30070.01). It may be necessary to adjust the percentage of CO₂ in the incubator depending on the NaHCO₃ level in the basal medium.

It is recommended to quickly thaw the frozen cells from liquid nitrogen in a 37°C water-bath, transfer to a tube containing 10 ml of growth medium without Geneticin and Hygromycin, spin down cells, and resuspend cells in pre-warmed growth medium without Geneticin and Hygromycin. Transfer resuspended cells to a T25 flask and culture in a CO₂ incubator at 37°C overnight. The next day, replace the growth medium with fresh growth medium without Geneticin and Hygromycin, and continue growing culture in a CO₂ incubator at 37°C until the cells are ready to be split. At first passage, switch to growth medium containing Geneticin and Hygromycin. Cells should be split before they reach complete confluence.

To passage the cells, rinse cells with phosphate buffered saline (PBS), and detach cells from the culture vessel with 0.05% Trypsin/EDTA. Add complete growth medium and transfer to a tube, spin down the cells, then resuspend cells and seed appropriate aliquots of cell suspension into new culture vessels. Subcultivation ration: 1:10 to 1:20 weekly or twice a week.

To freeze down the cells, rinse cells with phosphate buffered saline (PBS), and detach cells from culture vessel with Trypsin/EDTA. Add complete growth medium and transfer to a tube, spin down cells, and resuspend in freezing medium (10% DMSO + 90% FBS). Place at -80 °C overnight and place in liquid nitrogen the next day. Alternatively, vials may be placed directly in liquid nitrogen.

Functional Validation and Assay Performance

The following assays are designed for 96-well format. To perform the assay in different tissue culture formats, the cell number and reagent volume should be scaled appropriately.

Materials Required but Not Supplied

- DAPT (Selleckchem # S2215): inhibitor of Notch pathway (γ-secretase inhibitor). Prepare 10 mM of stock solution in DMSO.
- Assay medium: growth medium without Geneticin and Hygromycin
- 96-well tissue culture plate or 96-well tissue culture-treated white clear-bottom assay plate
- ONE-Step™ Luciferase Assay System (Cat. #60690)
- Luminometer

Mycoplasma testing

The cell line has been screened using the PCR-based Venor® GeM Mycoplasma Detection kit (Sigma-Aldrich) to confirm the absence of Mycoplasma species.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.


Inhibition of Notch reporter activity by inhibitor of Notch signaling pathway in Notch CSL Reporter – HEK293 cells

1. Harvest Notch CSL Reporter – HEK293 cells from culture in growth medium and seed cells at a density of 35,000 cells per well into white clear-bottom 96-well microplate in 45 μ l of assay medium.
2. Dilute the inhibitor (DAPT) stock in assay medium. Add 5 μ l of diluted inhibitor to the test wells. The final concentration of DMSO in assay medium can be up to 0.5%.
Add 5 μ l of assay medium with same concentration of DMSO without inhibitor to control wells.
Add 50 μ l of assay medium with DMSO to cell-free control wells (for determining background luminescence).
Set up each treatment in at least triplicate.
3. Incubate the plate at 37°C in a CO₂ incubator for 24 hours.
4. The next day, perform luciferase assay using the ONE-Step™ Luciferase Assay System according to the protocol provided: Add 100 μ l of ONE-Step™ Luciferase reagent per well and rock at room temperature for ~15 minutes. Measure luminescence using a luminometer.
If using other luciferase reagents from other vendors, follow the manufacturer's assay protocol.
5. Data Analysis: Obtain the background-subtracted luminescence by subtracting the average background luminescence (cell-free control wells) from the luminescence reading of all wells.


OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

Figure 1. Response of Notch CSL Reporter – HEK293 cells to Notch pathway inhibitor.

A. DAPT Blocked Notch Reporter Activity. The CSL Reporter-HEK293 cells contain the firefly luciferase gene under the control of Notch-response elements (CSL-responsive elements) stably integrated into HEK293 cells, but not the expression construct for Notch1 Δ E, so it is unable to express the Notch luciferase reporter.



B. DAPT Inhibition Dose Response Curve. The results were shown as percentage of luminescence. The background-subtracted luminescence of cells in the absence of DAPT was set at 100%.



References

Lu FM *et al.* (1996) Constitutively active human Notch1 binds to the transcription factor CBF1 and stimulates transcription through a promoter containing a CBF1-responsive element. *Proc. Natl. Acad. Sci. USA* **93(11)**: 5663-5667.

Kanungo *et al* (2008) The Notch signaling inhibitor DAPT down-regulates cdk5 activity and modulates the distribution of neuronal cytoskeletal proteins. *J. Neurochem.* **106**: 2236

License Disclosure

Purchase of this cell line grants you with a 10-year license to use this cell line in your immediate laboratory, for research use only. This license does not permit you to share, distribute, sell, sublicense, or otherwise make the cell line available for use to other laboratories, departments, research institutions, hospitals, universities, or biotech companies. The license does not permit the use of this cell line in humans or for therapeutic or drug use. Inappropriate use or distribution of this cell line will result in revocation of the license and result in an immediate cease of sales and distribution of our products to your laboratory. Publications using this cell line should reference BPS Bioscience, Inc., San Diego

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

Related Products

<u>Product name</u>	<u>Catalog#</u>	<u>Size</u>
Notch Pathway Reporter Kit	60509	500 rxns
TCF/LEF Reporter Kit (Wnt / β -catenin signaling Pathway)	60500	500 rxns
SRE Reporter Kit (MAPK/ERK signaling Pathway)	60511	500 rxns
CRE/CREB Kit (cAMP/PKA signaling Pathway)	60611	500 rxns
Wnt Signaling Pathway TCF/LEF (Luc) Reproter HEK293	60501	2 vials
Hedgehog Signaling PathwayGli Reporter-NIH3T3 Cell Line	60409	2 vials
ERK Signaling Pathway SRE Reporter-HEK293 Cell Line	60406	2 vials
JNK Signaling Pathway AP1 Reporter-HEK293 Cell Line	60405	2 vials
JAK/STAT Signaling Pathway ISRE Reporter-HEK293	60409	2 vials
ONE-Step™ Luciferase Assay System	60690-1	10 ml
ONE-Step™ Luciferase Assay System	60690-2	100 ml

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

150408

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85