

Product Specification Sheet

Product Name	Stemgent® Alkaline Phosphatase Staining Kit II
Description	The undifferentiated state of embryonic stem (ES) and induced pluripotent stem (iPS) cells can be characterized by a high level of alkaline phosphatase (AP) expression which, along with the expression of other surface markers, indicates undifferentiated cells with the potential to self-renew. AP is a hydrolase enzyme responsible for dephosphorylating molecules such as nucleotides, proteins, and alkaloids under alkaline conditions. When fixed ES or iPS cells are stained using the AP Staining Kit, undifferentiated cells appear red or purple, whereas differentiated cells appear colorless.
Catalog Number	AMS.00-0055
Size	50 assays
Kit Components	Fix Solution: 25 ml (Cat. No. 09-0042) AP Staining Solution A: 10 ml (Cat. No. 09-0043) AP Staining Solution B: 10 ml (Cat. No. 09-0065) AP Staining Solution C: 10 ml (Cat. No. 09-0066)
Storage and Stability	Store at 4°C. Do not freeze. Stable for 6 months from date of receipt when stored as directed.
Quality Control	The AP Staining Kit II is functionally tested on human and mouse ES cells to ensure product quality (Figure 1).




Figure 1. Phenotypic assessment of R1 mouse ES cells and H1 human ES cells using AP staining. Cells were fixed and stained following the AP Staining Kit protocol. [A] Undifferentiated R1 cells, [B] undifferentiated H1 cells, and [C] partially differentiated H1 cells.

- | | |
|-------------------|--|
| References | <ol style="list-style-type: none"> 1. Pease, S., Braghetta, P., Gearing, D., Grail, D., and Williams, R.L. (1990) Isolation of embryonic stem (ES) cells in media supplemented with recombinant leukemia inhibitory factor (LIF). <i>Dev. Biol.</i> 141: 344-352. 2. Smith, A.G., Nichols, J., Robertson, M., and Rathjen, P.D. (1992) Differentiation inhibiting activity (DIA/LIF) and mouse development. <i>Dev. Biol.</i> 151: 339-351. 3. Takahashi, K., Tanabe, K., Ohnuki, M., Narita, M., Ichisaka, T., Tomoda, K., and Yamanaka, S. (2007) Induction of pluripotent stem cells from adult human fibroblasts by defined factors. <i>Cell</i> 131: 861-872. 4. Yu, J., Vodyanik, M.A., Smuga-Otto, K., Antosiewicz-Bourget, J., Frane, J.L., Tian, S., Nie, J., Jonsdottir, G.A., Ruotti, V., Stewart, R., Slukvin, I.I., and Thomson, J.A. (2007) Induced pluripotent stem cell lines derived from human somatic cells. <i>Science</i> 318: 1917-1920. |
|-------------------|--|

For research use only. Not for use in diagnostic procedures.

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85