

Data Sheet

Fluorogenic HDAC Class2a Assay Kit

Catalog #: AMS.50041

DESCRIPTION: The *Fluorogenic HDAC Class2a Assay Kit* is a complete assay system designed to measure histone deacetylase (HDAC) class2a activity for screening and profiling applications. It comes in a convenient 96-well format, with all the reagents necessary for 100 fluorescent HDAC activity measurements. In addition, the kit includes purified HDAC4 and a potent HDAC inhibitor, Trichostatin A, for use as a positive and negative control. The *Fluorogenic HDAC Class 2a Assay Kit* is based on a unique fluorogenic substrate and developer combination. This assay method eliminates dealing with the radioactivity, extraction, and chromatography aspects of traditional assays. Using this kit, only two simple steps on a microtiter plate are needed to analyze the HDAC class2a activity level. First, the HDAC fluorometric substrate is incubated with purified HDAC4 enzyme. The deacetylation sensitizes the substrate so subsequent treatment with the Lysine Developer produces a fluorophore that can then be measured using a fluorescence reader.

COMPONENTS:

Catalog #	Component	Amount	Storage	
50004	HDAC4 human recombinant enzyme	1 µg	-80°C	Avoid freeze/ thaw cycles!
50040	Fluorogenic HDAC substrate class 2A (5 mM)	50 µl	-80°C	
50030	2x HDAC Developer (contains Trichostatin A) (50 µM)	6 ml	-80°C	
	Trichostatin A (1 mM) in DMSO	100 µl	-20°C	
50031	HDAC assay buffer	10 ml	-20°C	
	black, low binding NUNC black microtiter plate	1 plate	Room temp.	

MATERIALS OR INSTRUMENTS REQUIRED BUT NOT SUPPLIED:

0.1% solution (1 mg/ml) of bovine serum albumin (BSA) in water
 Fluorimeter capable of excitation at 350-380 nm and detection at 440-460 nm
 Adjustable micropipettor and sterile tips
 Rotating or rocker platform

APPLICATIONS: Great for studying enzyme kinetics and screening small molecular inhibitors of class2a HDACs for drug discovery and HTS applications. This kit is suitable for Class 2a and class 4 HDACs (HDACs 4, 5, 7, 9, 11). For class 1 and class 2b HDACs (HDACs HDACs 1, 2, 3, 6, 10), we recommend the Fluorogenic HDAC Assay Kits, Cat. #50033 or #50034, or for HDAC8, we recommend the Fluorogenic HDAC8 Assay Kit, Cat. # 50068.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

STABILITY: One year from date of receipt when stored as directed.

REFERENCE(S): Ontoria, J.M., *et al.*, *J. Med. Chem.* 2009 Nov 12;**52(21)**:6782-9.

ASSAY PROTOCOL:

Immediately prior to assay:

- 1) Dilute **Trichostatin A** 1 mM stock 10-fold with **HDAC Assay Buffer** to make a 100 μ M solution. (Make only sufficient quantity needed for the assay; store remaining 1 mM **Trichostatin A** stock solution in aliquots at -80°C.)
- 2) Dilute **HDAC substrate** 5 mM stock 25-fold with **HDAC Assay Buffer** to make a 200 μ M solution. (Make only sufficient quantity needed for the assay; store remaining 5 mM stock solution in aliquots at -80°C.)
- 3) Dilute **HDAC4** in **HDAC Assay Buffer** to 12 pg/ μ l (60 pg/reaction)*. Aliquot any remaining enzyme and store undiluted at -80°C. Keep diluted enzyme on ice. Discard any remaining diluted enzyme after use. **Note: optimal enzyme concentration may vary with the specific activity of the enzyme.*

Step 1:

In duplicate, add the reaction mixtures (below) to the microtiter black plate as follows:

- 1) Prepare the master mixture: N wells \times (5 μ l **HDAC substrate** (200 μ M) + 5 μ l BSA (1 mg/ml) + 30 μ l **HDAC Assay Buffer**). Add 40 μ l of master mixture to all wells.
- 2) Add 5 μ l of inhibitor solution of each well designated "Test Inhibitor". For the "Positive Control" and "Blank", add 5 μ l of the same solution without inhibitor (inhibitor buffer). Add 5 μ l of diluted **Trichostatin A** (100 μ M) to the wells designated "Inhibitor Control". Keep final DMSO concentration at or below 1%.
- 3) Add 5 μ l of **HDAC Assay Buffer** to the wells designated "Blank".
- 4) Initiate reaction by adding 5 μ l of diluted **HDAC4 enzyme** to the wells designated "Positive Control", "Test Inhibitor", and "Inhibitor Control". Incubate at 37°C for 30 min.

	“Blank”	Positive Control	Test Inhibitor	Inhibitor Control
HDAC substrate (200 μ M)	5 μ l	5 μ l	5 μ l	5 μ l
BSA (1 mg/ml)	5 μ l	5 μ l	5 μ l	5 μ l
HDAC Assay Buffer	35 μ l	30 μ l	30 μ l	30 μ l
Diluted Trichostatin A (100 μ M)	–	–	–	5 μ l
Test Inhibitor	–	–	5 μ l	–
Inhibitor buffer (no inhibitor)	5 μ l	5 μ l	–	–
Diluted HDAC4 (0.012 ng/ μ l)	–	5 μ l	5 μ l	5 μ l
Total	50 μl	50 μl	50 μl	50 μl


Step 2:

Add 50 μ l of undiluted **HDAC Developer (2x)** to each well. Incubate the plate at room temperature for 15 minutes.

Step 3:

Read sample in a microtiter plate-reading fluorimeter capable of excitation at a wavelength in the range of 350-380 nm and detection of emitted light in the range of 440-460 nm. “Blank” value is subtracted from all other values.

Example of Assay Results:



HDAC4 enzyme activity, measured using the *Fluorogenic HDAC4 Assay Kit*, #50064. Fluorescence was measured using a Bio-Tek fluorescent microplate reader. *Data shown is lot-specific.*

RELATED PRODUCTS:

<u>Product Name</u>	<u>Catalog #</u>	<u>Size</u>
HDAC1	50051	50 µg
HDAC2 (C-His)	50002	50 µg
HDAC2 (C-Flag)	50052	50 µg
HDAC3/NcoR2	50003	50 µg
HDAC4	50004	10 µg
HDAC5	50005	10 µg
HDAC6 (C-Flag)	50056	50 µg
HDAC6 (N-GST)	50006	50 µg
HDAC6 (H216A)	50046	50 µg
HDAC6 (H611A)	50066	50 µg
HDAC7	50007	10 µg
HDAC8	50008	50 µg
HDAC9	50009	10 µg
HDAC10	50010	50 µg
HDAC11	50011	50 µg
HDAC Assay Kit	50033	96 reactions
HDAC Assay Kit (Green)	50034	96 reactions
HDAC Class 2a Assay Kit	50041	96 reactions
HDAC2 Assay Kit	50062	96 reactions
HDAC3 Assay Kit	50073	96 reactions
HDAC5 Assay Kit	50065	96 reactions
HDAC6 Assay Kit	50076	96 reactions
HDAC7 Assay Kit	50077	96 reactions
HDAC8 Assay Kit	50068	96 reactions
HDAC9 Assay Kit	50069	96 reactions

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85