


Data Sheet

Fluorogenic HDAC1 Assay Kit

Catalog #: AMS.50061

DESCRIPTION: The Fluorogenic HDAC1 Assay Kit is a complete assay system designed to measure histone deacetylase 1 (HDAC1) activity for screening and profiling applications. It comes in a convenient 96-well format, with all the reagents necessary for 100 fluorescent HDAC1 activity measurements. In addition, the kit includes purified HDAC1 enzyme and a potent HDAC inhibitor, Trichostatin A, for use as a positive and negative control. The Fluorogenic HDAC1 Assay Kit is based on a unique fluorogenic substrate and developer combination. This assay method eliminates dealing with the radioactivity, extraction, and chromatography aspects of traditional assays. Using this kit, only two simple steps on a microtiter plate are needed to analyze the HDAC activity level. First, the HDAC fluorometric substrate, containing an acetylated lysine side chain, is incubated with purified HDAC1. The deacetylation sensitizes the substrate so subsequent treatment with the Lysine Developer produces a fluorophore that can then be measured using a fluorescence reader.

HDACs regulate cellular processes by catalyzing the hydrolysis of an acetyl group from acetyllysines in modified proteins. In the HDAC assay, fluorescent-dye molecules are attached to a peptide containing acetyllysine. Attachment to the peptide quenches the fluorescence of the dye. After treatment of the peptide with an HDAC, the reaction is mixed with a development solution that is specific for nonacetylated lysines. If the acetyl group has been removed from the lysine by the HDAC, this solution will release the dye allowing for fluorescence. Fluorescence is therefore directly related to HDAC activity.



OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

North America

23591 El Toro Rd, Suite #180
Lake Forest, CA 92630
Tel: +1 800 987 0985
Fax: +1 949 265 7703

amsbio

info@amsbio.com

www.amsbio.com
AMS Biotechnology

COMPONENTS:

Catalog #	Reagent	Amount	Storage	
50051	HDAC1 human recombinant enzyme	1 μ g	-80°C	Avoid Freeze/Thaw Cycles!
50037	Fluorogenic HDAC substrate (5 mM)	50 μ l	-80°C	
50030	2x HDAC Developer (contains Trichostatin A) (50 μ M)	6 ml	-80°C	
	Trichostatin A (200 μ M)	100 μ l	-20°C	
50031	HDAC assay buffer	10 ml	-20°C	
	black, low binding NUNC black microtiter plate	1 plate	Room temp.	

MATERIALS OR INSTRUMENTS REQUIRED BUT NOT SUPPLIED:

BSA (bovine serum albumin) (1 mg/ml)
Fluorescent microplate reader
Adjustable micropipettor and sterile tips
Rotating or rocker platform

APPLICATIONS: Great for studying enzyme kinetics and screening small molecular inhibitors for drug discovery and HTS applications.

STABILITY: One year from date of receipt when stored as directed.

REFERENCE(S):

1. Santo, L., et al., *Blood*. 2012 Mar 15; **119(11)**:2579-89.
2. Bradner, J.E., et al., *Nat Chem Biol*. 2010 Mar; **6(3)**: 238-243.

ASSAY PROTOCOL:**Immediately prior to assay:**

- 1) Dilute **HDAC substrate** 5 mM stock 25-fold with **HDAC assay buffer** to make a 200 μ M solution. (Make only sufficient quantity needed for the assay; store remaining 5 mM stock solution in aliquots at -80°C.)
- 2) Dilute **HDAC1** in **HDAC assay buffer** to 0.4 ng/ μ l (2 ng/reaction)*. Aliquot any remaining enzyme and store undiluted at -80°C. Keep diluted enzyme on ice. Discard any remaining diluted enzyme after use. *Note: optimal enzyme concentration may vary with the specific activity of the enzyme.

Step 1:

In duplicate, add the reaction mixtures (below) to the microtiter black plate as follows:

- 1) Prepare the master mixture: N wells \times (5 μ l HDAC substrate (200 μ M) + 5 μ l BSA (1 mg/ml) + 30 μ l HDAC assay buffer). Add 40 μ l of master mixture to all wells.
- 2) Add 5 μ l of inhibitor solution of each well designated "Test Inhibitor". For the "Positive Control", "Inhibitor Buffer Control" and "Blank", add 5 μ l of the same solution without inhibitor (inhibitor buffer). Add 5 μ l of Trichostatin A (200 μ M) to the well designated "TSA Inhibitor Control".

	"Blank"	Positive Control	Inhibitor Buffer Control	Test Inhibitor	TSA Inhibitor Control
HDAC substrate (200 μ M)	5 μ l				
BSA (1 mg/ml)	5 μ l				
HDAC assay buffer	35 μ l	30 μ l	30 μ l	30 μ l	30 μ l
Trichostatin A	-	-		-	5 μ l
Test Inhibitor	-	-	-	5 μ l	-
Inhibitor buffer (no inhibitor)	5 μ l	5 μ l	5 μ l	-	
HDAC1 (0.4 ng/ μ l)	-	5 μ l	5 μ l	5 μ l	5 μ l
Total	50 μl				

- 3) Add 5 μ l of HDAC assay buffer to the wells designated "Blank".
- 4) Initiate reaction by adding 5 μ l of diluted HDAC1 enzyme to the wells designated "Positive Control", "Inhibitor Buffer Control", "Test Inhibitor", and "TSA Inhibitor Control". Incubate at 37°C for 30 min.


Step 2:

Add 50 μ l of HDAC assay developer (2x) to each well. Incubate the plate at room temperature for 15 minutes.

Step 3:

Read sample in a microtiter plate-reading fluorimeter capable of excitation at a wavelength in the range of 350-380 nm and detection of emitted light in the range of 440-460 nm. "Blank" value is subtracted from all other values.

Example of Assay Results:



HDAC1 enzyme activity, measured using the HDAC1 Fluorogenic Assay Kit, #AMS.40061. Fluorescence was measured using a Bio-Tek fluorescent microplate reader. Data shown is lot-specific.

RELATED PRODUCTS:

<u>Product Name</u>	<u>Catalog #</u>	<u>Size</u>
HDAC1	50051	50 μg
HDAC2 (C-His)	50002	50 μg
HDAC2 (C-Flag)	50052	50 μg
HDAC3/NcoR2	50003	50 μg
HDAC4	50004	10 μg
HDAC5	50005	10 μg
HDAC6 (C-Flag)	50056	50 μg
HDAC6 (N-GST)	50006	50 μg
HDAC6 (H216A)	50046	50 μg
HDAC6 (H611A)	50066	50 μg
HDAC7	50007	10 μg
HDAC8	50008	50 μg
HDAC9	50009	10 μg
HDAC10	50010	50 μg
HDAC11	50011	50 μg
HDAC Assay Kit	50033	96 reactions

<i>HDAC Assay Kit (Green)</i>	50034	96 reactions
<i>HDAC Class 2a Assay Kit</i>	50041	96 reactions
<i>HDAC2 Assay Kit</i>	50062	96 reactions
<i>HDAC3 Assay Kit</i>	50073	96 reactions
<i>HDAC6 Assay Kit</i>	50076	96 reactions
<i>HDAC8 Assay Kit</i>	50068	96 reactions



UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

North America

23591 El Toro Rd, Suite #180
Lake Forest, CA 92630
Tel: + 1 800 987 0985
Fax: + 1 949 265 7703

amsbio

info@amsbio.com

www.amsbio.com
AMS Biotechnology