

Data Sheet

Fluorogenic SIRT3 Assay Kit Catalog #: AMS.50088

DESCRIPTION: The *Fluorogenic SIRT3 Assay Kit* is a complete assay system designed to measure Sirtuin 3 (SIRT3) activity for screening and profiling applications. It comes in a convenient 96-well format, with all the reagents necessary for **32 fluorescent SIRT3 activity measurements***. The kit includes purified SIRT3 enzyme and a SIRT inhibitor, Nicotinamide, for use as a positive and negative control, respectively. The *Fluorogenic SIRT3 Assay Kit* is based on a unique fluorogenic substrate and developer combination. This assay method eliminates dealing with the radioactivity, extraction, and chromatography aspects of traditional assays. Using this kit, only two simple steps on a microtiter plate are needed to analyze the SIRT3 activity level. First, the HDAC fluorometric substrate (HDAC substrate 1) is incubated with purified SIRT3 enzyme. The deacetylation sensitizes the substrate so subsequent treatment with the Lysine Developer produces a fluorophore that can then be measured using a fluorescence reader.

**Note: The kit includes sufficient HDAC substrate 1, buffers, NAD⁺, nicotinamide, and detection reagents for a full 96-well plate. Researchers can test their own SIRT3-containing samples in the remaining wells, or additional SIRT3 enzyme may be ordered separately (Cat. #50014).*

COMPONENTS:

Catalog #	Component	Amount	Storage	
AMS.50014	SIRT3 human recombinant enzyme	70 µg	-80°C	Avoid freeze/thaw cycles!
AMS.50032	Fluorogenic HDAC substrate 1 (5 mM)	50 µl	-80°C	
	Nicotinamide Adenine Dinucleotide (NAD ⁺) (50 mM)	50 µl	-80°C	
	Nicotinamide (10 mM)	500 µl	-80°C	
	2x SIRT Developer (contains 2 mM Nicotinamide)	6 ml	-80°C	
AMS.50090	SIRT assay buffer	10 ml	-20°C	
AMS.79685	Black, low binding NUNC black microtiter plate	1 plate	Room temp.	

APPLICATIONS: Great for studying enzyme kinetics and screening small molecular inhibitors for drug discovery and HTS applications.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, U.K.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Via Lisano 3
(CP.683)
6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

REFERENCES:

1. A. Ito *et al.* (2001) *EMBO J.* **20** 1331.
2. N.A. Barlev *et al.* (2001) *Mol. Cell* **8** 1243.
3. A. Ito *et al.* (2002) *EMBO J.* **21** 6236.

ASSAY PROTOCOL:

Immediately prior to assay:

- 1) Dilute HDAC substrate 5 mM stock 50-fold with SIRT assay buffer to make a 100 μ M solution. (Make only sufficient quantity needed for the assay; store remaining 5 mM stock solution in aliquots at -80°C.)
- 2) Dilute SIRT3 in SIRT assay buffer to 100 ng/ μ l (2000 ng/reaction)*. Aliquot any remaining enzyme and store undiluted at -80°C. Keep diluted enzyme on ice. Discard any remaining diluted enzyme after use. **Note: optimal enzyme concentration may vary with the specific activity of the enzyme.*

Step 1:

In duplicate, add the reaction mixtures (below) to the microtiter black plate. Incubate at 37°C for 60 min.

	Enzyme Positive Control	Inhibitor Negative Control	Test Inhibitor	“Blank” Negative Control
SIRT3 (100 ng/ μ l)	20 μ l	20 μ l	20 μ l	-
HDAC substrate (100 μ M)	5 μ l	5 μ l	5 μ l	5 μ l
BSA (1 mg/ml)	5 μ l	5 μ l	5 μ l	5 μ l
NAD ⁺ (50 mM)	0.5 μ l	0.5 μ l	0.5 μ l	0.5 μ l
Nicotinamide (10mM)	-	5 μ l	-	-
Test Inhibitor	-	-	X μ l	-
SIRT assay buffer	19.5 μ l	14.5 μ l	19.5 - X μ l	40 μ l
Total	50 μl	50 μl	50 μl	50 μl

Step2:

Add 50 μ l of SIRT assay developer (2x) to each well. Incubate the plate at room temperature for 15 minutes.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
 184 Park Drive, Milton Park
 Abingdon OX14 4SE, U.K.
 T: +44 (0) 1235 828 200
 F: +44 (0) 1235 820 482

 **North America**
 1035 Cambridge Street,
 Cambridge, MA 02141.
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

 **Germany**
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880

 **Switzerland**
 Via Lisano 3
 (CP.683)
 6900 Massagno
 T: +41 (0) 91 604 55 22
 F: +41 (0) 91 605 17 85

Step 3:

Read sample in a microtiter-plate reading fluorimeter capable of excitation at a wavelength in the range of 350-380 nm and detection of emitted light in the range of 440-460 nm. "Blank" value is subtracted from all other values.

RELATED PRODUCTS:

Fluorogenic SIRT1 Assay Kit	#AMS.50081	96 rxns.
Fluorogenic SIRT2 Assay Kit	#AMS.50082	96 rxns.
Chemiluminescent SIRT6 Assay Kit	#AMS.50086	96 rxns.
SIRT1 (Sir2) Enzyme	#AMS.50012	100 µg
SIRT2 Enzyme	#AMS.50013	100 µg
SIRT3 Enzyme	#AMS.50014	100 µg
SIRT4 Enzyme	#AMS.50015	100 µg
SIRT5 Enzyme	#AMS.50016	100 µg
SIRT6 Enzyme	#AMS.50017	100 µg
SIRT7 Enzyme	#AMS.50018	100 µg
SIRT Assay Developer	#AMS.50089	6 mL

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

AMSBIO | www.amsbio.com | info@amsbio.com

**UK & Rest of the World**

184 Park Drive, Milton Park
Abingdon OX14 4SE, U.K.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

**North America**

1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

**Germany**

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

**Switzerland**

Via Lisano 3
(CP.683)
6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85