

Data Sheet
ACE2 Inhibitor Screening Assay Kit
Catalog #AMS.79923
Size: 96 reactions

BACKGROUND: Angiotensin converting enzyme 2 (ACE2) is an exopeptidase that catalyzes the conversion of angiotensin II to angiotensin 1-7 and L-phenylalanine. Angiotensin II is part of the classical renin angiotensin system (RAS), a hormone system that regulates fluid balance, blood pressure and maintains vascular tone. ACE2 has been also proved to be the receptor for the human respiratory coronavirus NL63, the SARS-coronavirus (SARS-CoV) and the novel coronavirus 2019-nCoV/SARS-CoV-2.

DESCRIPTION: The *ACE2 Inhibitor Screening Assay Kit* is designed to measure the exopeptidase activity of ACE2 for screening and profiling applications. The ACE2 assay kit comes in a convenient 96-well format, with purified ACE2, its substrate, and ACE2 buffer for 96 reactions.

COMPONENTS:

Catalog #	Component	Amount	Storage	
11003	ACE2	2 µg	-80 °C	Avoid multiple freeze/thaw cycles!
	ACE2 Fluorogenic Substrate	3 ml	-80 °C	
	ACE2 Buffer	3 ml	-20 °C	
79685	96-well black microplate	1	Room Temp	

MATERIALS OR INSTRUMENTS REQUIRED BUT NOT SUPPLIED:

Adjustable micropipettor and sterile tips
 Fluorescent microplate reader

APPLICATIONS: Useful for studying enzyme kinetics and screening small molecule inhibitors and antibodies for drug discovery and HTS applications.

STABILITY: Up to 6 months from date of receipt, when stored as recommended.

REFERENCES:

Danilczyk, U., *et al. Circulation Research*. 2006. **98(4)**: 463-471.
 Ge, X. Y., *et al. Nature*. 2013. **503(7477)**: 535-538.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
 184 Park Drive, Milton Park
 Abingdon OX14 4SE, UK
 T: +44 (0)1235 828 200
 F: +44 (0) 1235 820 482



North America
 1035 Cambridge Street,
 Cambridge, MA 02141
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218



Germany
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880



Switzerland
 Centro Nord-Sud 2E
 CH-6934 Bioggio-Lugano
 T: +41(0) 91 604 55 22
 F: +41(0) 91 605 17 85

ASSAY PROTOCOL:***All samples and controls should be tested in duplicate.***

1. Thaw **ACE2** on ice. Upon first thaw, briefly spin tube containing enzyme to recover full contents of the tube. Note: **ACE2** is sensitive to freeze/thaw cycles. Avoid multiple freeze/thaw cycles. Do not re-use thawed aliquots or diluted enzyme.
2. Prepare **Enzyme solution** (0.5 ng/μl ACE2) by diluting ACE2 in **ACE2 Buffer**.
3. Add 20 μl of **Enzyme solution** (0.5 ng/μl ACE2) to each well designated "Positive Control" and "Test Inhibitor," and 20 μl of **ACE2 buffer** to each well designated "Blank."
4. Add 5 μl of **Test Inhibitor** solution to each well designated "Test Inhibitor." For the wells labeled "Positive Control" and "Blank," add 5 μl of 10% DMSO in water (Inhibitor buffer). Note: Keep the final concentration of DMSO in the reaction ≤1%.
5. Add 25 μl of **ACE2 Fluorogenic Substrate** to all wells. Protect from light and incubate reaction at room temperature for 60 minutes.

	Positive Control	Test Inhibitor	Blank
Enzyme solution (0.5 ng/μl ACE2)	20 μl	20 μl	-
ACE2 Buffer	-	-	20 μl
Test inhibitor	-	5 μl	-
10% DMSO in water (Inhibitor buffer)	5 μl	-	5 μl
ACE2 Fluorogenic Substrate	25 μl	25 μl	25 μl
Total	50 μl	50 μl	50 μl

6. Read fluorescence intensity of the samples ($\lambda_{\text{excitation}} = 555 \text{ nm}$; $\lambda_{\text{emission}} = 585 \text{ nm}$) in an appropriate microplate reader.

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218




Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85

Example of assay results:



ACE2 inhibition by DX600 (Cayman, #22186) and MLN-4760 (Sigma, #5306160001) measured using the *ACE2 Inhibitor Screening Assay Kit*, AMSBIO #AMS.79923. Fluorescence was measured using a Bio-Tek microplate reader. Data shown is lot-specific. For lot-specific information, please contact AMSBIO at info@amsbio.com.

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
 184 Park Drive, Milton Park
 Abingdon OX14 4SE, UK
 T: +44 (0)1235 828 200
 F: +44 (0) 1235 820 482

 **North America**
 1035 Cambridge Street,
 Cambridge, MA 02141
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

 **Germany**
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880

 **Switzerland**
 Centro Nord-Sud 2E
 CH-6934 Bioggio-Lugano
 T: +41(0) 91 604 55 22
 F: +41(0) 91 605 17 85