

Data Sheet

ACE2 – HEK293 Recombinant Cell Line

Catalog #AMS.79951

Description

Recombinant clonal stable HEK293 cell line constitutively expressing full length human ACE2, Genbank #NM_021804.3). Surface expression of ACE2 was confirmed by flow cytometry.

Background

Human Angiotensin converting enzyme 2 (ACE2), also known as ACEH, is an integral membrane protein found in the outer space of cells in the lungs, arteries, heart, kidney, and intestines. ACE2 serves as the entry point into cells for some coronaviruses, including the SARS-CoV-2 virus that is responsible for the COVID-19 pandemic.

Application

This cell line is useful for ACE2 binding assays, flow cytometry, or for screening ACE2 antibodies.

Format

Each vial contains ~ 2 x 10⁶ cells in 1 ml of 10% DMSO in FBS.

Storage

Store in liquid nitrogen immediately upon receipt.

Mycoplasma Testing

This cell line has been screened using the MycoAlert™ Mycoplasma Detection Kit (Lonza, #LT07-118) to confirm the absence of Mycoplasma contamination.

Cell Culture

Thaw Medium 1 (#AMS.60187): MEM medium (Hyclone #SH30024.01) supplemented with 10% FBS, 1% non-essential amino acids (Hyclone #SH30238.01), 1 mM Na pyruvate (Hyclone #SH30239.01), 1% Penicillin/Streptomycin (Hyclone #SV30010.01).

Growth Medium 1N (#AMS.79801): Thaw Medium 1 (#AMS.60187) and 0.5 µg/ml of Puromycin (InvivoGen, #ant-pr-1).

Cells should be grown at 37°C with 5% CO₂ using Growth Medium 1N.

To thaw the cells, it is recommended to quickly thaw the frozen cells from liquid nitrogen in a 37°C water-bath, and transfer to a tube containing 10 ml of Thaw Medium 1 (no puromycin). Spin down cells, resuspend cells in pre-warmed Thaw Medium 1 (no puromycin), transfer resuspended cells to a T25 flask and culture in 37°C CO₂ incubator. At first passage switch to

Growth Medium 1N (contains puromycin). Cells should be split before they reach complete confluence.

To passage the cells, rinse cells with phosphate buffered saline (PBS), detach cells from culture vessel with Trypsin/EDTA, add Growth Medium 1N and transfer to a tube. Spin down cells, resuspend cells and seed appropriate aliquots of cell suspension into new culture vessels. Subcultivation ration: 1:5 to 1:10 weekly or twice a week.

To freeze down the cells, rinse cells with phosphate buffered saline (PBS), and detach cells from culture vessel with Trypsin/EDTA. Add Growth Medium 1N and transfer to a tube, spin down cells, and resuspend in freezing medium (10% DMSO + 90% FBS) at $\sim 2 \times 10^6$ cells/ml. Dispense 1 ml of cell aliquots into cryogenic vials. Place vials in an insulated container for slow cooling and store at -80°C overnight. Transfer to liquid nitrogen the next day for storage. It is recommended to expand the cells and freeze down more than 10 vials of cells for future use at early passage.

Sequence

MSSSSWLLLSLVAVTAAQSTIEEQAKTFLDKFNHEAEDLFYQSSLASWNYNTNITEENVQNMN
NAGDKWSAFLKEQSTLAQMYPLQEIQNLTVKQLQLALQQNGSSVLSEDKSKRLNTILNTMSTIY
STGKVCNPDNPQECLLLEPGLNEIMANSLDYNERLWAWESWRSEVGKQLRPLYEEYVVLKNE
MARANHYEDYGDYWRGDYEVNGVDGYDYSRGLIEDVEHTFEEIKPLYEHLHAYVRAKLMNA
YPSYISPIGCLPAHLLGDMWGRFWTNLYSLTVFPGQKPNIDVTDAMVDQAWDAQRFKEAEKF
FVSVGLPNMTQGFWENSMLTDPGNVQKAVCHPTAWDLGKGDFRILMCTKVTMDDFLTAHHE
MGHIQYDMAYAAQPFLLRNGANEGFHEAVGEIMSLSAATPKHLKSIGLLSPDFQEDNETEINFL
LKQALTIVGTLPFTYMLEKWRWMVFKGEIPKDQWMKKWEMKREIVGVVEPVPHDETYCDPA
SLFHVSNDYSFIRYYTRTLYQFQFQEALCQAAKHEGPLHKCDISNSTEAGQKLFNMLRLGKSEP
WTLAENVVGAKNMNVRPLLNYFEPLFTWLKDQNKNSFVGWSTDWSPYADQSIKVRISLKSAL
GDKAYEWNNDNEMYLFRSSVAYAMRQYFLKVKNQMILFGEEDVRVANLKPRISFNFFVTAPKNV
SDIIPRTEVEKAIRMSRSRINDAFRLNDNSLEFLGIQPTLGPPNQPPVSIWLIVFGVVMGVIVVGIV
ILIFTGIRDRKKKNKARSGENPYASIDISKGENNPGFQNTDDVQTSF

Materials Required but Not Supplied

- HEK293 growth medium or use Thaw Medium 1 (#AMS.60187): MEM medium (Hyclone #SH30024.01) supplemented with 10% FBS, 1% non-essential amino acids (Hyclone #SH30238.01), 1mM Na pyruvate (Hyclone #SH30239.01), 1% Penicillin/Streptomycin (Hyclone#SV30010.01).
- Growth Medium 1N (#AMS.79801): Thaw medium 1 (#AMS.60187) and 0.5 µg/ml of Puromycin (InvivoGen, #ant-pr-1).
- Spike (SARS-CoV-2) Pseudotyped Lentivirus (Luc Reporter) (#AMS.79942)
- Bald lentiviral pseudovirion (Luciferase reporter) (#AMS.79943)
- Anti-SARS-CoV-2 Spike antibody (clone AM001414, #91361)
- Anti-ACE2 antibody (#AF933)
- Recombinant ACE2 protein (#AMS.11003)
- 96-well tissue culture treated, white clear-bottom assay plate (Corning, #3610)
- ONE-Step™ luciferase assay system (#AMS.60690)

Assay Protocol

The following protocol is a general guideline for transducing ACE2-HEK293 cells using SARS-CoV-2 Spike pseudotyped lentivirus (Luciferase reporter). The optimal transduction conditions (e.g. MOI, concentration of polybrene, time of assay development) should be optimized according to the cell type and the assay requirements.

1. Day 1: Harvest ACE2-HEK293 cells from culture and seed cells at a density of 5,000-10,000 cells per well into white opaque 96-well microplate in 50 µl of Thaw Medium 1 (#AMS.60187). Incubate cells at 37°C with 5% CO₂ overnight.
2. Day 2: prepare serial dilutions of anti-Spike or anti-ACE2 antibody in Thaw Medium 1.

To test anti-Spike antibody, preincubate 5 µl of the SARS-CoV-2 Spike pseudotyped lentivirus with 5 µl of diluted anti-Spike antibody for 30 minutes. After incubation, add 10 µl of virus/antibody mix into each well of the ACE2-HEK293 cells.

To test anti-ACE2 antibody, add 5 µl of diluted anti-ACE2 antibody into each well of ACE2-HEK293 cells and incubate for 30 minutes. At the end of the incubation, add 5 µl of SARS-CoV-2 Spike pseudotyped lentivirus into each well.

For control wells, the same number of ACE2-HEK293 cells were seeded, but no virus or antibody was added.

Incubate the plates at 37°C with 5% CO₂ overnight.

3. Day 3: Remove the medium containing the lentivirus from the wells. Add 50 µl of fresh Thaw Medium 1 to each well.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

If the tested antibody does not adversely affect the target cells, it is not necessary to change the medium on Day 3.

- Day 4, approximately 48-60 hours after transduction, prepare the ONE-Step™ Luciferase reagent per recommended protocol. Add 50 µl of ONE-Step™ Luciferase Assay reagent per well. Incubate at room temperature for ~15 to 30 minutes and measure luminescence using a luminometer. The transduction efficacy was determined by measuring the luciferase activity.




Figure 1. Expression of ACE2 validated by flow cytometry. ACE2-HEK293 cells (green) or parental HEK293 cells (blue) were stained with anti-human ACE2 polyclonal goat IgG primary antibody (#AF933) and Alexa Fluor 488 conjugated rabbit anti-goat IgG secondary antibody (Thermo Fisher #A-21222). The ACE2 expression was analyzed by FACS.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.




Figure 2. Transduction of ACE2-HEK293 Cells using SARS-CoV-2 Spike Pseudotyped Lentivirus. A. Approximately 10,000 cells/well of ACE2-HEK293 cells or HEK293 parental cells were transduced with 5 μ l/well of SARS-CoV-2-Spike pseudotyped lentivirus (Luc reporter) #AMS.79942). After 18 hours of transduction, the medium was changed to fresh HEK growth medium (Thaw Medium 1). After 48 hours of transduction, ONE-Step Luciferase reagent (#AMS.60690) was added to cells to measure the luciferase activity. The SARS-CoV-2 Spike pseudotyped lentivirus transduced ACE2-HEK293 cells with much greater efficiency compared with HEK293 parental cells, indicating the transduction is dependent upon ACE2 expression. **B.** Approximately 10,000 ACE2-HEK293 cells/well were transduced with 10 μ l/well of SARS-CoV-2 Spike pseudotyped lentivirus (Luc reporter) mixed with anti-Spike antibody (clone #AM001414, #91361) or recombinant ACE2 (#AMS.11003). After 18 hours of transduction, the medium was changed to fresh HEK growth medium (Thaw Medium 1). After 48 hours of transduction, ONE-Step Luciferase reagent (#AMS.60690) was added to cells to measure the luciferase activity.

License Disclosure: Purchase of this cell line grants you with a 10-year license to use this cell line in your immediate laboratory, for research use only. This license does not permit you to share, distribute, sell, sublicense, or otherwise make the cell line available for use to other laboratories, departments, research institutions, hospitals, universities, or biotech companies. The license does not permit the use of this cell line in humans or for therapeutic or drug use. The license does not permit modification of the cell line in any way. Inappropriate use or distribution of this cell line will result in revocation of the license and result in an immediate cease of sales and distribution of BPS products to your laboratory. BPS does not warrant the suitability of the cell line for any particular use, and does not accept any liability in connection with the handling or use of the cell line. Modifications of this cell line, transfer to another facility, or commercial use of the cells may require a separate license and additional fees; contact sales@bpsbioscience.com for details. Publications using this cell line should reference BPS Bioscience, Inc., San Diego.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

Related Products

Product

Product	Cat. #	Size
ACE2 CHO Recombinant Cell Line	AMS.79959	2 vials
ACE2 HeLa Recombinant Cell Line	AMS.79958	2 vials
ACE2 Lentivirus	AMS.79944	2 vials
SARS-CoV-2 Spike Pseudotyped Lentivirus (Luciferase Reporter)	AMS.79942	500 µl x2
Bald Lentiviral Pseudovirion (Luciferase Reporter)	AMS.79943	500 µl x2
ACE2, His-tag	AMS.11003-2	100 µg
Thaw Medium 1	AMS.60187	100 ml
Growth Medium 1N	AMS.79801	500 ml

AMS BIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85