

# iMatrix-511MG

Product No. AMS.892 005 1,050 µg

Version 003  
Store at 2-15°C

**Product description:** iMatrix-511MG is a recombinant human laminin-511 E8 fragment protein expressed in Chinese Hamster Ovary (CHO)-S cells. iMatrix-511MG contains the integrin-binding site of the laminin-511 molecule. iMatrix-511MG is a useful cell culture substrate for feeder-free culture and single-cell passage of ES cells and iPS cells, facilitating stable culture expansion. iMatrix-511MG is also useful for the culture of other cells adhering to laminin-511.

**Content:** Recombinant human laminin-511 E8 fragment protein in PBS(-)

**Concentration:** 0.5 mg/mL

**Amount:** Product No. 892 005 1,050 µg / 6 tubes  
175 µg / 0.35 mL / tube

**Storage:** Store at 2°C to 15°C, protect from light.

**Expiration date:** The shelf life is 2 years from the date of manufacture. The expiration date is printed on the outer carton.

**Activity:** The dissociation constant (Kd) for the binding with integrin  $\alpha 6\beta 1$  is 1 nM or less.

**Methods of use:** By either of the following methods, iMatrix-511MG can be coated onto a culture vessel. **The optimum coating density may differ by cell-type, cell-line, medium selected, or purpose.** Insufficient coating density may result in the detachment of cells and varied cell conditions while the excessive coating density may lead to difficulty in detaching cells for passage.

## A. Pre-coating method

Determine the optimal coating density. 0.5 µg/cm<sup>2</sup> is a standard but test between 0.1 and 1.5 µg/cm<sup>2</sup>.

- 1) Dilute iMatrix-511MG with PBS(-). Use the diluted iMatrix-511MG immediately. To coat with 0.5 µg/cm<sup>2</sup> onto a 6-well plate with 9.6 cm<sup>2</sup>/well, dilute 9.6 µL of iMatrix-511MG with 2 mL of PBS(-) per well.
- 2) Place the diluted iMatrix-511MG into a culture vessel and incubate either at 37°C for 1 h, or at room temperature for 3 h, or at 4°C overnight.
- 3) Aspirate the coating solution. Then, immediately seed your cells. **Do not allow the coated surface to dry.**

## B. Pre-mixing method

Determine the optimal coating density for cell culture. The standard density is 0.25 µg/cm<sup>2</sup> but test between 0.1 and 1.5 µg/cm<sup>2</sup>. The optimal coating density may be affected by the medium and cell density of the cell suspension.

- 1) Add iMatrix-511MG to the cell suspension. To coat with 0.25 µg/cm<sup>2</sup> onto a 6-well plate with 9.6 cm<sup>2</sup>/well, add 4.8 µL of iMatrix-511MG to 2 mL of the cell suspension per well.
- 2) Place the cell suspension containing iMatrix-511MG into a culture vessel.

\*If you face difficulties in detaching cells for passage, re-adjust the conditions (e.g., reduce the coating density).

## References:

- Taniguchi Y. *et al.* (2009), *J. Biol. Chem.* **284** (12): 7820-31  
Miyazaki T. *et al.* (2012), *Nat. Commun.* **3**: 1236  
Nakagawa M. *et al.* (2014), *Sci. Rep.* **4**: 3594  
Takashima Y. *et al.* (2014), *Cell* **158** (6): 1254-69  
Miyazaki T. *et al.* (2017), *Sci. Rep.* **7**: 41165  
Kikuchi T. *et al.* (2017), *Nature* **548** (7669): 592-6  
Goparaju S.K. *et al.* (2017), *Sci. Rep.* **7**: 42367  
Hayashi R. *et al.* (2017), *Nat. Protoc.* **12** (4): 683-96  
Ishii K. *et al.* (2018), *Stem Cell Reports* **10** (2): 568-82

**Caution:** In the event of accidental ingestion or contact with the eyes, immediately wash the affected area and seek medical attention.

**Product information:** Current information including references and FAQs are available online at [www.amsbio.com/imatrix-recombinant-laminin-511/](http://www.amsbio.com/imatrix-recombinant-laminin-511/).

**Designed by:** Matrixome Inc.

3-2 Yamadaoka, Suita, Osaka 565-0871, Japan  
Institute for Protein Research, Osaka University

**Manufactured by:** Nippi, Incorporated

1-1-1 Senju Midori-cho, Adachi, Tokyo 120-8601, Japan

AMS BIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

AMS BIO LLC  
USA & Canada   
1035 Cambridge Street,  
Cambridge, MA 02141  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218

AMS BIO Europe BV  
EU   
Berenkoog 41,  
1822 BH Alkmaar,  
Netherlands  
T: +31 (0) 72 8080244  
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd  
UK & Rest of the World   
184 Park Drive, Milton Park  
Abingdon OX14 4SE  
T: +44 (0) 1235 828 200  
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd  
Switzerland   
Via Lisano 3,  
(CP.683)  
CH-6900 Massagno  
T: +41 (0) 91 604 55 22  
F: +41 (0) 91 605 17 85

AMS BIO Europe BV  
Deutschland   
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880