

For research use only  
Not for use in diagnostic procedures

# iMatrix-332

Product No. AMS.892 031 350 µg  
Product No. AMS.892 032 1,050 µg

Version 001  
Store at 2-15°C

**Product description:** iMatrix-332 is a recombinant human laminin-332 E8 fragment protein expressed in Chinese Hamster Ovary (CHO)-S cells. iMatrix-332 contains the integrin-binding site of the laminin-332 molecule. iMatrix-332 is useful for differentiation of iPSC and ESCs into corneal epithelial cells. iMatrix-332 also supports the culture of the other cells adhering to laminin-332 such as keratinocytes.

**Content:** Recombinant human laminin-332 E8 fragment protein in 20 mM phosphate buffer, 500 mM NaCl

**Protein concentration:** 0.5 mg/mL

**Amount:** 175 µg / 0.35 mL / tube  
Product No. 892 031 350 µg / 2 tubes  
Product No. 892 032 1,050 µg / 6 tubes

**Storage:** Store at 2°C to 15°C, protect from light.

**Expiration date:** The shelf life is 2 years from the date of manufacture. The expiration date is printed on the carton.

**Methods of use:** By the following method iMatrix-332 can be coated onto a culture vessel. **The optimum coating density may differ by cell-type, cell-line, medium selected, or purpose.** Insufficient coating density may result in the detachment of cells and varied cell conditions while the excessive coating density may lead to difficulty in detaching cells for passage.

## Coating protocol

Determine the optimal coating density. 0.5 µg/cm<sup>2</sup> is a standard but test between 0.1 and 1.5 µg/cm<sup>2</sup>.

- 1) Dilute iMatrix-332 with PBS(-). Use the diluted iMatrix-332 immediately. To coat with 0.5 µg/cm<sup>2</sup> onto a 6-well plate with 9.6 cm<sup>2</sup>/well, dilute 9.6 µL of iMatrix 332 with 2 mL of PBS( ) per well.
- 2) Place the diluted iMatrix-332 into a culture vessel and incubate either at 37°C for 1 h, or at room temperature for 3 h, or at 4°C overnight.
- 3) Aspirate the coating solution. Then, immediately seed your cells. **Do not allow the coated surface to dry.**

\*If you face difficulties in detaching cells for passage, re-adjust the conditions (e.g., reduce the coating

density).

\*Coating protocol is illustrative only.

## References:

Nishiuchi R. *et al.* (2006) *Matrix Biol.* **25** (3), 189-97  
Miyazaki T. *et al.* (2012) *Nat. Commun.* **3**, 1236  
Shibata S. *et al.* (2018) *Cell Rep.* **25** (6), 1668-79  
Shibata S. *et al.* (2020) *Stem Cell Reports* **14** (4), 663-76

**Caution:** For research use only. Not intended for human use. In the event of accidental ingestion or contact with the eyes, immediately wash the affected area and seek medical attention.

**Product information:** Current information including references and Q&A are available on the website of Matrixome, Inc. Please use the URL or QR code below.

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

UK & Rest of the World   
184 Park Drive, Milton Park  
Abingdon OX14 4SE.  
T: +44 (0) 1235 828 200  
F: +44 (0) 1235 820 482

North America   
1035 Cambridge Street,  
Cambridge, MA 02141.  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218

Europe   
Berenkoog 41,  
1822 BH Alkmaar,  
Netherlands  
T: +31 (0) 72 8080244  
F: +31 (0) 72 8080142

Switzerland   
Via Lisano 3,  
(CP.683)  
CH-6900  
T: +41 (0) 91 604 55 22  
F: +41 (0) 91 605 17 85

Germany   
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880