

CellScrub™

amsbio

Washing Buffer

Cat. #	Contents	Quantity
B100001	CellScrub Washing Buffer	1 x 100 ml


Shipping	Shipped on blue ice
Storage	Store CellScrub Washing Buffer at 4°C.

Related Products	Catalog #
GenePORTER® Transfection Reagent, 75 rxns.	T201007
GenePORTER® 2 Transfection Reagent, 75 rxns.	T202007
GeneSilencer® siRNA Transfection Reagent, 200 rxns.	T55500750
BioPORTER® Protein Delivery Reagent, 24 rxns.	BP502401
Cytofectin™ Transfection Reagent, 1 ml	T610001
Perfectin™ Transfection Reagent, 100 rxns.	T303007
BoosterExpress Reagent, 3 x 1.5 ml	T20100B

Introduction: CellScrub Buffer is a unique washing reagent designed to remove all cationic lipid/DNA complexes associated with cell surfaces after transfection. The CellScrub Buffer is non-toxic for the cells and allows discrimination between extracellular and intracellular plasmid DNA (see Figure 1 on the right). The use of the CellScrub Buffer as part of a cell washing procedure allows more precise quantification of the amount of DNA delivered into the cells for DNA stability or cytometry studies.

Figure 1. Removal of Extracellular Cationic lipid/DNA Complexes By the CellScrub Washing Buffer.

Adherent cells were incubated with cationic lipid/DNA complexes at 4°C for 1-2 hours. Several washing procedures were used to remove cell-associated complexes. After washes, cells were lysed, the DNA was extracted and analyzed by agarose gel electrophoresis. DNA was stained with ethidium bromide. Arrows show plasmid DNA. **Lane 1:** Trypsinisation and PBS washes. **Lane 2:** Washing procedure including the CellScrub Washing Buffer.



METHODS AND PROCEDURES

A. Example Protocol: DNA Stability Study

1. Transfect adherent cells cultured in 6-well plates with 5µg of DNA complexed with cationic lipids per well as described*.
2. After the desired times of incubation, remove culture medium and wash cells once with PBS containing calcium chloride and magnesium chloride.
3. Treat the cells at room temperature for 10-15 minutes with the CellScrub Washing Buffer (1-2 ml/well) to remove all extracellular cationic lipid/DNA complexes.
4. After two more washes with PBS (without calcium chloride and magnesium chloride), trypsinize the cells, then wash them once again with PBS.
5. Centrifuge and lyse cells as desired. Perform DNA extractions, agarose gel electrophoresis, or Southern blot analysis according to standard procedures.

NOTE: Intracellular plasmid delivery and transgene expression can be measured simultaneously by flow cytometry using the **pGeneGrip™ vector** (for example, Rhodamine/GFP; Cat # G101045 where rhodamine is the intracellular DNA marker and green fluorescent protein the transgene expression marker). Cell transfection, fixation and flow cytometry analysis can be done according to Tseng *et al.* **
In order to insure that the fluorescence is derived from intracellular plasmid, the cell washing procedure described above for the DNA stability study must be applied.

*Felgner, J.H. *et al.* Cationic lipid-mediated transfection in mammalian cells: "Lipofection". *J. Tiss. Cult. Meth.* **15**, 63-68 (1993).

Tseng, W.-C. *et al.* Transfection by cationic liposomes using simultaneous single cell measurements of plasmid delivery and transgene expression. *J. Biol. Chem.* **272, 25641-25647 (1997).

LIMITED LICENSE: The purchase price paid for the CellScrub™ Washing Buffer grants end users a non-transferable, non-exclusive license to use this buffer for **internal research use only** as described in this manual; in particular, research use only excludes and without limitation, resale, repackaging, or use for the making or selling of any commercial product or service without the written approval of Genlantis, a division of Gene Therapy Systems, Inc. (GTS) -- separate licenses may be available for non-research use or applications. CellScrub Washing Buffer is not to be used for human diagnostic or included/used in any drug intended for human use. Care and attention should be exercised in handling this reagent by following appropriate research laboratory practices. Purchasers may refuse this license by returning the enclosed materials unused. By keeping or using the enclosed materials, you agree to be bound by the terms of this license. The laws of the State of California shall govern the interpretation and enforcement of the terms of this License.

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada 
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland 
Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland 
T: +49 (0) 69 779099
F: +49 (0) 69 13376880