

BONE RESORPTION ASSAY PLATE 24

BONE RESORPTION ASSAY PLATE 48

This product is a calcium phosphate (CaP)-coated 24-well or 48-well plate used to measure the bone resorption activity of osteoclasts. The plate is coated with a synthetic CaP (carbonate apatite), similar to that of natural apatite, which is able to be used as an alternative to dentin discs.

1) Specifications

(A) BONE RESORPTION ASSAY PLATE 24:

A 24-well plate coated with a synthetic CaP (carbonate apatite) (1 plate, gamma sterilized)

(B) BONE RESORPTION ASSAY PLATE 48:

A 48-well plate coated with a synthetic CaP (carbonate apatite) (1 plate, gamma sterilized)



2) Example of use

- (1) **A:** Wash each well of the 24-well plate with 1 mL of culture medium.
B: Wash each well of the 48-well plate with 0.5mL of culture medium.
- (2) Inoculate RAW264 or RAW264.7 cells into each well in culture medium (DMEM/F-12 or α MEM containing 10% FBS). Add an inducer of osteoclastic differentiation, such as RANKL (100 ng/mL) and the test substances to be evaluated.
A: 24-well plate 1×10^4 cells/mL; 1 mL/well
B: 48-well plate 5×10^3 cells/mL; 0.5mL/well
- (3) On day 3, change the medium with freshly made medium (containing RANKL and drugs). This step may be eliminated. However, the induction of osteoclastic differentiation by RANKL would be reduced.
- (4) On day 5 or 6, remove the conditioned medium from each well and treat the wells with 5% sodium hypochlorite for five minutes. After washing the plate with water and then drying it, photograph the regions in each well using a microscope and measure the pit area with image analyzing software.

3) Assay precautions

- (1) To stimulate cells, we recommend a RANKL concentration ≥ 100 ng/mL.
- (2) This product is for research use only, and not for use in diagnostic or therapeutic procedures.

4) Expected Results



UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880


North America

23591 El Toro Rd, Suite #180
Lake Forest, CA 92630
Tel: +1 800 987 0985
Fax: +1 949 265 7703

amsbio

info@amsbio.com


www.amsbio.com
AMS Biotechnology



- A. A phase-contrast micrograph of RAW264 cells (day 6) cultured in CaP-coated plates stimulated with RANKL (Oriental Yeast Co., Ltd., Tokyo, Japan; 100 ng/mL). Osteoclast-like cells (OC) were observed.
- B. Photograph of the plate after removing cells. Pits can be observed macroscopically (Left: without RANKL; Right: with RANKL).
- C. Micrograph of the pits in a CaP-coated plate (with RANKL).
- D. RANKL-dependent increase of the pit area (mean \pm S.D., $n = 3$).
- E. The inhibitory effect of the bisphosphonate, Pamidronate, on CaP resorption induced by RANKL (100 ng/mL).

#20131009

For research use only. Not for clinical diagnosis.



UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

North America

23591 El Toro Rd, Suite #180
Lake Forest, CA 92630
Tel: +1 800 987 0985
Fax: +1 949 265 7703

amsbio

info@amsbio.com

www.amsbio.com
AMS Biotechnology