

FITC-dextran

(Fluorescein isothiocyanate dextran)

Chemical names:

- Dextran(3',6'-dihydroxy-3-Oxospiro (isobenzofuran-1- (3H), 9'-[9Hxanthen]-5 (or 6)-ylcarbamo)thiate.
- Fluoresceinisothiocyanate- dextran
- Fluoresceinyl thiocarbamoyl- dextran

CAS number: 60842-46-8

Structure:




Fig.1. Structural representation of fragment of FITC-dextran molecule.

Properties

FITC-dextran is supplied as a yellow/orange powder which dissolves freely in water or salt solutions giving a yellow solution. The product also dissolves in DMSO, formamide and other polar organic solvents but is insoluble in lower aliphatic alcohols, acetone, chloroform and dimethylformamide.

Spectral data

Excitation is best performed at 493 nm and fluorescence measured at 518 nm (Fig. 2). Since the charge status of the fluorescein moiety is dependent on the pH and ionic strength of the medium, the fluorescence intensity will also vary with these parameters. The maximal intensity is observed at pH > 8. Measurements in biological media will significantly affect the fluorescence intensity which may be enhanced or depressed.

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada 
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland 
Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland 
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

©AMSBIO / AMS Biotechnology (Europe) Ltd. AMSBIO is the global source of all the products listed here.
All trademarks mentioned here are the property of their respective owners.




Fig. 2. Fluorescence scan of FITC-dextran 70 in phosphate buffered saline pH 9.




Fig. 3. Fluorescence (Emission 520nm) of FITC-dextran in the range pH 4-9.

Storage and stability

The stability of FITC-dextran has been investigated in various media and at various temperatures. From these studies, it is concluded that the stability of FITC-dextran in vitro and in vivo is excellent. Only at elevated pH (>9) and elevated temperatures is there a risk for hydrolysis of the fluorescein label. Studies at 37°C in rabbit plasma, muscle homogenate, liver homogenate and urine established that FITC-dextran is stable for at least 3 days. No changes in the mol. wt and no release of fluorescein moieties was noted. FITC-dextran was stable in 6 % trichloroacetic acid at room temperature for 3 days. The hydrolysis shows specific catalysis by hydroxide ions in the pH range 10-10.75 (1).

Hydrolysis of the thiocarbamoyl linkage gives rise to 4- or 5-aminofluorescein which is readily determined by HPLC. In unpublished studies, the stability of an autoclaved FITC-dextran 70 solution was studied at temperatures from 8 to 50°C over a period of 5 months. Only at 50°C could a slight increase (1%) in free aminofluorescein be noted. Autoclaving alone gives a 2.7 % release of free aminofluorescein. In other unpublished studies, FITC-dextran was found to be stable in solution at pH 4 for up to 1 month at temperatures up to 35°C. At 80°C and pH 4, the thiocarbamoyl linkage was stable for 30 min. However, the dextran may degrade. At pH 9, considerable (24 %) decrease in fluorescence took place at 35 °C over 1 month. Several studies have confirmed the in vivo stability of FITC-dextran during the duration of the experiments (2).

Toxicity

In studies in mice, FITC-dextran was found to be tolerated well when injected intravenously or intraperitoneally in doses up to 6 g/kg body-weight. Their toxicity patterns follow those of parent dextrans. Clinical dextran fractions have been employed for over 50 years as plasma volume expanders. Dextran-induced anaphylactoid reactions (DIARs) have been observed

in humans after injection of clinical dextran solutions (3,4). FITC-dextrans are also likely to display this type of behavior but few reports of problems with experimental animals have appeared.

Synthesis

Selected dextran fractions prepared from native Dextran B512F are labelled with fluorescein (5). The fluorescein moiety is attached by a stable thiocarbamoyl linkage and the labelling procedure does not lead to any depolymerization of the dextran. The FITC-dextrans have from 0.002-0.008 mol FITC per glucose unit and at these low levels of substitution confer minimal charges to the dextran, which is an essential requirement for permeability studies. Physical chemical properties of FITC-dextran. The dextran molecule at molecular weights greater than 5000 Daltons behaves as a flexible and extended coil in solution. Table 1. (below) shows the molecule dimensions at various molecular weights.

Dextran MW	Stokes radius (Å)	Radius of gyration (Å)
2 x 10 ⁶	270	380
1 x 10 ⁶	199	275
500 000	147	200
200 000	130	130
100 000	69	95
70 000	58	80
40 000	44.5	62
10 000	23.6	-

Table 1. Molecular dimensions of dextran

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada 
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland 
Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland 
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

Dextrans and FITC-dextrans will exhibit Newtonian flow characteristics i.e. the viscosity is independent of shear rate (Fig. 4). Studies in the range pH 4-10 establish that the viscosity is independent of pH. The isoelectric point of FITC-dextran lies in the range 8-9 and between pH 6.5 and 9.5 they show no migration on electrophoresis (1).




Fig. 4. The viscosities of dextran fractions at various concentrations.

Applications

FITC-dextrans are primarily used for studying permeability and transport in cells and tissues. An added benefit is that measurements of the fluorescence provide quantitative data on the permeability of healthy and diseased tissues. Such studies can be performed in real time by intravital fluorescence microscopy. The technique offers high sensitivity and concentrations down to $1\mu\text{g/ml}$ can be detected in tissue fluids. FITC-dextrans have also been used as a pH probe in cells (6,7). It may also be noted from polarization experiments that the rotational freedom of fluorescein conjugated to dextran remains high and fluorescent lifetime of the excited state is similar to that before conjugation (6).

1. General procedures

The microvasculature of the hamster cheek pouch has proved to be a useful model for studying plasma leakage in different experimental conditions, e.g. following ischemia/reperfusion, or topical application of a whole range of inflammatory mediators, parasites and bacteria. With this technique, vascular permeability changes can be studied in real time and be related to other microvascular events such as leukocyte adhesion and activation. The cheek pouches are examined by intravital fluorescence microscopy using suitable filters (490/520 nm) and images are captured with a digital camera. A 5 % FITC-dextran 150 solution (approx. 100 mg/kg body-weight in normal saline) is administered i.v. (8-10).

Permeability studies using combined fluorescence stereomicroscopy, fluorescence light microscopy was reported by Thorball (2). This paper also includes tissue fixation techniques in the presence of

FITC-dextran and details of the microscopy set-up (filters, illumination). Regenerative titanium ear chambers (rabbits) have been used to study the blood/lymph systems in the microcirculation with FITC-dextran. Lymph ingrowth is seen after 4-8 weeks of implantation (11).




Fig.5. Images taken from cheek pouch after infusion of FITC-dextran 150. The second image shows the leakage of the microvasculature after subjection to histamine. (By kind permission of E. Svensjö).

2. Permeability studies on intestinal tissues

Permeability of intestinal epithelial monolayer during inflammation was studied using a FITC-dextran 150 (12). The action of a protease inhibitor on mucosal erosions and epithelial dysfunction in the GI tract was studied using FITC-dextran 4 (13). Thorball made extensive studies on tissue fixation in connection with FITC-dextran in the GI tract (2). FITC-dextran (4000-70000) were used to study permeability changes following cutaneous thermal injury in vitro using modified Ussing chambers (14). See also useful references (15,16).

3. Permeability studies of brain and nervous system.

For studies of FITC-dextran in the nervous system, techniques are required that immobilize the tracers enabling good optical resolution of features such as neuron, glial cells, sciatic endoneurium, ganglion etc. A method employing freeze-drying of the dissected tissue samples, fixation with formaldehyde vapour at 80 °C and embedding in paraffin wax or plastic proved to give good results (17). The sections from these preparations were dipped in xylene before mounting and examination.

The distribution of intravenously injected FITC-dextran (Mw from 3000 to 150 000) was studied in hamsters and mice. The normal perineural diffusion barrier of the cerebral cortex was not permeable to any FITC-dextran in the above range (19). The hamsters received 5mg FITC-dextran/10g bodyweight in 0.5ml saline (18). Permeability characteristics of the hippocampus were investigated with FITC-dextran. The blood-brain barrier was subjected to focused ultrasound and the passage of FITC-dextran was analysed (19).

The time-dependent permeability of the cerebral vessels following ischemia-reperfusion injury was examined using iv. FITC-dextran 150 (20). Brain edema following acute liver failure has been studied by monitoring protein related-permeability changing using FITC-dextran (21).

3. Permeability studies on neoplastic tissues.

Gerlowski (22) describes the use of transmitted light fluorescence TV microscopy for real time studies of microvascular permeability in neoplastic tissues. The studies were performed in the rabbit ear chamber model using FITC-dextran 150. The authors also examined the quenching of FITC-dextran concentrations from 0.6 mg to 30 mg/100 ml at intervals up to 90 min. Only the very highest

concentrations showed indications of quenching effects. T. Li and co-workers (23) studied membrane permeability of Sarcoma cells following focused ultrasound using FITC-dextran 500. FITC-dextran 500 was used to study the endocytosis of cells in experiments with ovarian cancer cells (24).

4. Permeability studies within the ocular chamber

E. Mannerma and co-workers (25) studied drug permeation process in the retinal pigment epithelium using FITC-dextran 40 and other probes. FITC-dextran 40 was used to assess barrier dysfunction in response to TNF-alpha (26). S. Lightman and co-workers (27) were able to study in detail the leakage of the retinal vessels in an inflammation model. Repeated doses of FITC-dextran were well tolerated. Studies on the pathways by which FITC-dextran and fluorescein leave the vitreous body have been described (28). Uveoscleral outflow in normal inflamed eyes has also been examined by C.B.Toris and co-workers (29) using FITC-dextran. An in vitro model of the outer blood retinal barrier has been described (30).

5. Permeability studies of renal tissues

FITC-dextran 4 is very rapidly excreted into the urine with 70 % of the dose excreted within the first hour (rabbits, rats) (unpublished data). Lencer and coworkers (31) reported studies with FITC-dextran 10 as a probe for endosome function and localization in the kidney. Frozen sections of excised tissue were examined by epifluorescent microscopy.

Diverse applications

FITC-dextran has been used to explore the permeability of nasal mucosa in vitro (32). In this study charged FITC-dextran derivatives were also employed. Diffusion of drugs vertically into the skin was studied using FITC-dextran 20 loaded dissolving micro-needles (33).

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada 
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland 
Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland 
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

©AMSBIO / AMS Biotechnology (Europe) Ltd. AMSBIO is the global source of all the products listed here.
All trademarks mentioned here are the property of their respective owners.

References

1. P.Kurtzhals, C.Larsen, and M.Johansen. High performance size-exclusion chromatographic procedure for the determination of fluoresceinyl isothiocyanate dextrans of various molecular masses in biological media. *J Chromatogr.* 1989;491:117-127.
2. N.Thorball. FITC-dextran tracers in microcirculatory and permeability studies using combined fluorescence Stereo Microscopy, Fluorescence Stereo microscopy and electron microscopy. *Histochemistry.* 1981;71:209-233.
3. K.G.Ljungström, H.Rench, K.Strandberg et al. Adverse reactions to dextran in Sweden 1970-1979 *Acta Chir Scan.* 1983; 149:253-262
4. H.Hedin, W.Richter and J.Ring. Dextran-induced anaphylactoid reactions in man: role of dextran reactive antibodies *Int. Arch Allergy Appl Immun.* 1976; 52:145-159.
5. A.N de Belder and K.Granath. Preparation and properties of fluorescein labelled dextrans. *Carbohydr Res.* 1973;30:375-378.
6. M.J.Geisow, P.D'Arcy Hart and M.R.Young. Temporal changes of lysosome and phagosome pH during phagolysosome formation in macrophages; studies by fluorescence spectroscopy. *J Cell Biol.* 1981;89: 645-652.
7. S.Ohkuma and B.Poole. Fluorescence probe measurement of the intralysosomal pH in living cells and the perturbation of pH of by various agents. *Proc Natl Acad Sci USA.* 1978;75:3327-3331.
8. E.Svensjö, E.M Saraiva, M.T.Bozza, et al. Salivary gland homogenates of *Lutzomyia longipalpis* and its vasodilatory peptide maxadilan cause plasma leakage via PAC 1 receptor activation. *J. Vasc Res.* 2009;46:435-446.
9. E.Svensjö. The hamster cheek pouch as a model in microcirculation research. *Eur Respir J Suppl* 12. 1990; 595s-600s; discussion 600s-601s.
10. E.Svensjö. The hamster cheek pouch as a research model in inflammation. Chapter 30 in David Shepro (Editor), *Microvascular Research-Biology and Pathology*, p.195-207, 2006, Elsevier Academic Press.
11. J.Jonsson, K.E.Arfor and H.C.Hint. Studies on relationships between blood and lymphatic systems within the microcirculation. 6th Europ Conf Microcirculation, Aalborg. 170;214-218 (Karger; Basel 1971).
12. T.Mochizuki, H.Satsu, M.Totsuka and M.Shimizu. Transepithelial transport of macromolecular substances in IL-4 treated human intestinal T84 cell monolayers. *Biosci Biotechnol Biochem.* 2009;73:2422-6.
13. B.Lei, W.Zha, Y.Wang, et al. Development of a novel self-microemulsifying drug delivery system for reducing HIV protease inhibitor-induced intestinal epithelial barrier dysfunction. *Mol Pharm.* 2010; 7:844-53.
14. F.Bertiaume, R.M.Ezzell, M. Toner et al. Transport of fluorescent dextrans across the rat ileum after cutaneous thermal injury. *Critical Care Medicine.* 1994;22:455-64.
15. N.Pantzar, B.R.Weström, A.Luts and S.Lundin. Regional small-intestinal permeability in vitro to different-sized dextrans and proteins in the rat. *Scand J Gastroenterol.* 1993;28:205-11.
16. C.Tagesson, R.Sjödahl, B.Thorén, Passage of molecules through the wall of the gastrointestinal tract.. A simple experimental model. *Scand J Gastroenterol.* 1978;13:519-24.
17. D.Hultström. FITC-dextrans in neurobiological research. *Acta Universitatis Upsaliensis*, 438, Almquist and Wiksell, Uppsala, 1982 and references cited therein.
18. Y.Olsson, E.Svensjö, K.E.Arfor. Fluorescein labelled dextrans as tracers for vascular permeability studies in the nervous system. *Acta Neuropathol.* 1975; 33:45-50.

19. W.Shougang, B.Baseri, J.Choi et al. Delivery of fluorescent dextrans through the ultra-sound induced blood-brain barrier opening in mice. Ultrasonics Symposium, 2008; IUS2008:1702-1705.
20. T.Ishi, T.Asai, T.Urakami et al. Accumulation of macromolecules in brain parenchyma in acute phase of cerebral infarction/reperfusion. Brain Res. 2010;19:164-8.
21. F.Chen, N.Ohashi, W.Li et al. Disruptions of occludin and claudin-5 in brain endothelial cells in vitro and in brains of mice with acute liver failure. Hepatology. 2009;50:1914-23.
22. L.E.Gerlowski and R.K.Jain. Microvascular permeability of normal and neoplastic tissues. Microvasc. Res. 1986; 31:288-305.
23. T.Li, Q.Hao, and X.Wang et al. The effect of focused ultrasound on the physiochemical properties of Sarcoma 180 cell membrane. Sheng Wu Yi Xue Gong Cheng Xue Za Zhi, 2009 Oct;26(5);941-6. Chinese
24. F.Chen, M.Hou, F.Ye et al. Ovarian cancer cells induce peripheral mature dendritic cells to differentiate into macrophage-like cells in vitro. Int J Gynecol Cancer. 2009;19:1487-93.
25. E.Mannermaa, M.Reinisaalo, V.P.Ranta et al. Filter cultured ARPE-10 cells as outer blood-retinal barrier model. Eur J Pharm Sci. 2010;40:289-296.
26. M.Shivanna, G.Rajashekhar and S.P.Srinivas. Barrier dysfunction of the corneal endothelium in response to TNF alpha. Invest Ophthalmol Vis Sci. 2010;51:1575-82.
27. S.L.Lightman, L.E.Caspers-Velu, S.Hirose et al. Angiography with fluorescein-labelled dextran in a primate model of uveitis. Arch Ophthalmol. 1987;105:844-848.
28. M.Arare and D.M.Maurice. The loss of fluorescein, fluorescein glucuronide and fluorescein isothiocyanate dextran from the vitreous by the anterior and retinal pathways. Exp Eye Res. 1991;52;27-39.
29. C.B.Toris, D.S.Gregerson and J.E.Pederson. Uveoscleral outflow using different-sized fluorescent tracers in normal and inflamed eyes. Exp Eye Res. 1987;45;525-32.
30. R.D.Hamilton, A.J.Foss and L.Leach. Establishment of a human in vitro model of the outer blood-retinal brain barrier. J. Anat. 2007;211:707-16.
31. W. I. Lencer, P.Weyer, A.S. Verkman et al. FITC-dextran as a probe for endosome function and localization in kidney. Amer J Physiol. 1990;258;C309-317.
32. N.Uchida, Y.Maitani Y.Machida et al. Influence of bile salts on the permeability of insulin through the nasal mucosa of rabbits in comparison with dextran derivatives. Int J Pharm. 1991;74;95-104.
33. K.Fukushima, A.Isa, H.Morita et al. Two-layered dissolving microneedles for percutaneous delivery of peptide/protein drugs in rats. Pharm Res. 2010 (E-pub ahead of print).

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada 
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland 
Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland 
T: +49 (0) 69 779099
F: +49 (0) 69 13376880