### <u>amsbio</u>

# Multiple-Choice™ First Strand cDNA

#### **OVERVIEW**

Amplification and cloning of cDNA by PCR is an essential tool for today's molecular biologist. This technique can be used to clone members of a multi-gene family, either homologues across species or previously identified. PCR amplification of cDNA and its subsequent subcloning can often save time and reduce costs by eliminating the need for the construction and screening of a cDNA library.

## PRODUCT FORMAT AND QUALITY CONTROL


Multiple-Choice™ First-Strand cDNA Kits are formatted for the rapid and efficient screening of a panel of tissues for gene expression. The cDNA panels provided in these kits are also ideal for the characterization of alternatively spliced mRNA either within a single tissue or across multiple organism tissues.

Each step in the production of Multiple-Choice First Strand cDNA is carefully monitored to ensure product quality. First strand cDNA is synthesized from poly A+RNA using an oligo(dT) primer and a reverse transcriptase that favors the production of long product. Each panel of cDNAs is examined by PCR using standard conditions for the detection of "house-keeping" cDNAs, such as b-actin and cyclophilin. Small aliquots of these primers are included as controls in each kit.

#### PRODUCT APPLICATION

Multiple-Choice First Strand cDNA is PCR\* ready and available as individual cDNA (30 reactions) or in sets of six cDNAs (10 reaction each).

When Multiple-Choice cDNA is used as a template for PCR, it results in specific amplification of the desired sequence. As shown the figure below, an actin signal is detectable in each cDNA sample of the rat Multiple-Choice cDNA Set 1, whereas the liver-specific cytochrome P450 can only be amplified from the liver cDNA sample.



#### **KIT COMPONENTS**

- First strand cDNAs derived from either human, rat or mouse tissues (2-10ng/ul)
- Control primers for B-Actin
- Control primers for Cyclophilin (Multiple Choice™ cDNA Set only)







## Multiple-Choice™ First Strand cDNA

#### **ORDERING INFORMATION & PRICING**

MULTIPLE-CHOICE FIRST STANDARD cDNA	TISSUES AVAILABLE	CATALOG NUMBER
Human cDNA from individual tissues (30 reactions)	Brain, Heart, Kidney, Liver, Placenta, PBL, Lung, Muscle, Ovary, Prostate, Small Intestine, Spleen, Stomach, Testis, Fetal Brain, Fetal Kidney, Fetal Liver, Fetal Muscle, Fetal Spleen, Fetal Testis	CH-10XX
Mouse cDNA from individual tissues (30 reactions)	Brain, Heart, Kidney, Liver, Lung, Muscle, Skin, Small Intestine, Spleen, Stomach, Testis, Thymus	CM-10XX
Rat cDNA from individual tissues (30 reactions)	Brain, Heart, Kidney, Liver, Lung, Muscle, Skin, Small Intestine, Spleen, Stomach, Testis, Thymus	CR-10XX
Human cDNA Set 1	Brain, Heart, Kidney, Spleen, Liver, PBL	CH-1101
Human cDNA Set 2	Lung, Muscle, Ovary, Prostate, Small Intestine, Testis	CH-1102
Human cDNA Set 3	Fetal Brain, Fetal Kidney, Fetal Liver, Fetal Muscle, Fetal Spleen, Fetal Testis	CH-1103
Mouse cDNA Set 1	Brain, Heart, Kidney, Liver, Spleen, Thymus	CM-1101
Mouse cDNA Set 2	Lung, Muscle, Skin, Small Intestine, Stomach, Testis	CM-1002
Rat cDNA Set 1	Brain, Heart, Kidney, Liver, Spleen, Thymus	CR-1001
Rat cDNA Set 2	Lung, Muscle, Skin, Small Intestine, Stomach, Testis	CR-1002

#### **PRODUCT DELIVERY AND STORAGE**

Multiple-Choice First Stand cDNA is shipped on dry ice. Product storage at −20°C is recommended.

\*PCR is covered by patents issued to Cetus Corporation and is owned and licensed by Hoffman-LaRoche Molecular Systems, Inc. Purchase of any of OriGene's PCR-related products does not convey a license to use the PCR process covered by these patents. Purchasers of these products must obtain a license to use the PCR process before performing PCR.





Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22

Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19 60325 Frankfurt/Main Tel: +49 (0) 69 779099 Fax: +49 (0) 69 13376880 United States

23591 El Toro Rd, Suite #167 Lake Forest, CA 92630 Tel: + 1 855 267 2464 Fax: + 1 949 768 8613