

Human Adipogenesis Assay Kit

PPAR γ Agonists

Cat# DIF-AG; DIF-AG-NC

INSTRUCTION MANUAL ZBM0005.01

STORAGE CONDITIONS

Frozen subcutaneous preadipocytes

Store in liquid nitrogen IMMEDIATELY upon receipt. No expiration date is applicable; however, the cells must be plated within 1 week of receiving the kit to account for the expiration of the kit components.

Media, Reagents A & B, Buffers:

Store at 2 - 8°C.

Glycerol Standard

-20°C

For Research Use Only

Not For Use In Diagnostic Procedures

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, expressed or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Zen-Bio, Inc. Zen-Bio, Inc. shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

The product accompanying this document is intended for research use only and is not intended for diagnostic purposes or for use in humans.



UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

United States

23591 El Toro Rd, Suite #167
Lake Forest, CA 92630
Tel: + 1 800 987 0985
Fax: + 1 949 265 7703

INTRODUCTION

The differentiation assay kits provide the tools to study the compounds that stimulate cultured human adipocyte differentiation or lipogenesis. Such compounds may be PPAR γ agonists or a combination of thiazolidinediones and glucocorticoids that are potentially useful in the treatment of diabetes.

This kit is designed to test compounds as potential PPAR γ agonists. It is our experience that the PPAR γ agonist used as the Positive Control sufficiently stimulates human adipocyte differentiation after 7 days of treatment. The protocol is designed so that the test compounds are also used in a 7-day treatment regimen. This kit contains sufficient reagents to assay 100 assay points in a 96 well format.

ITEMS INCLUDED IN THE KIT

Item	Description	Unit	QTY	Item	Storage
Human SQ Preadipocytes	Human subcutaneous preadipocytes, $\geq 2.0 \times 10^6$ cells/vial, cryopreserved	VIAL	1	1	Liquid nitrogen
PM-1	Preadipocyte medium (See Appendix A)	BOTTLE	50ml	2	4°C
NC	Negative control: (See Appendix A)	VIAL	500 μ l	2	4°C
PC	Positive control: (See Appendix A)	VIAL	500 μ l	2	4°C
VC	Vehicle control (See Appendix A)	VIAL	500 μ l	2	4°C
IM	Initiation medium (See Appendix A)	BOTTLE	300ml	3	4°C
MM	Maintenance medium (See Appendix A)	BOTTLE	100ml	2	4°C
Wash buffer		BOTTLE	50ml	2	4°C
Lysis buffer		BOTTLE	25ml	2	4°C
Glycerol Reagent A	Reconstitute w/ 11 ml deionized water prior to use	BOTTLE	11ml	2	4°C
Reagent B	Reconstitute w/ 2.5 ml deionized water prior to use	BOTTLE	2.5 ml	2	4°C
Glycerol standard	Glycerol @ 1mM [Reconstitute with 200 μ l Standards Diluent to make the 200 μ M glycerol standard; see page 5 for recommended dilution scheme]	VIAL	50 μ l	2	-20°C
Standards diluent		BOTTLE	2 ml		4°C
Tray	Clear polyvinyl tray for multi-channel pipetters	EACH	3	-----	-----
Data sheet	Certificate of Analysis and protocol	EACH	1	-----	-----
Plate A	96-well plate for plating and differentiating	PLATE	1	---	
Assay Plates, blank	96-well assay plate, blank	PLATE	2	---	-----

Other equipment/reagents required but not provided with the kit:

- ▶ Additional Maintenance Medium, if necessary (see background information)
- ▶ Single-channel pipetter
- ▶ Multi-channel pipetter
- ▶ Plate reader with a filter of 540 nm
- ▶ Tubes to dilute glycerol standards

ASSAY PROCEDURE

A. DIFFERENTIATION PROCEDURE

On each day of the procedure, the appropriate medium must be warmed to 37° C prior to use.

Note: This protocol is designed to accommodate a weekday work schedule if started on a Monday-Thursday. Any deviation may require weekend work. We strongly recommend testing all compounds in triplicate.

Day 1: This is the day the cells are plated.

1. Remove cells from liquid nitrogen and place immediately into a 37° C water bath and agitate while in bath. Be careful not to submerge the cap of the vial into water. Do not leave the vials in water bath after most of the content has thawed. Rinse the vials with 70% ethanol before taking them to the culture hood.
2. Upon thawing, transfer the cells to a sterile conical bottom centrifuge tube containing 10 ml of Preadipocyte medium (cat # PM-1).
3. Centrifuge: 1,200 rpm (282Xg) / 20°C / 5 minutes. Aspirate the supernatant. **TAKE CARE TO NOT ASPIRATE ANY OF THE CELL PELLETT.**
4. The cell vial contains a minimum of 2.0×10^6 viable cells; however, we recommend performing a cell count to determine a more exact number of cells. Resuspend the cell pellet in 2 ml Preadipocyte medium; dilute an aliquot in 0.4% trypan blue solution. We suggest withdrawing an aliquot of 50 μ l of cells and mixing with 100 μ l of the trypan blue solution, resulting in a dilution factor of 3. Count live (unstained) cells on a hemacytometer. The cell concentration required for approximately 40,000 cells / cm^2 in the 96 well format with 150 μ l /well is 1.3×10^6 cells in 15 ml Preadipocyte medium.
5. Plate cells in one of the 96 well plates provided in the kit. Do not agitate the plate, as cells will not plate evenly.
6. Place plate in 37°C incubator, 5% CO₂, 97% humidity. The cells will be maintained in the incubator after each manipulation until Day 14.

NOTE: This kit contains a sufficient volume of Initiation medium (IM) to use 10 ml of medium per compound dilution for a maximum of 29 compounds tested in triplicate (87 wells) on the 96-well plate, leaving 9 wells for controls. If a compound stock is too concentrated to accomplish the desired dilution, use an appropriate solution (not supplied) to prepare an intermediate concentration that would allow for a final volume of 10 ml.

Also the positive control in this kit, the PPAR γ agonist, has a final solvent concentration of 0.1% DMSO. This is included in the vehicle control (VC). If the concentration of any solvent for the compounds used is high enough to potentially alter differentiation, please include that solvent

concentration as an additional treatment. We do not recommend treating the cells with solutions exceeding 1% of any solvent, as higher concentrations may be toxic to the cells.

Day 2:

7. Twenty-four hours later, check the cells for confluence.
8. Using the Initiation Medium (IM), prepare treatments. Plan to do all treatments in triplicate. A blank plate map is included in these instructions to record the well treatments.
9. When all treatments are prepared, remove Preadipocyte medium from control wells. We recommend doing the treatments in small groups so the cells do not dry out.
10. Pipet 150 μ l Positive Control, (PC), 150 μ l Negative Control (NC), and 150 μ l Vehicle Control (VC) into appropriate wells.
11. Remove media from experimental wells and pipet 150 μ l each treatment media into appropriate wells. Incubate the plate for 7 days.

Day 8:

12. Using a multichannel pipetter remove 100 μ l media from all wells. Gently feed all wells with 100 μ l of the Maintenance Medium (MM) that is provided with this kit. Incubate the plate for 7 days.

Day 15:


13. Cells are now mature. Proceed to part B.
14. The positive control wells should exhibit significantly greater lipid accumulation than the negative control wells or the vehicle control wells. Refer to page 7 for a picture of a typical positive control when the adipocytes are mature.

B. TRIGLYCERIDE ASSAY

1. Warm the Wash buffer and Lysis buffer in a 37°C water bath.
2. Prepare the Reagent B by adding 2.5ml deionized water per bottle and gently invert. DO NOT VORTEX! Use a pipet to ensure that the powder is completely dissolved. Keep at room temperature. Store in a light protected bottle. Reconstituted Reagent B is stable for 60 days refrigerated (2-8°C); store any remaining solution refrigerated (2-8°C). Bring Reagent B to room temperature.
3. Remove all media. Using about 15 ml of the wash buffer, wash the cells one time with 150 μ l wash buffer. Label the disposable tray "wash buffer" and retain for later use.
4. Remove all Wash buffer. Using a new tray, add 15 μ l Lysis buffer. Incubate at 37°C – 50°C for 20 minutes.
5. After the incubation is complete, visually confirm cell lysis by checking the wells under a microscope. If cells are not fully lysed, incubate another 10 minutes.
6. Add 135 μ l warm Wash Buffer and mix the lysates by pipetting up and down three times.
7. Add 20 μ l Reagent B to each well. It is not necessary to mix at this time, however, gently tap the plate to help mix the reagents. Incubate the plate at 37°C for 2 hours.

8. Bring Reagent A and the glycerol standards to room temperature during this time. The Wash Buffer can also be kept at room temperature at this point. Warm the Standards Diluent to 37°C. Prepare the standard curve as follows:


Pipette 200 μl of the Standards Diluent into the 1 mM glycerol standard tube provided and mix well by vortexing. This produces a diluted stock glycerol standard of 200 μM . Pipette 125 μl of diluent into 6 tubes (not provided). Using the newly diluted stock glycerol solution, prepare a dilution series as depicted below. Mix each new dilution thoroughly before proceeding to the next. The 200 μM stock dilution serves as the highest standard, and the diluent serves as the zero standard.



9. Also at this time prepare the Reagent A by adding 11ml deionized water per bottle and gently invert. DO NOT VORTEX! Use a pipette to ensure that the powder is completely dissolved. Keep at room temperature. If using a Reagent A solution previously prepared and stored at 2-8°C, also bring to room temperature. Make sure there is enough Reagent A from one solution to treat all the points in the assay. It may be necessary to combine solutions. Store in a light protected bottle. Reconstituted Glycerol Reagent A is stable for 60 days refrigerated (2-8°C); store any remaining solution refrigerated (2-8°C).
10. To a blank 96 well plate, add 80 μl wash buffer to each well needed for the assay (NOTE: do not add Wash Buffer to the wells used for the standard curve).
11. Working with one row or column at a time, mix the lysates very well using a multi-channel pipette. Immediately transfer 20 μl per well of the lysates to the corresponding well of the plate containing the wash buffer. **This results in a Dilution Factor of 5.**
12. Prepare the standard curve. Pipet 100 μl of each standard into a well. (NOTE: Eight wells are necessary for the curve. If there are remaining wells on the assay plate, you can utilize the remaining wells. If not, a second plate is included in this kit).
13. Using the third tray, add 100 μl Reagent A to samples and standards. Mix by pipetting up and down one time. Incubate at room temperature for 15 minutes.
14. Read at 540 nm using a microtiter plate reader.

GLYCEROL STANDARD CURVE

This kit is designed to show relative lipid accumulation of experimental treatments compared to controls. The assay is based on the equation



The reagent measures the concentration of glycerol released after lysing the cells and hydrolyzing the triglyceride molecules. The triglyceride concentration can then be determined from the glycerol values.

Generate standard curve: see example below


[DO NOT use this standard curve to generate your data. This is an example.]

Subtract the OD value of the 0 μ M standard from all OD values including the standard curve.

Zero (blank) = .040

μ M Glycerol	OD	OD - blank
3.125	0.054	0.014
6.25	0.066	0.026
12.5	0.082	0.042
25	0.138	0.098
50	0.214	0.174
100	0.402	0.362
200	0.711	0.671

slope =	0.0034
intercept=	0.0075
r^2 =	0.9985



y = observed O.D. minus the blank

x = concentration of glycerol in μ M

To calculate x for each y, (i.e. to change the observed O.D. into glycerol concentration) use the following equation: $y = (\text{slope}) \times (x) + \text{intercept}$ $y = mx + b$ so $x = (y - b) / m$

$x = (y - 0.0075) / 0.003$ where 0.003 = slope of the line and 0.0075 = y intercept. Be careful to enter the proper sign for the y intercept value as it may be a negative number.

The R^2 value should be equal or greater than 0.98 for the standard curve to be valid. Any R^2 values below 0.98, must have the standard curve run again.

Solve for the Total Glycerol concentration (i.e. total triglyceride concentration) for each OD. **Remember to include the Dilution Factor in the equation.** Data is expressed as μ M Glycerol.






NOTE: Any OD values that are negative after the blank is subtracted should be considered to be 0 for the OD value. Also any OD values greater than the highest standard (200 μ M), should be suspect. The compound should be re-assayed using a lower dose of the compound at treatment OR a more dilute solution of the condition medium at the time of the assay.

TROUBLESHOOTING

Problem	Suggestions
High background or the triglyceride reagent turns a darker color before the assay begins.	<ul style="list-style-type: none"> • Use clean tray and tips • Change pipet tips frequently
Edge effects	<ul style="list-style-type: none"> • Ensure a saturated humidity in the incubator to prevent evaporation from the outside wells
Inconsistent OD reading	<ul style="list-style-type: none"> • Be careful when pipetting to avoid bubbles. If bubbles persist, burst the bubbles using a large gauge needle prior to reading and read the plate again. • Mix the lysates well before transferring the 20μl to the Wash buffer plate.
Cells appear dead after 7 days treatment with my compound	<ul style="list-style-type: none"> • An acute treatment for 3 days in Initiation Medium followed by additional feedings each 7 days should yield suitable positive control signal to complete the assay.

APPENDIX A. DIFFERENTIATION PICTURES

PREADIPOCYTE \longrightarrow MATURE ADIPOCYTE



APPENDIX B: COMPOSITIONS OF KIT REAGENTS



Reagents	Vehicle Control (VC) IBMX, 0.1% DMSO	Negative Control (NC) IBMX, PPAR γ agonist, TNF α	Positive Control (PC) IBMX, PPAR γ agonist
Components	<ul style="list-style-type: none"> • DMEM / Ham's F-12 medium • HEPES • Fetal bovine serum • Biotin • Pantothenate • Human insulin • Dexamethasone • Penicillin • Streptomycin • Amphotericin B • Isobutylmethylxanthine (IBMX) • DMSO 	<ul style="list-style-type: none"> • DMEM / Ham's F-12 medium • HEPES • Fetal bovine serum • Biotin • Pantothenate • Human insulin • Dexamethasone • Penicillin • Streptomycin • Amphotericin • Isobutylmethylxanthine (IBMX) • PPARγ agonist • TNFα 	<ul style="list-style-type: none"> • DMEM / Ham's F-12 medium • HEPES • Fetal bovine serum • Biotin • Pantothenate • Human insulin • Dexamethasone • Penicillin • Streptomycin • Amphotericin B • Isobutylmethylxanthine (IBMX) • PPARγ agonist
Reagents	Preadipocyte Medium		
Components	<ul style="list-style-type: none"> • DMEM / Ham's F-12 medium • HEPES • Fetal bovine serum • Penicillin • Streptomycin • Amphotericin B 	Maintenance Medium (MM) <ul style="list-style-type: none"> • DMEM / Ham's F-12 medium • HEPES • Fetal bovine serum • Biotin • Pantothenate • Human insulin • Dexamethasone • Penicillin • Streptomycin • Amphotericin B 	Initiation Medium (IM) <ul style="list-style-type: none"> • DMEM / Ham's F-12 medium • HEPES • Fetal bovine serum • Biotin • Pantothenate • Human insulin • Dexamethasone • Penicillin • Streptomycin • Amphotericin B • Isobutylmethylxanthine (IBMX)

APPENDIX C: PLATE LAYOUT

H	G	F	E	D	C	B	A	
								1
								2
								3
								4
								5
								6
								7
								8
								9
								10
								11
								12

APPENDIX D: DIFFERENTIATION FLOWCHART _____

DAY 1
PLATE CELLS. INCUBATE 24 HOURS 37°C



DAY 2
APPLY TREATMENTS IN INITIATION MEDIUM
INCUBATE 7 DAYS 37°C




DAY 8
CHANGE CELLS TO MAINTENANCE MEDIUM
INCUBATE 7 DAYS 37°C



DAY 15
CELLS ARE MATURE. MOVE ON TO
TRIGLYCERIDE ASSAY PROTOCOL

APPENDIX E: TRIGLYCERIDE ASSAY FLOWCHART



UK & Rest of World
184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

United States
23591 El Toro Rd, Suite #167
Lake Forest, CA 92630
Tel: +1 800 987 0985
Fax: +1 949 265 7703

amsbio
info@amsbio.com