

Avidin/Biotin Blocking Kit

Avidin Solution and Biotin Solution for Endogenous Avidin or Biotin

Storage: 2-8°C

Catalog No: E08-18

18 ml

Intended Use:

This kit is composed of 2 reagents: Avidin Solution and Biotin Solution. This product is designed to eliminate the non-specific binding caused by the endogenous avidin or biotin like substances in immunohistochemistry and cytochemistry, when either avidin or Streptavidin/biotin is used as the signal amplification system.

Kit Components:

- | | |
|---|------|
| 1) Reagent 1: one dropper bottle of Avidin solution | 18ml |
| 2) Reagent 2: one dropper bottle of Biotin solution | 18ml |

Recommended Protocol:

These blocking steps should be finished prior the primary antibody incubation step.

After incubation with normal goat serum:

- 1) Apply 2 drops (50 µL/ drop) of reagent 1 (Avidin Solution) to cover the tissue and incubate for 10 minutes at room temperature.
- 2) Rinse well with phosphate buffered saline (PBS).
- 3) Apply 2 drops (50 µL/ drop) of reagent 2 (Biotin Solution) to cover the tissue and incubate for 10 minutes at room temperature
- 4) Rinse well with phosphate buffered saline (PBS).

The endogenous avidin/ biotin activity is now blocked and the tissue is ready for standard immunostaining procedure

Storage:

Store at 2-8°C. Return reagent back to 2-8°C after each use.

Remarks:

For research use only.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85