

Luciferase (Firefly and Renilla) Expression Stable Cell Lines

Catalog #	Product	Amount
HEK293 Host cell line		1 vial of cells (2 x 10 ⁶ cells) in 90% complete medium, 10% DMSO
<u>SC002-Bsd</u>	Luciferase (firefly), HEK293 stable cells with blasticidin antibiotic marker	
<u>SC002-Puro</u>	Luciferase (firefly), HEK293 stable cells with Puromycin antibiotic marker	
<u>SC002-Neo</u>	Luciferase (firefly), HEK293 stable cells with Neomycin antibiotic marker	
<u>SC002-GB</u>	Luciferase (firefly), HEK293 stable cells with GFP fluorescent and Blasticidin antibiotic dual marker	
<u>SC002-GP</u>	Luciferase (firefly), HEK293 stable cells with GFP fluorescent and Puromycin antibiotic dual marker	
<u>SC002-RB</u>	Luciferase (firefly), HEK293 stable cells with RFP fluorescent and Blasticidin antibiotic dual marker	
<u>SC002-RP</u>	Luciferase (firefly), HEK293 stable cells with RFP fluorescent and Puromycin antibiotic dual marker	
<u>SC020-Puro</u>	Luciferase (Renilla), HEK293 stable cells with Puromycin antibiotic marker	
<u>SC020-RP</u>	Luciferase (Renilla), HEK293 stable cells with RFP fluorescent and Puromycin antibiotic dual marker	
<u>SC021-Puro</u>	Luciferase (firefly) and CRE recombinase, co-expression stable cell line in HEH293 cells with Puromycin antibiotic marker	
<u>SC021-RP</u>	Luciferase (firefly) and CRE recombinase, co-expression stable cell line in HEK293 cells with RFP fluorescent and Puromycin antibiotic dual marker	
<u>SC021-GB</u>	Luciferase (firefly) and CRE recombinase,	

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85

	co-expression stable cell line in HEK293 cells with GFP fluorescent and Blastcidin antibiotic dual marker
HeLa Host cell line	
<u>SC032-Bsd</u>	Luciferase (firefly), HeLa stable cells with blastcidin antibiotic marker
<u>SC032-Puro</u>	Luciferase (firefly), HeLa stable cells with Puromycin antibiotic marker
<u>SC032-GB</u>	Luciferase (firefly) and GFP co-expression stable cell line in HeLa cells with Blastcidin antibiotic dual marker
<u>SC032-GP</u>	Luciferase (firefly) and GFP co-expression stable cell line in HeLa cells with Puromycin antibiotic dual marker
<u>SC032-GN</u>	Luciferase (firefly) and GFP co-expression stable cell line in HeLa cells with Neomycin antibiotic dual marker
<u>SC032-RN</u>	Luciferase (firefly) and RFP co-expression stable cell line in HeLa cells with Neomycin antibiotic dual marker
<u>SC032-RB</u>	Luciferase (firefly) and RFP co-expression stable cell line in HeLa cells with Blastcidin antibiotic dual marker
<u>SC032-RP</u>	Luciferase (firefly) and RFP co-expression stable cell line in HeLa cells with Puromycin antibiotic dual marker
MDA-MB-231 Host cell line	
<u>SC041</u>	MDA-MB-231 / Luciferase-2A- RFP Stable Cell Line (Blastcidin)
<u>SC044</u>	MDA-MB-231 / Luciferase-2A- GFP Stable Cell Line (Puromycin)
A549 Host cell line	
<u>SC043-Luc</u>	A549 / Luciferase stable cell line (Puromycin)





<u>SC043-LG</u>	A549 / Luciferase-2A- GFP stable cell line (Puromycin)	
Jurkat T Cells		
<u>SC048</u>	Jurkat T Cell /Firefly Luciferase (Puromycin) Stable Cell Line	
<u>SC050-L</u>	MCF7 / Firefly Luciferase (Puromycin) Stable Cell Line	
<u>SC051-L</u>	ZR-75-1 / Firefly Luciferase (Puromycin) Stable Cell Line	

Storage: Liquid Nitrogen

Product Description

The **HEK293** Cell Line is a permanent line established from primary embryonic human kidney transformed with sheared human adenovirus type 5 DNA. The expressed E1A adenovirus gene allows these cells to produce very high levels of protein.

The **HeLa** cell line is derived from human cervical cancer cells, and widely used as the host cells for all kinds of biological tests, as it divides an unlimited number of times in a laboratory cell culture plate as long as fundamental cell survival conditions are met.

The **MDA-MB-231** cell line is one of the common human breast cancer cell line, derived from metastatic breast cancer, mammary gland epithelial cells.

The **A549** cell line a human lung epithelial cancer cell line, first developed through culture of lung carcinomatous tissue from a 58-year-old Caucasian male.

The **Jurkat** cells are an immortalized line of human T lymphocyte cells that are used to study acute T cell leukemia, T cell signaling, and the expression of various chemokine receptors susceptible to viral entry, particularly HIV. They are also widely used for the drug screening, signal transduction and immunotherapy model cells.

The **MCF7** is human breast cancer cell line, an epithelial cancer cell line derived from breast adenocarcinoma. It is widely used for *in vitro* breast cancer studies. MCF-7 cells are able to process estrogen in the form of estradiol via estrogen receptors in the cell cytoplasm. This results in the MCF-7 cell line being an estrogen receptor (ER) positive cell line. MCF-7 is also progesterone receptor positive and HER2 negative.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880




Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85



The **ZR-75-1** is an adhesive, human breast epithelial cancer cell line, derived from a malignant ascitic effusion in a 63 year old female Caucasian with infiltrating ductal carcinoma. The cells are reported to express both the wild type and variant estrogen receptors, progesterone receptor and other steroid hormones. It is widely used for *in vitro* breast cancer research.

Luciferase stable cell lines are transformed from the different host cell lines listed above. Each cell line stably expresses a luciferase gene (**firefly luciferase** or **Renilla luciferase**). The cell lines are established by transduction of luciferase expressing lentivirus with different antibiotic selection marker. The luciferase is constitutively expressed at high level under either a super strong CMV promoter (suCMV) or an enhanced EF1a promoter. When included, a fluorescent protein marker (GFP or RFP) is co-expressed under the same CMV promoter (Note: the fluorescent protein and luciferase are expressed as individual protein, not fusion, mediated by a 2A element). The antibiotic selection marker is expressed under RSV promoter. AMSBIO also provides the stable cell lines expressing both a luciferase and a CRE recombinase with different antibiotic markers. Please see schemes below for the **expression cassette structure** which was integrated into the stable cell line's genome.



Culture procedures

1. Thaw the frozen vial of cells quickly in a 37°C water bath (1-3min), decontaminate the outside of the vial with 70% ethanol.
2. Transfer the entire contents of the cryovial into a T-75 cm² flask containing 15 ml of pre-warmed complete medium. Incubate the cells overnight in a 37°C incubator, 5% CO₂.
3. (For adhesive cell type: on the following day, carefully replace the medium with 15 ml of pre-warmed, **complete medium**.)
4. Incubate the cells and monitor cell density.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85



5. Passage cells (1:10 dilution) when the culture reaches 80-90% confluency (or 1×10^6 cells/ml for suspension cells).
6. Freeze cells at a density of $2-5 \times 10^6$ cells/ml using 90% complete medium with 10% DMSO.

Complete medium

Depends on host cell type. Please refer to ATCC website for medium details of each host cell line.

Quality Control

Each vial contains greater than 2×10^6 cells with >95% viability before freeze. Cells are tested free of bacteria, viruses, mycoplasma.

Warranty and user terms

1. This product is warranted to perform as described when used in accordance with this manual. AMSBIO MAKES NO REPRESENTATIONS AND EXTENDS NO WARRANTIES OF ANY KIND, EITHER EXPRESS OR IMPLIED. AMSBIO's sole remedy for breach of warranty should be, at the option of AMSBIO, to repair or replace the product if this product does not meet the stated quality standard.
2. By paying the purchase price, the buyer is granted a non-transferable, non-exclusive license to use the product. This product is sold **for research and development purposes only**.
3. This product is limited to the laboratory that the product is delivered to. This Product is not for resale, distribution, or transfer for any purpose, including transfer of the Product as a component of any products
4. This Product should be used only for non-profit purposes, including any products and services usages. Furthermore, **research use only** means that this product is excluded, without limitation, from resale, repackaging, or modification for the purpose of making or selling of any commercial products or services without the written approval of AMSBIO.
5. AMSBIO is not liable, and does not have any responsibility or liability, whatsoever for any direct and indirect, consequential, or other damages resulting from using this Product.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85