

## MMLV-RT enzyme

**Storage Conditions:** Store at  $-20^{\circ}\text{C}$ . It is stable a year from the date of purchase if stored and handled properly.

### Quality Control:

1. Nonspecific Endonuclease/Nickase Activity: None detected after 14-16 hour incubation of 100 units of enzyme with 300 ng of supercoiled plasmid DNA, analysed on agarose gel.
2. Exonuclease Activity: None detected after 14-16 hour incubation of 100 units of enzyme with 40 ng of  $^{32}\text{P}$ -labeled Sau3A fragments of pUC19, analysed by PAGE.
3. Protease Activity: None detected after 14-16 hour incubation of 1200 units of enzyme with 1  $\mu\text{g}$  of protease substrate, analysed by fluorescence.

**Source:** *E.coli*

**Size and concentration:** 16,000 units or 100,000 units with 1 ml or 2.5 ml 5x reaction buffer, respectively at 200 units/ $\mu\text{l}$  concentration.

**Storage buffer:** 20 mM Tris-HCl (pH 7.5), 1 mM DTT, 0.05% (v/v) Triton X-100, 0.1 mM EDTA, 0.1 M NaCl and 50% (v/v) glycerol.

**Reaction buffer (5x):** 250 mM Tris-HCl (pH 8.3), 375 mM KCl, 15 mM MgCl<sub>2</sub>, and 50 mM DTT

**Unit definition:** One unit of the enzyme incorporates 1 nmole of dTTP into acid-precipitable material in 10 minutes at  $37^{\circ}\text{C}$  using poly (A): oligo (dT)<sub>25</sub> as template-primer.

## PROTOCOL

### First-Strand cDNA Synthesis

#### Materials to Be Supplied by the User:

Recombinant RNasin Ribonuclease Inhibitor (Cat.# Z5040001)

- dATP, 10mM (Cat.# K6011101-400 , 100mM)
- dCTP, 10mM (Cat.# K6011103-400 , 100mM)
- dGTP, 10mM (Cat.# K6011102-400 , 100mM)
- dTTP, 10mM (Cat.# K6011104-400 , 100mM)
- Nuclease-Free Water

1. The following procedure uses 2  $\mu\text{g}$  of RNA. In a sterile RNase-free microcentrifuge tube, add 0.5  $\mu\text{g}$  of the primer or primer-adaptor per microgram of the mRNA sample in a total volume of 15  $\mu\text{l}$  in water. Heat the tube to  $70^{\circ}\text{C}$  for 5 minutes to melt secondary structure within the

**Intended use of the products:** All products sold by AMSBIO are intended for research use only and must not be used for medical diagnostic or drug purposes. AMSBIO products may not resell, modified for resale or used to manufacture commercial products without written approval from AMSBIO.

**Limited product warranty:** AMSBIO has a total quality assurance policy. If an y product does not meet the specification stated in the product data sheet will be replaced at no charge. This warranty limits AMSBIO liability to the replacement of the product only. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by AMSBIO. AMSBIO shall not be liable for any incidental, direct, special or consequential damages, even if it is aware of the possibility of such damages.



#### UK & Rest of World

184 Milton Park, Abingdon  
OX14 4SE, Oxon, UK  
Tel: +44 (0) 1235 828 200  
Fax: +44 (0) 1235 820 482

#### Switzerland

Centro Nord-Sud 2E  
CH-6934 Bioggio-Lugano  
Tel: +41 (0) 91 604 55 22  
Fax: +41 (0) 91 605 17 85

#### Deutschland

Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
Tel: +49 (0) 69 779099  
Fax: +49 (0) 69 13376880

#### North America

23591 El Toro Rd, Suite #180  
Lake Forest, CA 92630  
Tel: + 1 800 987 0985  
Fax: + 1 949 265 7703

amsbio

info@amsbio.com

www.amsbio.com  
AMS Biotechnology

template. Cool the tube immediately on ice to prevent secondary structure from reforming and then spin briefly to collect the solution at the bottom of the tube.

2. Add the following components to the annealed primer/template in the order shown.

Note: Do not alter the ratio of primer to mRNA.

M-MLV 5X Reaction Buffer 5 µl; dATP, 10mM 1.25 µl; dCTP, 10mM 1.25 µl; dGTP, 10mM 1.25 µl; dTTP, 10mM 1.25 µl; RNasin.

Ribonuclease Inhibitor 25 units; M-MLV RT 200 units; Nuclease-Free Water to final volume 25 µl.

3. Mix gently by flicking the tube and incubate for 60 minutes at 37°C for random primers or 42°C for other primers or primer- adaptors. The extension temperature may be optimized between 37 and 42°C.
4. Perform second-strand synthesis using a protocol of your choice. Standard protocols for second-strand synthesis may be found in reference 2.

Note: The M-MLV RT Reaction Buffer is compatible with enzymes used in a number of downstream applications. Phenol extractions and ethanol precipitations typically are not necessary before performing second-strand synthesis and amplification.

#### IV. References

1. Roth, M.J., Tanese, N. and Goff, S.P. (1985) Purification and characterization of murine retroviral reverse transcriptase expressed in Escherichia coli. J. Biol. Chem. 260, 9326–35.

2. Sambrook, J., Fritsch, E.F. and Maniatis, T. (1989) In: Molecular Cloning: A Laboratory Manual, Cold Spring Harbor Laboratory, Cold Spring Harbor, New York, 8.64.

Note: M-MLV Reverse Transcriptase is less processive than AMV Reverse Transcriptase, and therefore, more units of the M-MLV enzyme are required to generate the same amount of cDNA as in the AMV reaction.

**Intended use of the products:** All products sold by AMSBIO are intended for research use only and must not be used for medical diagnostic or drug purposes. AMSBIO products may not resell, modified for resale or used to manufacture commercial products without written approval from AMSBIO.

**Limited product warranty:** AMSBIO has a total quality assurance policy. If any product does not meet the specification stated in the product data sheet will be replaced at no charge. This warranty limits AMSBIO liability to the replacement of the product only. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by AMSBIO. AMSBIO shall not be liable for any incidental, direct, special or consequential damages, even if it is aware of the possibility of such damages.



#### UK & Rest of World

184 Milton Park, Abingdon  
OX14 4SE, Oxon, UK  
Tel: +44 (0) 1235 828 200  
Fax: +44 (0) 1235 820 482

#### Switzerland

Centro Nord-Sud 2E  
CH-6934 Bioggio-Lugano  
Tel: +41 (0) 91 604 55 22  
Fax: +41 (0) 91 605 17 85

#### Deutschland

Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
Tel: +49 (0) 69 779099  
Fax: +49 (0) 69 13376880

#### North America

23591 El Toro Rd, Suite #180  
Lake Forest, CA 92630  
Tel: + 1 800 987 0985  
Fax: + 1 949 265 7703

amsbio

info@amsbio.com

www.amsbio.com  
AMS Biotechnology