

PRODUCT INFORMATION

Neutral Protease NB

General Neutral protease from *Clostridium histolyticum* is a proteolytic enzyme and therefore can be used for isolation of various cell types by tissue dissociation.

Description Neutral Protease NB is chromatographically purified and largely free from collagenolytic activity (PZ U).

Specification Neutral protease activity ≥ 0.50 U/mg (DMC)

Application Neutral Protease NB is, mostly in combination with Collagenase NB 1 Premium Grade (Cat. No. S1745503) or another Collagenase NB product, suitable for cell isolation from various tissue types.

If a GMP conforming product is required, Neutral Protease NB GMP Grade (Cat. No. N0002936) is recommended. This product has comparable enzymatic activities to Neutral Protease NB.

Storage conditions Neutral Protease NB is provided as a lyophilized powder. It should be stored at +2 to +8 °C in a dry environment. Under these conditions the product is stable until the minimum shelf life stated on the certificate of analysis if repeated opening and closing of the vial is avoided.

For storage of solutions please refer to "Stock solution".

Documents For each lot a specific certificate of analysis is provided. A certificate of origin is available.

Product size

Product	Cat.No.	Size (U)
Neutral Protease NB	S3030111	50
	S3030112	100

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
 184 Park Drive, Milton Park
 Abingdon OX14 4SE, U.K.
 T: +44 (0) 1235 828 200
 F: +44 (0) 1235 820 482



North America
 1035 Cambridge Street,
 Cambridge, MA 02141.
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218



Germany
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880



Switzerland
 Via Lisano 3
 (CP.683)
 6900 Massagno
 T: +41 (0) 91 604 55 22
 F: +41 (0) 91 605 17 85

PRODUCT INFORMATION

Neutral Protease NB

Instructions for use:

General

Neutral Protease NB is, in combination with Collagenase NB 1 Premium Grade, especially suited for isolation of islets of Langerhans from human pancreas, and also for cell isolation from several other human or animal tissues.

Tissue dissociation

Recommended starting concentrations for isolation of islets of Langerhans from human pancreas:

Collagenase NB 1 Premium Grade: 15 – 20 PZ U/g tissue

Neutral Protease NB: 0.8 – 1.5 DMC U/g tissue

Neutral Protease NB is provided in vials with 50 DMC U or 100 DMC U.

In general, the appropriate neutral protease concentration depends on tissue type and origin as well as on the isolation procedure.

Neutral protease activity is at an optimum at 37 °C.

Stock solution

A stock solution of Neutral Protease NB is prepared by dissolving the enzyme in water. The enzyme solution must be constantly stored on ice.

Since neutral protease depends on calcium, absolutely no calcium chelating agents (e.g. EDTA) should be present at all.

Reconstituted Neutral Protease NB can be 0.22 µm filtered, aliquoted and stored at -20 °C. Aliquots are stable for 1 year if repeated freezing and thawing is avoided.

For 0.22 µm filtration filters with low protein-binding properties (e.g. cellulose acetate, PVDF, or PES) are recommended.

Working solution

To prepare a working solution, the stock solution is diluted with water to achieve the required neutral protease concentration. The working solution must be constantly stored on ice until use.

If Neutral Protease NB solution is mixed with Collagenase NB solution, the blend should be used immediately.

Inactivation and inhibitors

The dissociation process can be reduced, e.g. by cooling down or dilution of the enzyme solution.

Inhibitors of neutral protease include chelating agents like EDTA.

Important note

Neutral Protease NB is not intended for direct application in humans.

N0056794 001

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World

184 Park Drive, Milton Park
 Abingdon OX14 4SE, U.K.
 T: +44 (0) 1235 828 200
 F: +44 (0) 1235 820 482



North America

1035 Cambridge Street,
 Cambridge, MA 02141.
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218



Germany

Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880



Switzerland

Via Lisano 3
 (CP.683)
 6900 Massagno
 T: +41 (0) 91 604 55 22
 F: +41 (0) 91 605 17 85