

ATP Assay Kit

Rapid bioluminescent determination of ATP

DESCRIPTION

Adenosine 5'-triphosphate (ATP) is the chemical energy for cellular metabolism and is often referred to as "energy currency" of the cell. ATP is produced only in living cells during photosynthesis and cellular respiration and consumed in cellular processes including biosynthetic reactions, motility and cell division. It is a key indicator of cellular activity and has been utilized as a measure of cell viability and cytotoxicity in research and drug discovery.

Biochain's ATP Assay Kit provides a rapid method to measure intracellular ATP. The single working reagent lyses cells to release ATP, which, in the presence of *luciferase*, immediately reacts with the Substrate *D-luciferin* to produce light. The light intensity is a direct measure of intracellular ATP concentration.

Luciferase

ATP + *D-luciferin* + O₂ → oxyluciferin + AMP + PP_i + CO₂ + light

This non-radioactive, homogeneous cell-based assay is performed in microplates. The reagent is compatible with all culture media and with all liquid handling systems for high-throughput screening applications in 96well and 384-well plates.

APPLICATIONS

ATP determination in cells and other biological samples.

KEY FEATURES

Safe. Non-radioactive assay.

Sensitive and accurate. As low as 0.1 μM ATP or 40 cells can be quantified.

Homogeneous and convenient. "Mix-incubate-measure" type assay. No wash and reagent transfer steps are involved.

Robust and amenable to HTS: Z' factors of > 0.5 are routinely observed in 96-well and 384-well plates. Can be readily automated on HTS liquid handling systems for processing thousands of samples per day.

KIT CONTENTS (100 tests in 96-well plates)

Assay Buffer: 10 mL

Substrate: 120 μL

ATP Enzyme: 120 μL

Standard: 100 μL 3 mM ATP

Storage conditions: store all reagents at -20°C. This kit is shipped on dry ice. Shelf life of at least 6 months.

Precautions: reagents are for research use only. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to Material Safety Data Sheet for detailed information.

PROCEDURES

1. **Standard Curve.** Prepare 1000 μL 30 μM ATP Premix by mixing 10 μL 3 mM Standard and 990 μL distilled water (for cell culture samples dilute ATP in culture media). Dilute standard as follows. Transfer 100 μL standards into wells of a white opaque 96-well plate.

No	Premix + H ₂ O/media	Vol (μL)	ATP (μM)
1	150 μL + 0 μL	150	30
2	120 μL + 30 μL	150	24
3	90 μL + 60 μL	150	18
4	60 μL + 90 μL	150	12
5	45 μL + 105 μL	150	9
6	30 μL + 120 μL	150	6
7	15 μL + 135 μL	150	3
8	0 μL + 150 μL	150	0

Samples. Add 100 μL sample per well in separate wells. For tissue samples, homogenize 20 mg sample in 200 μL of cold phosphate-buffered saline, spin at 12,000 g for 5 min to pellet any debris. Transfer 1-100 μL supernatant to each well and bring the volume to 100 μL with PBS. Test several doses of the sample and choose the readings that are within the standard curve range for ATP calculation.

For cell cultures, plate cells (100 μL/96 well plate, 25 μL/384well plate) in white opaque tissue culture plates. If desired, add 5 μL test compounds and controls dissolved in PBS or culture medium per well. Rock plate lightly to mix and incubate for desired period of time (e.g. overnight).

2. **Assay.** Bring Assay Buffer and Substrate to room temperature. Thaw enzyme on ice or at 4°C. Fresh Reconstitution is recommended. Store unused reagents including the enzyme at -20°C.


For each 96-well, mix 95 μL Assay Buffer with 1 μL Substrate and 1 μL ATP Enzyme. Add 90 μL Reconstituted Reagent to each well.

For each 384-well, mix 30 μL Assay Buffer with 0.3 μL Substrate and 0.3 μL ATP Enzyme. Add 25 μL Reconstituted Reagent to each well. Mix by tapping the plate. Incubate for 10 minutes at room temperature.

3. **Read luminescence** on a luminometer. For most luminometers (Berthold Luminometer, L.J.L Analyst, Top Count, MicroBeta Counters, CLIPR and LeadSeeker), integration time of 0.1 to 5 sec is appropriate.

GENERAL CONSIDERATIONS

Signal stability. After adding the Reconstituted Reagent, the luminescence signal is stable for about 20 min and decreases slowly thereafter.



ATP Standard Curve in Water

LITERATURE

- [1]. Kangas L, et al. (1984). Bioluminescence of cellular ATP: a new method for evaluating agents in vitro. *Medical Biology*, 62: 338 - 343.
- [2]. Zhelev Z, et al (2004). Phenothiazines suppress proliferation and induce apoptosis in cultured leukemic cells without any influence on the viability of normal lymphocytes. Phenothiazines and leukemia. *Cancer Chemother Pharmacol*. 53(3):267-75.
- [3]. Ingram PR, et al (2004). A comparison of the effects of ocular preservatives on mammalian and microbial ATP and glutathione levels. *Free Radic Res*. 38(7):739-50.



UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

United States

23591 El Toro Rd, Suite #167
Lake Forest, CA 92630
Tel: +1 800 987 0985
Fax: +1 949 265 7703