

Preparation for 3D cell culture on alvetex[®]

1. TERA2.cl.SP12 cells^{1,2} were routinely maintained in T75 flasks.




Figure 1. Phase contrast micrographs of TERA2.cl.SP12 cells grown in conventional 2D culture plates. Images show cells at low (left) and high (right) confluency. Scale bars: 100 μ m

2. Complete media consisted of: Dulbeccos Modified Eagles Medium (DMEM) supplemented with 10% v/v heat-treated FBS, 2 mM L-glutamine and 100 U/ml Penicillin/Streptomycin.
3. Cells were harvested by trypsinisation and centrifuged for 5 minutes (1000 rpm). The supernatant was discarded and the cell pellet was re-suspended in an appropriate volume of media for cell counting by Trypan Blue.
4. Cells were re-suspended at a concentration of 4×10^6 cells/ml for seeding.
5. Alvetex[®] 6-well inserts were prepared for seeding by dipping in 70% ethanol and washed twice with 7 ml of media per well.
6. 125 μ l of the cell suspension was added to the centre of the alvetex[®] disc, which was equivalent to 500,000 cells per well.
7. The plate was incubated 30 min at 37 °C with 5% CO₂ to allow the cells to settle into the scaffold.
8. 10 ml of media was added to each well taking care not to dislodge cells from alvetex[®].
9. Plates were re-incubated and maintained by complete media exchange after every 2-3 days.

¹ Przyborski S.A. (2001). Isolation of human embryonal carcinoma stem cells by immuno-magnetic sorting. *Stem Cells*, **19**, 500-504.

² Stewart R., Christie V. & Przyborski S.A. (2003). Manipulation of human pluripotent embryonal carcinoma stem cells and the development of neural subtypes. *Stem Cells*, **21**, 248-256.



UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centre Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

United States

23591 El Toro Rd, Suite #147
Lake Forest, CA 92630
Tel: +1 800 987 0985
Fax: +1 949 265 7703




Figure 2. Brightfield micrographs showing the structure of TERA2.cl.SP12 cells cultured for 7 days on 22 mm diameter alvetex® discs presented in 6-well inserts in 6-well plates. Cells were fixed, embedded in paraffin wax, sectioned (10 µm) and counterstained with haematoxylin and eosin.





Figure 3. Biochemical analysis of cell viability using a standard MTT assay. Data from 3 sample replicates of TERA2.cl.SP12 cells are shown ($n=3$, mean \pm SD). Cells were cultured for 3 days on 22 mm alvetex® discs presented in the 6-well inserts in 6-well-plate format.