

Introduction:

The following protocol outlines a method for the retrieval of cells cultured in Alvetex® Scaffold. Example data were obtained using this protocol to extract HepG2 hepatocyte cells cultured in Alvetex® Scaffold for 7 days in 6-well inserts (AMS.AVP004-32) in the Well Insert Holder in a Deep Petri Dish (AMS.AVP015-2) format. Two methods of cell extraction may be employed; a standard method and a syringe-based method which produces higher average cell yields. Both methods demonstrate that partial retrieval of cells from Alvetex® Scaffold is possible.

Methods:

1.0. Standard Enzymatic Method

- 1.1. Unclip inserts and carefully remove the Alvetex® Scaffold discs using flat-ended forceps.
- 1.2. Gently wash each disc in PBS, and transfer to a new 6-well plate.
- 1.3. 3 ml of 0.25 % Trypsin-EDTA (Sigma, T4049).
- 1.4. Incubate plate at 37 °C, 5 % CO₂ on a shaking platform set to 100 rpm, for 15 minutes.
- 1.5. Transfer the resulting cell suspension to a 15 ml centrifuge tube.
- 1.6. Add 3 ml medium to the well containing the Alvetex® Scaffold membrane and gently triturate to remove residual detached cells.
- 1.7. Combine the two solutions in the 15 ml centrifuge tube to neutralise the trypsin.
- 1.8. Centrifuge at 1000 rpm for 5 minutes, to pellet the cells.
- 1.9. Resuspend the cell pellet in an appropriate volume of medium for downstream processes.

2.0. Enzyme and Syringe Method

- 2.1. Unclip inserts and carefully remove the Alvetex® Scaffold discs using flat-ended forceps.
- 2.2. Gently wash each disc in PBS, and transfer to a new 6-well plate.
- 2.3. 3 ml of 0.25 % Trypsin-EDTA (Sigma, T4049).
- 2.4. Incubate plate at 37 °C, 5 % CO₂ on a reciprocating shaker set to 100 strokes/min, for 15 minutes.

- 2.5. Transfer the resulting cell suspension to a 15 ml centrifuge tube.
- 2.6. Transfer the Alvetex® Scaffold membrane to the bottom of a 30 ml plastic syringe barrel. Add 3 ml medium and depress the plunger slowly to force residual detached/loose cells from the scaffold.
Note: Care should be taken not to force the solution through the membrane too quickly as this is liable to induce considerable shear stress upon the cells.
- 2.7. Combine the two solutions in the 15 ml centrifuge tube to neutralise the Typsin.
- 2.8. Centrifuge at 1000 rpm for 5 minutes, to collect the cells.
- 2.9. Resuspend the cell pellet in an appropriate volume of medium required for the downstream processes.

Example: HepG2 cells grown for 7 days in Alvetex® Scaffold 6-well insert in a Petri dish format. HepG2 cells (ATCC, HB-8065) were routinely maintained in T-75 flasks. HepG2 complete media consisted of: MEM media (Gibco, 21090) supplemented with 10 % v/v FBS, 2 mM L-glutamine and 100 U/ml Penicillin/Streptomycin. Cells were seeded onto Alvetex® Scaffold discs in 6-well inserts (AMS.AVP004-32) in Well Insert Holders in Deep Petri dishes (AMS.AVP015-2), at a density of 1×10^6 cells in 100 μ l media per insert. After settling for 1 hour in an incubator (5 % CO₂, 37 °C), complete media was carefully added (70 ml per Petri dish). Cultures were maintained for 7 days, with media changes on days 2 and 4. After 7 days, cultures were processed according to the protocol described above.

The total number, and number of viable cells retrieved were counted using Trypan Blue exclusion method. Alvetex® Scaffold discs were analysed post-retrieval by MTT assay (Figure 1.) and histological staining to estimate cell numbers remaining in the 3D scaffold (Figure 2.).




Figure 1. Biochemical analysis of cell viability using a standard MTT assay on Alvetex® Scaffold discs post-cell extraction using trypsin-EDTA. Data from 3 sample replicates of HepG2 cultures are shown ($n=3$, mean \pm SE).

On the basis of residual in-scaffold MTT assay data, 12 % of cells were retrieved (compared to untreated control discs) using the standard protocol described above. When the syringe-based method was employed, cell recovery increased to 24 %. Cell viability post-retrieval was 99 % for both methods, as determined by cell counting in conjunction with Trypan Blue exclusion (data not shown).




Figure 2. Bright field micrographs showing the morphology of HepG2 cells cultured for 7 days on Alvetex[®]Scaffold in 6-well insert within a Petri dish format (x10 and x40 objective). **A:** untreated HepG2 cells and **B:** HepG2 cells after Trypsin-EDTA treatment for 15 minutes. Cells were fixed, embedded in paraffin wax, sectioned (10 μ m) and counterstained with Haematoxylin and Eosin. The effect of Trypsin action can clearly be seen on residual cells which lose contact with each other and adopt a rounded morphology.

www.amsbio.com - info@amsbio.com



UK and rest of the world

Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482



DE

Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880



CH

Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85



US

Tol Free: +1 800 987 0985
Tel: +1 949 265 7717
Fax: +1 949 265 7703

amsbio