

Introduction

Immunofluorescence uses the recognition of cellular targets by fluorescent dyes or antigen-specific antibodies coupled to fluorophores. Depending on the antibody or dye used, proteins, lipids and DNA can be visualised within individual cells and tissues. Alvetex® Scaffold can easily be processed like a standard tissue sample, allowing established immunofluorescence methods to be followed with excellent results. An example protocol is outlined below.

Materials and Methods

1.0 Materials Required for Fixing, Embedding and Sectioning

- 4 % Paraformaldehyde (PFA) fixative (for full method refer to 'Histology Series Part 1. Choosing the Right Fixative to Preserve 3D Cell Cultures' located at www.amsbio.com/alvetex.aspx).
- Dehydration ethanols (30 %, 50 %, 70 %, 80 %, 90 %, 95 %, 100 %)
- HistoClear (*Note: the use of HistoClear is recommended over the use of Xylene*).
- Paraffin wax
- Basic equipment: spade forceps to handle 3D Alvetex® Scaffold cultures, scalpel (to trim / cut scaffold if required), disposable pipettes, embedding moulds, convection oven, microtome.

2.0 Materials Required for Immunofluorescence

- 1x citrate buffer pH 6.0
- Blocking buffer (5 % Normal goat serum; 1 % bovine serum albumin; 0.2 % Triton X-100 in PBS)
- Permeabilisation solution (1 % w/v Triton X-100 in PBS)
- Primary and secondary antibodies (specific to experimental analysis)
- Microwave oven (if antigen retrieval is required)
- Humidified chamber
- Vectashield/DAPI mountant

3.0 4 % PFA Fixation of Scaffolds

- 3.1. Aspirate off the medium and carefully wash the 3D culture in PBS twice.
- 3.2. If using well inserts in Petri dishes, remove inserts from holder and place in a conventional 6-well plate. To fix, add in 5 ml 4 % PFA fixative at 4 °C. Fix specimens at 4 °C for a minimum of 12 hours but no longer than 24 hours.
- 3.3. Aspirate off the fixative, wash 3 times using 5 ml PBS to thoroughly remove excess fixative, discarding the waste liquid.
- 3.4. If using inserts, remove Alvetex® Scaffold from the well insert either by separating the two parts of the device or using a scalpel to cut out the scaffold. At this time samples can be transferred to a tissue processing cassette to minimize direct handling and damage to the 3D culture.
- 3.5. Aspirate off the PBS and add 5 ml of 30 % ethanol. Leave to equilibrate for at least 15 minutes. Aspirate off the ethanol and discard.
- 3.6. Repeat with 50 %, 70 %, 80 %, 90 % and then 95 % ethanol. A gradual dehydration of the sample will result in less tissue shrinkage. Material can be stored in 95 % ethanol prior to paraffin embedding.

4.0 Paraffin Embedding and Sectioning

- 4.1. Fully dehydrate specimens stored in 95 % ethanol by replacing with 100 % absolute ethanol for at least 30 minutes.
- 4.2. Aspirate off the ethanol and replace with HistoClear for at least 30 minutes. *(Note: HistoClear will dissolve certain types of plastic, therefore it is best to process samples in glass tissue processing cassettes).*
- 4.3. Replace the HistoClear with a 50:50 solution of HistoClear and molten paraffin wax (60 °C) mix and incubate in a convection oven at 60 °C for 30 minutes.
- 4.4. Replace the HistoClear:wax mix with 100 % molten wax and incubate at 60 °C for a further 60 minutes.
- 4.5. Transfer the scaffold to plastic embedding moulds and orientate into the required embedding position, with plane of section in mind. Embed in molten wax.
- 4.6. Allow wax to cool and set at room temperature for 1-2 hours or overnight.
- 4.7. Once the wax has hardened, remove the wax embedded block from the plastic mould. The sample is now ready for sectioning on a suitable microtome (e.g. Leica RM2125).
- 4.8. Following the microtome manufacturers' instructions throughout, align the block correctly with the microtome blade and proceed to cut 10 µm sections of the sample block.

- 4.9. Transfer sections to a slide water bath (40 °C), floating them on the surface of the water to enable them to flatten out.
- 4.10. Transfer selected sections to slides by flotation. Superfrost Plus slides (Thermo, 4951PLUS4) are recommended since sections adhere to these slides well.
- 4.11. Place on a slide drier and leave overnight. The sections should now be ready for immunofluorescence.

5.0 Immunofluorescent Analysis

- 5.1. Deparaffinise sections in HistoClear for 10 minutes. Handle samples carefully to avoid loss of section from slide.
- 5.2. Hydrate specimen through a graded series of ethanols (100 %, 90 %, 70 %) with 5 minutes incubation in each solution.
- 5.3. If antigen retrieval is required, place slides in 1x citrate buffer (pH 6) and microwave (800 W) for 6 minutes. If antigen retrieval is not required go to point 5.6.
- 5.4. Leave to stand outside microwave for 1 minute then microwave for a further 3 minutes at 800 W.
- 5.5. Allow to cool for 20 minutes and wash in PBS for 10 minutes. Repeat three times.
- 5.6. Treat cells with permeabilisation solution for 15 minutes.
- 5.7. Aspirate off permeabilisation solution, replace with blocking buffer and incubate for 15 minutes.
- 5.8. Incubate with primary antibody diluted in blocking buffer (this will be specific to the antigen chosen and will be used at a pre-determined concentration) overnight at 4 °C in a humidified chamber, or as stated in the datasheet for the specific antibody used.
- 5.9. Wash 3 times for 10 minutes in PBS.
- 5.10. Incubate with secondary antibody diluted in blocking buffer (this will be specific to the antigen chosen and will be used at a pre-determined concentration) for 2 hours in the dark at room temperature.
- 5.11. Wash 3 times for 10 minutes with PBS.
- 5.12. Mount in Vectashield/DAPI. This solution simultaneously mounts the specimen and stains cell nuclei.

- 5.13. Seal the perimeter of the cover slip with nail varnish and dry in the dark.
- 5.14. Store slides in the dark at 4 °C until ready for inspection using a fluorescence microscope equipped with the appropriate filters.

Example results

Figure 1. below shows an example of a 3D Alvetex® Scaffold culture fixed with 4 % PFA, embedded in paraffin wax, sectioned and probed with Ki67 antibody using immunocytochemistry and standard fluorescence microscopy.




Figure 1 Protein expression in 3D cell cultures can be localised using immunocytochemistry and fluorescence microscopy. Three corresponding images from the same region of Alvetex® Scaffold show: **(A)** cell morphology by phase microscopy; **(B)** expression of cell nucleus marker DAPI; **(C)** expression of cell proliferation marker Ki67.

AMSBIO is the global source for alvetex®.
alvetex® is a registered trade mark of and manufactured by Reinnervate.

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
AMS Biotechnology (Europe) Ltd
184 Park Drive, Milton Park
Abingdon, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
amsbio LLC
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
AMS Biotechnology (Europe) Ltd
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
AMS Biotechnology (Europe) Ltd
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85