


Extragel Matrix Manual



AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada 

1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU 

Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World 

184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482






AMS Biotechnology (Europe) Ltd
Switzerland 

Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland 

T: +49 (0) 69 779099
F: +49 (0) 69 13376880

CONTENTS

-  01# **Product Introduction**
-  01# **Product Characteristics**
-  02# **Operation Method**
-  05# **Practical Application Cases**
-  08# **Reference**

Product Introduction

The basement membrane is a matrix under the basal surface of epithelial cells of animals. AMSBIO's Extragel Matrix is a reconstituted matrix hydrogel formed by basement membrane components extracted from mouse tumor tissues. This matrix hydrogel is mainly composed of laminin, collagen IV, and heparan sulfate proteoglycans (Kleinman et al. 1986). Besides, it contains various growth factors, such as epidermal growth factor (EGF), platelet-derived growth factor (PDGF), nerve growth factor (NGF), basic fibroblast growth factor (FGF-2), transforming growth factor- β (TGF- β), and insulin-like growth factor (IGF) (Vu-kicevic et al. 1992).

Product Characteristics

Extragel Matrix is liquid at 4°C but gelled when heated to 37°C. This transformation phenomenon is reversible. It can be liquefied again when it is stored at 4°C overnight. (Tip: It is recommended to store the Extragel Matrix in an ice box in a refrigerator at 4°C to realize the full liquefaction of the reconstituted matrix hydrogel.)

Dispense Extragel Matrix into appropriate aliquots when first use. Stable for 2 years when stored at -80°C

PRODUCT APPLICATION

Extragel Matrix can be applied to the growth, differentiation, metabolism and toxicology of organoids, in vivo and in vitro angiogenesis experiments, and tumor formation experiment of immunodeficient mice.

PRECAUTIONS

Extragel Matrix would start solidifying after the temperature is higher than 10°C, so the operation should be performed on ice. The matrix hydrogel can be dissolved in basic culture medium pre-chilled at 4°C, and the organoid can be released from the Extragel Matrix.

Operation Method

ORGANOID CULTURE (1 HOUR)

1. Thaw the Extragel Matrix (-OM-1 or -OM-1-ph or -OM-3 or -OM-3-ph) in refrigerator at 4°C overnight.
 2. Preheat the 24-well plate in cell culture incubator.
 3. Prepare aliquots of Extragel Matrix using pre-chilled tips.
 4. Prepare the single cell pellet with 1×10^6 cells derived from patients or animal tissue, centrifuged at 300 g for 5 minutes.
 5. Mix the cell pellet with 50 μ L Extragel Matrix thoroughly.
 6. Add the mixture into the well of plate (50 μ L per well).
 7. Keep the plate in incubator for 10 minutes, flip and keep after another 5 minutes.
 8. Add 500 μ L culture medium to the well with matrix and cells.
 9. Change the medium every 3 days.
-


TUMOR FORMATION EXPERIMENT IN IMMUNODEFICIENT MICE (1 HOUR)

1. Mix 1×10^6 cells with Extragel Matrix (volume ratio=1:1) (OM-2 or OM-2-ph or OM-5 or OM-5-ph). The total volume of mixture for each injection should be 100 μ L at least.
2. Inject the mixture of cells and Extragel Matrix with a 1 mL syringe mounted with a 16 g needle and inject the mixture into subcutaneous tissues or thigh muscles.

Operation Method

● ANGIOGENESIS EXPERIMENT (1 HOUR)

1. Starving cells: replace the complement medium with starving medium: DMEM medium containing 0.2% FBS, 2 mM L-glutamine, 1 mM sodium pyruvate, 100 U/mL penicillin, and 100 µg/mL streptomycin for the incubation for 24h.
2. Uniformly spread Extrigel Matrix (OM-1 or OM-1-ph) on the 96-wellplate.
(Note: The pipette should be pre-cooled for 30 minutes. The operation shall be performed on ice to avoid the gelation of Extrigel Matrix and the generation of bubbles.)
Incubate the 96-well plate in a cell incubator for 30 minutes to solidify the
3. Extrigel Matrix.
Passage HUVECs and count them.
4. Add 200 µL of HUVEC suspension (containing 5×10^4 cells) to a 96-well plate
5. containing Extrigel Matrix, and put the 96-well plate in an incubator.
The vascular network structure will be formed in 3-12 hours.
6. At the optimal time for the vascular network formation, remove the culture
7. medium properly, stain it with the culture medium containing living cell dye 1/1000 Calcein AM (green), and take photos for records with a microscope.



NEURITE OUTGROWTH ASSAY



1. Prepare the rat neural stem cell (NSC) pellet with 1×10^5 cells, centrifuged at 300 g for 5 minutes.
2. Mix the cell pellet with 50 μ L Extragen Matrix (OM-1 or OM-1-ph), thoroughly.
3. Keep the plate in incubator for 10 minutes, flip and keep after another 5 minutes.
4. Add 500 μ L culture medium to the well with matrix and cells.
5. Change the medium every 3 days.

Notes: NSC culture medium contains Neurobasal medium, 2% B27, 2mM L-glutamine, 5% FBS, 20 ng/mL EGF, 20 ng/mL bFGF, 100 U/mL penicillin and 100 μ g/mL streptomycin. Neurites will be observed from day 2.

hESC AND iPSC CULTURE WITHOUT FEEDER LAYER (1 HOUR)

1. Thaw the Extragen Matrix (OM-4 or OM-4-ph) in refrigerator at 4°C overnight.
2. Mix the Extragen Matrix up&down with Pipette for 5 times.
3. Shortly spin the tube with benchtop centrifuge to remove any bubbles.
4. Preheat the well plate in cell culture incubator.
5. Prepare the diluted Extragen Matrix using basal medium at a ratio of 1:100.
6. Coating the well plate with diluted Extragen Matrix with volume of 300 μ L/cm² for 30 minutes at room temperature.
7. Remove diluted Extragen Matrix after coating and seed the hESC or iPSC solution.
8. Transfer the plate to incubator.

Practical Application Cases



Practical Application Cases




Figure 3. Growth of human colorectal cancer organoids derived from patients in the matrix hydrogel of Company A and Extragal Matrix for 13 days, respectively.







Figure 4. Characterization of hepatocellular carcinoma organoids derived from patients in the matrix hydrogel of Company A and Extragal Matrix, with immunostaining by DAPI (nuclus), HNF4a (a hepatic biomarker) and GPC3 (hepatocellular carcinoma biomarker).

Practical Application Cases



Practical Application Cases



Reference

1. Kleinman HK, et al, Basement membrane complexes with biological activity. *Biochemistry* 25: 312 (1986).
2. Vukicevic, Slobodan, et al. Identification of multiple active growth factors in basement membrane Matrigel suggests caution in interpretation of cellular activity related to extracellular matrix components. *Experimental cell research* 202: 1 (1992).
3. Guillen, K P, et al. A human breast cancer-derived xenograft and organoid platform for drug discovery and precision oncology. *Nature Cancer* 3: 232 (2022).