

## PROTOCOLS FOR BIOTIN-HABP (BOVINE NASAL CARTILAGE)

### Immunohistochemistry for Hyaluronic acid using B-HABP

#### Reagents

Product	Source	Cat.No.
<b>Biotinylated Hyaluronic Acid Binding Protein</b>	<b>Seikagaku</b>	<b>400763</b>
<b>Hyaluronidase (<i>Streptomyces hyalurolyticus</i>)</b>	<b>Seikagaku</b>	<b>100740-1</b>
<b>Chondroitinase ABC Protease Free (<i>Proteus vulgaris</i>) *</b>	<b>Seikagaku</b>	<b>100332-1A</b>
Trypsin *	Sigma	
Avidin solution/ Biotin solution	Vector	
Fluorophore conjugated Streptavidin	Jackson	
* For proteoglycan digestion when HA is masked.		

#### Experimental procedure

##### 1. Pretreatment with Hyaluronidase (Negative control)

Treat with reaction buffer of Hyaluronidase (100mM Sodium acetate buffer, pH6.0) for 15min at 37°C.

No wash

Treat with Hyaluronidase (200TRU/mL; 100mM Sodium acetate buffer, pH6.0) for 2hrs at 60°C.

Wash with PBS

##### 2. Blocking endogenous avidin biotin activity

Treat with avidin solution for 20min at RT.

Wash

Treat with biotin solution for 20min at RT.

Wash

##### 3. Treat with 0.1% BSA solution for 1hr at RT.

Wash


##### 4. Treat with Biotinylated HABP (2µg/mL) for 1-2hrs at RT.

Wash

##### 5. Treat with Fluorophore conjugated Streptavidin for 15min at RT.

**Example: Histochemical staining of hyaluronan (Texas red) in rabbit cornea.**

**N:** Normal cornea; **1, 3, 14** and **28 days:** corneas 1, 3, 14 and 28 days after wounding



### ELISA Procedure for Hyaluronic Acid using B-HABP

1. Coat Costar plate with HA (1 ug) or use BSA-HA plates
2. Wash with T-PBS buffer (200ul) 3 times
3. Add bHABP (e.g. 50 ng/mL) [T-PBS/1% BSA as a negative control]
4. Incubate covered by Parafilm at 37 C for 1 hour
5. Discard excess bHABP
6. Wash with T-PBS buffer (200ul) 3 times
7. Add 100ul HRP-sAv solution (ImmunoPure® Streptavidin, Horseradish Peroxidase Conjugated, 1mg)
8. Incubate covered by Parafilm at 37 C for 1 hour
9. Discard excess HRP-sAv
10. Wash with T-PBS buffer (200ul) 3 times
11. Add 100 ul TMB substrate
12. Incubate covered by Parafilm at 37 C for 15 min (prevent light exposure)
13. Add 100ul 1N HCl
14. Measure the endpoint at A450 minus A630nm