



Mimetix[®] Protocol: Cell Viability Assays

CellTiter Blue assay

1. Add 10 μ l of CellTiter Blue reagent per well of a 96-well plate.
 2. Incubate for 90 minutes at 37°C.
 3. Place plate into plate reader and shake vigorously for several seconds.
 4. Measure fluorescence at 580-620 nm (excite at 530-570 nm) according to the manufacturer's instructions.
- Please note: The Mimetix scaffold does not cause background fluorescence.

CellTiter Glo assay

1. Add 100 μ l of CellTiter Glo reagent per well of a 96-well plate.
 2. Incubate for 10 minutes at 37°C.
 3. Place plate into plate reader and shake vigorously for several seconds.
 4. Measure luminescence according to the manufacturer's instructions.
- Please note: It is not necessary to use the "CellTiter Glo 3D" assay kit. The scaffold thickness is only 50 μ m and thereby compatible with the standard CellTiter Glo assay kit.

MTT/MTS/WST-1 assay (various suppliers)

1. Add 10-20 μ l of the MTT/MTS/WST-1 assay reagent and incubate for 2-4 hours (refer to the respective instructions for further guidance).
2. MTT only: Carefully aspirate the supernatant and add 200 μ l DMSO or other recommended solvent.
3. Shake vigorously for 2-5 min until a homogeneous solution is obtained and all crystals have dissolved (MTT only).
4. Transfer 100-200 μ l of the content of each well into a fresh, clear plastic tissue culture plate: the scaffold gives rise to a background signal in absorbance-based assays!
5. Measure absorbances at 540 nm (MTT), 490 nm (MTS) and 450 nm (WST-1), respectively.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
AMS Biotechnology (Europe) Ltd
184 Park Drive, Milton Park
Abingdon, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
amsbio LLC
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
AMS Biotechnology (Europe) Ltd
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
AMS Biotechnology (Europe) Ltd
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85