



Mimetix[®] Protocol: Cell Visualisation

Imaging cells within the Mimetix scaffold with fluorescent or confocal microscopy

Confocal microscopy is a useful tool to obtain high-resolution images of cells within the Mimetix[®] scaffold and to investigate their depth profile. But cells can be easily observed using a “classic” fluorescence microscope too. In the example protocol in the next paragraph, nuclei are stained with TO-PRO[®] 3 and F-actin with Alexa Fluor[®] 488 Phalloidin. The Mimetix scaffold used here is a bespoke product and has Rhodamine 6G incorporated into the fibres. Please [contact us](#) if this could also be of interest for your research.

For a typical experiment using Mimetix[®] in a 12-well plate format, prepare a cell suspension at a concentration of 10^5 cells/mL for short experiments (1-3 days) or 5×10^4 cells/mL for longer experiments (up to 7 days) and add 1 mL into each well. Allow cells to adhere for 24 hours or overnight.

- To fix the cells in the scaffolds, aspirate medium and add 0.5 mL pre-warmed 4% paraformaldehyde in PBS to each well.
- Incubate for 15-30 min, then aspirate paraformaldehyde, wash 2x with PBS and add 0.5 mL 0.1% Triton X-100 in PBS to each well to permeabilize cell membranes (only necessary when performing Alexa Fluor[®] 488 Phalloidin staining).
- Incubate for 5 min, then aspirate Triton X-100, wash 2x with PBS and add 0.5 mL 0.5% BSA in PBS to each well to act as a blocking agent (only necessary when performing Alexa Fluor[®] 488 Phalloidin staining).
- Incubate for 30-60 min, then aspirate 0.5% BSA in PBS, wash 2x with PBS and add 200 μ L Alexa Fluor[®] 488 Phalloidin working solution into each sample.
- Incubate for 30-60 min, then aspirate Alexa Fluor[®] 488 Phalloidin, wash 2x with PBS and add 1 mL TO-PRO[®] 3 working solution (1:10,000 in PBS) to each dish.
- Incubate for 10 min, then aspirate TO-PRO[®] 3, wash 2x with PBS and 1x with water to prevent the formation of salt crystals.
- Remove scaffold discs from wells and place on a glass slide, cell-seeded side facing up. Add a drop of mounting medium, place cover slip on top (dry off excess mounting medium with paper towel) and seal edges with nail varnish.
- Image stained cells within Mimetix[®] by exciting Alexa Fluor[®] 488 Phalloidin at 488 nm using an Argon laser and detecting at 500-550 nm, Rhodamin 6G-labelled fibres at 543 nm using a He-Ne laser and detecting at 550-615 nm, and TO-PRO[®] 3 at 633 nm using a He-Ne laser and detecting at 650-710 nm.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
AMS Biotechnology (Europe) Ltd
184 Park Drive, Milton Park
Abingdon, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
amsbio LLC
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
AMS Biotechnology (Europe) Ltd
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
AMS Biotechnology (Europe) Ltd
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85