



Application Note

Isolation of Porcine Hepatocytes with Collagenase NB 8 Broad Range

Product information Collagenase NB 8 Broad Range

Collagenase NB 8 Broad Range is designed for dissociation of different tissues to isolate various cell types, e.g. porcine hepatocytes. Purified Collagenase NB 8 Broad Range is characterized by an increased collagenase activity and reduced levels of other proteases and bacterial endotoxins.

Specification

 $\begin{array}{ll} \mbox{Collagenase activity} & \geq 0.9 \ \mbox{U/mg (PZ*)} \\ \mbox{Neutral protease activity} & \leq 0.2 \ \mbox{U/mg (DMC)} \\ \mbox{Trypsin-like activities} & \leq 0.5 \ \mbox{U/mg (BAEE)} \\ \mbox{Endotoxin} & \leq 100 \ \mbox{EU/mg} \\ \end{array}$

* PZ units acc. to Wünsch

The lot-specific activities of Collagenase NB 8 Broad Range are stated on the certificate of analysis which is attached to every delivery.

Storage conditions

Collagenase NB 8 Broad Range is provided as a lyophilized powder and should be stored in a dry environment at +2 to +8 °C.

Stock solutions can be 0.2 μm filtered and stored at -20 °C for 1 year without loss of activity.

Isolation of Porcine Hepatocytes

| Recommended Starting Conditions | |
|--|------------------------------------|
| Concentration of Collagenase NB 8 Broad Range | 0.2 – 0.3 PZ U/ml* |
| Isolation Conditions | perfusion and subsequent digestion |

*P7 units acc. to Wünsch

In general the required collagenase concentration depends on tissue type and origin as well as isolation procedure. Therefore Collagenase NB concentrations stated above should be considered as starting points and progress of the digestion process should be monitored visually in order to determine the optimal conditions.

The digestion process can be stopped by adding EDTA, cooling down or diluting the digestion solution.

It is recommended to measure the concentration of the collagenase solution in PZ U/ml and not in mg/ml.

This ensures optimal reproducibility of isolations with collagenases from different batches, although Nordmark guarantees a high lot-to-lot consistency.

Ordering Information

| Enzyme | Cat. No. | Pack size |
|------------------------------|----------|-----------|
| Collagenase NB 8 Broad Range | S1745601 | 250 mg |
| Collagenase NB 8 Broad Range | S1745602 | 1 g |

Collagenase NB products

- > High cell yields and viability
- > Reliable lot-to-lot consistency
- > GMP grades available

AMSBIO | www.amsbio.com | info@amsbio.com



F: +44 (0) 1235 820 482

1035 Cambridge Street, Cambridge, MA 02141. T: +1 (617) 945-5033 or T: +1 (800) 987-0985

F: +1 (617) 945-8218



Germany Bockenh

Bockenheimer Landstr. 17/19 60325 Frankfurt/Main T: +49 (0) 69 779099 F: +49 (0) 69 13376880

