



● Cloning efficiency with coating-free protocol and pluripotency analysis of isolated hiPSC clones



1231A3 and 1210B2 (Human episomal iPSC lines established by CiRA) were cloned on iMatrix-511 or Matrigel® by the precoating or coating-free method. (A) Comparison of the cloning efficiency of hiPSCs by two different methods. Bars represent the means ± S.D. (n=3). (B, C) Analysis of pluripotent markers in isolated hiPSC clones. Expression levels of pluripotent markers were evaluated by (B) FACS and (C) ICC (Blue: DNA, Green: Tra1-60, Red: Oct3/4). Scale bars: 100 µm.

Feeder-free medium for ES/iPS cells



StemFit Technical tips

Key Points for *single-cell cloning* with *coating-free method*

Benefit 1

Superior colony-forming efficiency

Enables efficient single-cell cloning

Benefit 2

Coating-free protocol

No coating process, No incubation

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142


AMS Biotechnology (Europe) Ltd
UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland
Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

amsbio

Brief protocol for single-cell cloning by the coating-free method



1

Detach cells and resuspend in culture medium supplemented with 10 μ M Y-27632.

* Also see our Technical tips: Key points for successful single-cell passage

2

Prepare 10 mL of 10 cell/mL cell suspension by serial dilution with culture medium with 10 μ M Y-27632.

3

Add iMatrix-511 or Matrigel® to the prepared cell suspension and mix thoroughly.

Point-1

Concentration of ECMs !

ECM	Amount	Final conc.
iMatrix-511 (0.5 mg/ml)	35 μ l	1.75 μ g/ml
Matrigel®	100 μ l (x1/100)	10 μ l/ml

4

Plate 100 μ L (= 1 cell) in each well of the 96-well plate immediately.

5

Replace medium with fresh culture medium without Y-27632 at least every three days. Around day 8, select single colonies to be passaged to a 24-well plate.

<Medium Change Schedule Example>



Point-2

Extend or shorten culture period for colonies to be appropriate sizes

6

After washing the colonies with 100 μ l of PBS, detach the cells with 50 μ l of cell detaching solution and incubate at 37 $^{\circ}$ C for 10 min.

* Accutase or TrypLE™ can be used
* Incubation times may vary

7

Carefully remove cell detaching solution.

Point-3

Remove very carefully because the colonies easily detach

8

Dissociate colonies by pipetting with 100 μ l of culture medium with 10 μ M Y-27632. Transfer resuspended cells to ECM-coated 24-well plate with 400 μ l of culture medium containing 10 μ M Y-27632 immediately.

Point-4

Detach the colonies one by one as reattachment can occur soon after adding the culture medium

9

Change the medium to fresh culture medium without Y-27632 at least every three days.

amsbio