

Lipopolysaccharide from *E. Coli* 0111:B4

Catalog # 9028

For Research Use Only - Not Human or Therapeutic Use

DESCRIPTION:	Lipopolysaccharide (LPS) from <i>E. Coli</i> 0111:B4 in PBS
APPLICATION:	<p>LPS has two independent unique activities: one is triggering inflammatory reactions and the other is stimulating B-cells to produce antibodies. Use these unique activities for triggering, synchronizing, or enhancing the development of arthritis in both collagen-induced arthritis (CIA) and collagen antibody-induced arthritis (CAIA) models.</p> <p>A) In CAIA model: Use LPS to trigger arthritis in mice under the presence of sub-arthritis levels of autoantibodies to type II collagen.</p> <ol style="list-style-type: none"> 1) Inject sub-arthritis dose of polyclonal antibodies or a cocktail of monoclonal antibodies into mice on day 0, and then inject 25-50 µg of LPS by IP to trigger the development of arthritis. 2) An additional IP injection of LPS on day 7-14 re-activates arthritis. <p>B) In CIA model: Use to stimulate autoantibody production and development of arthritis in mice immunized with type II collagen.</p> <ol style="list-style-type: none"> 1) To synchronize arthritis development: Immunize mice with type II collagen/CFA on day 0 and then inject of 25-50 µg of LPS by IP on day 25-28, just before the expected onset of arthritis. 2) To stimulate antibody production and development of arthritis: Immunize mice with type II collagen/CFA on day 0 and then inject 20-50 µg of LPS by IP on day 21 or 28 to increase antibody production and the following development of arthritis. 3) To reactivate arthritis: Inject 25-50 µg of LPS by IP in mice on day 50-80, the late stage of CIA. <p>C) LPS is also used at low levels (0.01-1 µg/ml) as a positive control in vitro for stimulating the production and release of cytokines and chemokines.</p>
QUANTITY:	5 ml
FORM:	0.5 mg/ml solution
SOURCE:	<i>E. coli</i> 0111:B4
STORAGE TEMPERATURE:	-20°C
STABILITY:	1 year

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
 184 Park Drive, Milton Park
 Abingdon OX14 4SE, UK
 T: +44 (0)1235 828 200
 F: +44 (0) 1235 820 482

 **North America**
 1035 Cambridge Street,
 Cambridge, MA 02141
 T: +1 (617) 945-5033 or
 T: +1 (800) 987-0985
 F: +1 (617) 945-8218

 **Germany**
 Bockenheimer Landstr. 17/19
 60325 Frankfurt/Main
 T: +49 (0) 69 779099
 F: +49 (0) 69 13376880

 **Switzerland**
 Centro Nord-Sud 2E
 CH-6934 Bioggio-Lugano
 T: +41(0) 91 604 55 22
 F: +41(0) 91 605 17 85