

Certificate of Analysis

Osteoblast Differentiation Medium was designed to promote differentiation of osteoblasts in culture. This medium is fully supplemented and contains no fetal bovine serum. Osteoblast Differentiation Medium is sterile filtered and bicarbonate buffered, and ready to be used for the differentiation of osteoblast.

Product Name: Human Osteoblast Differentiation Medium

Quantity: 250 ml

Cat. No.: AMS.SMH417D250

Quality Control:

Bioassay: Each lot was tested for the ability of osteoblasts to mineralize their extracellular matrix within 20 days when cells are cultured in this differentiation medium.

Sterility Test: Negative for bacteria, yeast and fungi.

Storage & Stability: Avoid exposure to light. Stable for 2 months at 4°C.

For research purposes only. Not for diagnostic or in vivo use.

AMSBIO | www.amsbio.com | info@amsbio.com



UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85