

Attoglow Western Blot Enhancing Kit - with Millennium Enhancer Certificate of Analysis

Product Name: Attoglow Western Blot Enhancing Kit - with Millennium Enhancer

Catalog No.: K3172120/K3172250

Lot No.: _____

Shipping and Storage Condition: Ship and store at Room Temperature.

Shelf Life: 1 year from the date of receipt under proper storage conditions

Description

Components of Attoglow Western Blot Enhancing Kit are prepared with pure chemicals according to the proprietary technology. They are used for ultra-sensitive chemiluminescent analysis of western blot. This kit provides a complete system for western blot analysis, including 1) powerful, novel proprietary technology, Millennium Enhancer, which could significantly boost western blot signals, 2) protein interaction system and 3) Attoglow blocking agent. Two different sizes are available, which are good for 1,200 cm² and 2,500 cm² of membrane, respectively.

Quality Control

- Run two samples, which contain 10 fg and 10 pg of purified rabbit muscle GAPDH protein in duplicate, on SDS-PAGE. Transfer the gel to a PVDF membrane. Cut the membrane into two pieces, each contains one set of samples. Treat one piece with Millennium Enhancer, another as an untreated control, and both blotted with anti-GAPDH antibody clone 6C5 and HRP-conjugated anti-Mouse IgG. The 10 fg-sample should give one visible band only on the treated membrane, and the band should be almost as dark as that shown by the 10 pg-sample on the untreated membrane.
- Components in this kit are free of precipitate. When testing by visible spectrophotometry, the color of the solution is comparable to a control (previous) Lot. Results from this Lot of kit are comparable to those obtained with control (previous) lot of the solution.

Components

1. 1 Attoglow Western Blot Enhancing Kit
2. Certificate of Analysis
3. One complete Manual
4. MSDS

FOR RESEARCH USE ONLY

APPROVED BY: _____



UK & Rest of World

184 Milton Park, Abingdon
OX14 4SE, Oxon, UK
Tel: +44 (0) 1235 828 200
Fax: +44 (0) 1235 820 482

Switzerland

Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
Tel: +41 (0) 91 604 55 22
Fax: +41 (0) 91 605 17 85

Deutschland

Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
Tel: +49 (0) 69 779099
Fax: +49 (0) 69 13376880

United States

23591 El Toro Rd, Suite #167
Lake Forest, CA 92630
Tel: + 1 800 987 0985
Fax: + 1 949 265 7703