

IGF-1R (pY1161)

Insulin-like growth factor 1 receptor; IGF-1R

| | |
|-----------------------------|---|
| Catalogue #: | 500-12304 |
| Lot #: | See the label |
| Size: | 100ug |
| Host: | Rabbit |
| Isotyping: | Rabbit IgG |
| Applications: | IHC, ELISA |
| Reactivity: | Hu, Ms, Rt |
| Product Preparation: | N/A |
| Background: | N/A |
| Purification: | N/A |
| Specificity: | N/A |
| Formulation: | This affinity purified antibody is supplied in sterile Tris-buffered saline (pH7.2) containing antibody stabilizer. |
| Storage: | The antibodies are stable for 12 months from date of receipt when stored at -20°C. The antibodies can be stored at 2°C-8°C for one month without detectable loss of activity. Avoid repeated freezing-thawing cycles. |
| Gene ID: | 16001 |

Applications/Suggested Working Dilutions

| | |
|-----------------------------|---------------|
| Western Blot: | 0.1-1ug/ml |
| ELISA: | 0.01-0.1ug/ml |
| Immunoprecipitation: | 2-5ug/ml |
| IHC: | 2-10ug/ml |
| Flow Cytometry: | 5-10ug/ml |

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World
184 Park Drive, Milton Park
Abingdon, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482



North America
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218



Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880



Switzerland
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85