

IgM

Immunoglobulin Type M

Catalog #: 500-1764
Size: 100 µg
Host: Rabbit
Isotyping: Rabbit IgG
Applications: IHC(P)
Reactivity: Hu
Product Preparation: N/A
Background: N/A

PURIFICATION

The Rabbit IgG is purified by Epitope Affinity Purification

SPECIFICITY

N/A

FORMULATION

This affinity purified antibody is supplied in sterile Phosphate buffered saline (pH7.2) containing antibody stabilizer.

STORAGE

The antibodies are stable for 12 months from date of receipt when stored at -20oC. The antibodies can be stored at 2oC-8oC for one month without detectable loss of activity. Avoid repeated freezing-thawing cycles.

Gene ID

N/A

Applications/Suggested Working Dilutions

Western Blot: 0.1-1 µg/ml

ELISA: 0.01-0.1 µg/ml

Immunoprecipitation: 2-5 µg/ml

IHC: 2-10 µg/ml

Flow cytometry: 5-10 µg/ml

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85