



Oct4 mRNA

Code: 7001-101

Description: Human Oct4 mRNA – length 1154 nt + 120 A (polyA tail)

Concentration: ~ 1 µg/µl (final concentration will be provided in product insert)

RNA concentration will be determined by spectrophotometer (1 A₂₆₀ = 40 µg/ml RNA)

Acceptable ratio of A₂₆₀/A₂₈₀ : 1.8 – 2.3 (in TE buffer, pH 8)

The RNA is transcribed with T7 RNA polymerase and ARCA-NTP-Mix (containing 100% m5CTP and 100% ΨTP) from a linear template generated by PCR. Purified mRNA is treated with Antarctic phosphatase and re-purified.

RNA will be analysed on Agilent Bioanalyzer 2100.

ISO 13485 certified product

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85