

Mouse Anti-Ovalbumin IgG1 Monoclonal Antibody, Clone L71

Catalog # 7093

For Research Use Only - Not Human or Therapeutic Use

DESCRIPTION:	Mouse anti-ovalbumin (OVA) IgG1 monoclonal antibody, clone L71, unconjugated
APPLICATION:	Use for ELISA
QUANTITY:	0.1 ml, 1 mg/ml
FORM:	Dissolved in PBS, pH 7.2
SOURCE:	Mouse
REACTIVITY:	Reacts to native OVA, denatured OVA, and OVA peptides
IMMUNOGEN:	Native OVA
EPITOPE:	OVA 363-374: HPFLFCIKHIAT
PURITY:	Protein G purified
SUBTYPE:	Mouse IgG1
USAGE:	Dilute the antibody in buffer containing 2% heterologous serum, such as normal goat serum, to prevent non-specific reactions in ELISA.
STORAGE TEMPERATURE:	-20°C
STABILITY:	1 year

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85