

Human alpha Thrombin

Catalog #: 80111

Lot#: 150731

Description: The Human Thrombin was activated from homogeneous human prothrombin by activation with Factor Xa, Factor Va, and phospholipid. These activating enzymes were removed after activation. Complete activation was observed on SDS-PAGE. The above protein was purified from human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease). MW = 37 kDa.

Concentration: 3.41 mg/ml

Aliquots: 1,000 units

Volume/aliquot: 0.097 ml

Specific Activity: 3027 NIH units/mg

Total Protein: 0.33 mg

Buffer: 50 mM sodium citrate, 0.2 M NaCl, 0.1% PEG-8000, pH 6.5

Storage: < -70°C

Expiration Date: 2/2017

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

North America 
1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

Europe 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0)72 8080244
F: +31 (0)72 8080142

Switzerland 
Via Lisano 3,
(CP.683)
CH-6900
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

Germany 
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880