

## Data Sheet

### BINDING AGENT

**CATALOG #:** AMS.60390

**DESCRIPTION:** Extra binding agent for use with the PDE assay kit (Cat. # AMS.60300) or related PDE assay kits for the study of phosphodiesterase activity.

**APPLICATIONS:** Ideal for detection of phosphodiesterase activity, and screening of small molecule inhibitors for drug discovery and HTS applications.

**SUPPLIED AS:** Aqueous solution

**QUANTITY:** 100 µL

**STORAGE:** Store in aliquots at +4°C.

**STABILITY:** >6 months at +4°C

#### REFERENCES:

1. Jeon, Y., et al. (2005). *Cell Mol. Life Sci.* **62** (11): 1198–220.
2. Weiss, B. (1975). *Adv. Cycl. Nucl. Res.* **5**:195-211.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**UK & Rest of the World**  
AMS Biotechnology (Europe) Ltd  
184 Park Drive, Milton Park  
Abingdon, OX14 4SE, UK  
T: +44 (0)1235 828 200  
F: +44 (0)1235 820 482



**North America**  
amsbio LLC  
1035 Cambridge Street,  
Cambridge, MA 02141  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218



**Germany**  
AMS Biotechnology (Europe) Ltd  
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880



**Switzerland**  
AMS Biotechnology (Europe) Ltd  
Centro Nord-Sud 2E  
CH-6934 Bioggio-Lugano  
T: +41(0) 91 604 55 22  
F: +41(0) 91 605 17 85