

Anti-Chondroitin Sulfate D Monoclonal Antibody

Product Information

Product No.: AMS.A2872

Volume/Quantity: 0.1 mg

Product Form: Purified mouse monoclonal antibody/IgM, κ -chain

Buffer: Phosphate buffered-saline(-) (PBS(-))

Preparation: The antibody is purified from mouse ascites fluid.

Preservative: 0.05% sodium azide

Stabilizer: None

Immunogen: Chick embryo limb bud proteoglycan (PG-M/versican)

Specificity: The monoclonal antibodies preferentially recognize chondroitin sulfate containing the GlcUA(2-O-sulfate) β 1-3GalNAc(6-O-sulfate) disaccharide unit (D-unit). No cross-react with hyaluronic acid, dermatan sulfate, keratan sulfate and heparan sulfate.

Clone Name: MO-225

Applications

Immunohistochemistry, ELISA, Western blotting

Storage

Store in freezer (-20°C to -80°C). It is not recommended storing at 4°C. Avoid repeated freezing and thawing

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85

Anti-Chondroitin Sulfate D Monoclonal Antibody (Product No: AMS.A2872)

Note:

For research use only. Not for Diagnostic or Therapeutic Use.

References

A monoclonal antibody that specifically recognizes a glucuronic acid 2-sulfate-containing determinant in intact chondroitin sulfate chain

M. Yamagata, K. Kimata, Y. Oike, K. Tani, N. Maeda, K. Yoshida, Y. Shimomura, M. Yoneda, S. Suzuki, J. Biol. Chem. 1987, 262, 4146-4152.

Selective distributions of proteoglycans and their ligands in pericellular matrix of cultured fibroblasts

M. Yamagata, S. Saga, M. Kato, M. Bernfield, K. Kimata, J. Cell. Sci. 1993, 106, 55-65.

FAB CID-MS/MS characterization of tetrasaccharide tri- and tetrasulfate derived from the antigenic determinant recognized by the anti-chondroitin sulfate monoclonal antibody MO-225

T. Ii, M. Kubota, T. Hirano, M. Ohashi, K. Yoshida, S. Suzuki, Glycoconj. J. 1995, 12, 282-289.

Distribution of chondroitin sulphate proteoglycans and peanut agglutinin-binding molecules during bovine fetal palatine ridge formation

M. Takanosu, H. Amasaki, S. Matsumoto, K. Kimata, J. Anat. 1996, 189, 109-115.

Structural characterization of the epitopes of the monoclonal antibodies 473HD, CS-56, and MO-225 specific for chondroitin sulfate D-type using the oligosaccharide library

Y. Ito, M. Hikino, Y. Yajima, T. Mikami, S. Sirko, A. V. Holst, A. Faissner, S. Fukui, K. Sugahara, Glycobiology 2005, 15, 593-603.

Influences of retinal axons on the cultural substrate containing biotin-conjugated chondroitin sulfate in vitro

S. Ando, N. Sugiura, K. Kimata, H. Ichijo, Anat. Sci. Int. 2010, 85, 189-193.

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85