

Anti-Chondroitin Sulfate D Monoclonal Antibody

Product Information

Product No.: AMS.A2872

Volume/Quantity: 0.1 mg

Product Form: Purified mouse monoclonal antibody/IgM, κ -chain

Buffer: Phosphate buffered-saline(-) (PBS(-))

Preparation: The antibody is purified from mouse ascites fluid.

Preservative: 0.05% sodium azide

Stabilizer: None

Immunogen: Chick embryo limb bud proteoglycan (PG-M/versican)

Specificity: The monoclonal antibodies preferentially recognize chondroitin sulfate containing the GlcUA(2-O-sulfate) β 1-3GalNAc(6-O-sulfate) disaccharide unit (D-unit). No cross-react with hyaluronic acid, dermatan sulfate, keratan sulfate and heparan sulfate.

Clone Name: MO-225

Applications

Immunohistochemistry, ELISA, Western blotting

Storage

Store in freezer (-20°C to -80°C). It is not recommended storing at 4°C. Avoid repeated freezing and thawing

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85

Anti-Chondroitin Sulfate D Monoclonal Antibody (Product No: AMS.A2872)

Note:

For research use only. Not for Diagnostic or Therapeutic Use.

References

A monoclonal antibody that specifically recognizes a glucuronic acid 2-sulfate-containing determinant in intact chondroitin sulfate chain

M. Yamagata, K. Kimata, Y. Oike, K. Tani, N. Maeda, K. Yoshida, Y. Shimomura, M. Yoneda, S. Suzuki, *J. Biol. Chem.* 1987, 262, 4146-4152.

Selective distributions of proteoglycans and their ligands in pericellular matrix of cultured fibroblasts

M. Yamagata, S. Saga, M. Kato, M. Bernfield, K. Kimata, *J. Cell. Sci.* 1993, 106, 55-65.

FAB CID-MS/MS characterization of tetrasaccharide tri- and tetrasulfate derived from the antigenic determinant recognized by the anti-chondroitin sulfate monoclonal antibody MO-225

T. Ii, M. Kubota, T. Hirano, M. Ohashi, K. Yoshida, S. Suzuki, *Glycoconj. J.* 1995, 12, 282-289.

Distribution of chondroitin sulphate proteoglycans and peanut agglutinin-binding molecules during bovine fetal palatine ridge formation

M. Takanosu, H. Amasaki, S. Matsumoto, K. Kimata, *J. Anat.* 1996, 189, 109-115.

Structural characterization of the epitopes of the monoclonal antibodies 473HD, CS-56, and MO-225 specific for chondroitin sulfate D-type using the oligosaccharide library

Y. Ito, M. Hikino, Y. Yajima, T. Mikami, S. Sirko, A. V. Holst, A. Faissner, S. Fukui, K. Sugahara, *Glycobiology* 2005, 15, 593-603.

Influences of retinal axons on the cultural substrate containing biotin-conjugated chondroitin sulfate in vitro

S. Ando, N. Sugiura, K. Kimata, H. Ichijo, *Anat. Sci. Int.* 2010, 85, 189-193.