

HUMAN ADIPOCYTE DIFFERENTIATION MEDIUM

Human Adipocyte Differentiation Medium was designed to induce the differentiation of preadipocytes in culture. This medium is fully supplemented and contains 10% Fetal Bovine Serum. Human Adipocyte Differentiation Medium is sterile filtered and bicarbonate buffered, and ready to be used for differentiation of preadipocytes to adipocytes.

Catalog No.: AMS.SMH811D250

Quantity: 250ml

Quality Control:

Bioassay: Each lot was tested for the ability to induce differentiation of preadipocytes to adipocytes that accumulate lipid drops in cytoplasm within 2 weeks in this differentiation medium.

Sterility Test: Negative for bacteria, yeast and fungi.

Storage & Stability: Avoid exposure to light. Stable for 6 months at 4°C.

AMSBIO | www.amsbio.com | info@amsbio.com

AMSBIO LLC
USA & Canada 
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

AMSBIO Europe BV
EU 
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 8080244
F: +31 (0) 72 8080142

AMS Biotechnology (Europe) Ltd
UK & Rest of the World 
184 Park Drive, Milton Park
Abingdon OX14 4SE
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

AMS Biotechnology (Europe) Ltd
Switzerland 
Via Lisano 3,
(CP.683)
CH-6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

AMSBIO Europe BV
Deutschland 
T: +49 (0) 69 779099
F: +49 (0) 69 13376880