



Product Information: Recombinant Human SCF, 50µg

1. Introduction

SCF (stem cell factor) is a growth factor secreted from fibroblasts and endothelial cells. SCF binds to a receptor known as c-Kit and is thought to play an important role in the survival, proliferation and differentiation of hematopoietic stem cells. Recombinant Human SCF is an 18.5 kDa polypeptide containing 164 amino acid residues, which corresponds to the sequence of the secreted soluble form of SCF.

Source: *Corynebacterium glutamicum*

Amino Acids Sequence:

EGICRNRVTN NVKDVTKLVA NLPKDYMITL KYVPGMDVLP SHCWISEMNV QLSDSLTDLL

DKFSNISEGL SNYSIIDKLV NIVDDLVECV KENSSKDLKK SFKSPEPRLF TPEEFFRIFN

RSIDAFKDFV VASETSDCVV SSTLSPEKDS RSVTKPFML PPVA



2. Materials Provided

Component	Volume	Storage
Recombinant Human SCF, 50µg (0.1 mg/ml)	500 µl	Store at -80°C

3. Preparation for Use

Recombinant Human SCF is provided frozen and can be stored at -80°C until use. Before use, this frozen product should be thawed on ice. Avoid repeated freeze-thaw cycles.

4. Precaution and Disclaimer

Recombinant Human SCF is for research use only and is not for diagnostic use, therapeutic use (including clinical research), or as a food. SCF is dissolved in 20 mM HEPES (0.2 µm filtered).

For Research Use Only. Not Intended for Human or Animal Diagnostic or Therapeutic Uses

Version 2 (Jan 2022)

AMSBIO | www.amsbio.com | info@amsbio.com

UK & Rest of the World
184 Park Drive, Milton Park
Abingdon OX14 4SE.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

North America
1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

Europe
Berenkoog 41,
1822 BH Alkmaar,
Netherlands
T: +31 (0) 72 808 0244
F: +31 (0) 72 808 0142

Switzerland
Via Lisano 3,
(CP.683)
CH-6900
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85

Germany
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880