

PRODUCT SPECIFICATION

Product Code: AU.1000.SLC_10MM

Product Name: Silicon wafer chips (10mm x 10mm) coated with 100nm gold over a Titanium adhesion layer

Appearance (Color) : Gold

Appearance (Form) : Silicon Wafer Chips

Vendor Information:

Silicon Wafer, 10mm x 10mm × 500 μm

Nominal Film Thickness Settings: 50 Angstrom Ti, 1000 Angstrom Au

Packaging : 20 chips in gel box

Assay: 99.99% Ti, 99.999% Au

AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, UK
T: +44 (0)1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Centro Nord-Sud 2E
CH-6934 Bioggio-Lugano
T: +41(0) 91 604 55 22
F: +41(0) 91 605 17 85