

Animal adult normal multiple frozen tissue section panels

Quantity: 5 slides.

Storage Conditions: Keep at -80 °C. Half year from the date of receipt under proper storage condition.

Applications: Animal adult normal frozen tissue sections are ideal for:

- Rapid screening of genes or proteins against an extensive panel of tissues.
- Gene or protein expression pattern analysis in animal adult normal tissues.
- Comparison of locations and expression levels of novel genes or proteins in different specific animals tissues.

Quality Control: At least one of the tissue slides from each lot was stained with H & E to ensure the quality.

Description: Animal adult normal tissues were snap in liquid nitrogen immediately after excision and embedded in OCT. The tissue sections were 5-10 µm in thickness, mounted on positively charged glass slides and fixed by cold acetone. Frozen tissue section panels contain 8 different tissues on one slide.

Shipping: Slides are delivered on dry ice.

Intended use of the products: All products sold by AMSBIO are intended for research use only and must not be used for medical diagnostic or drug purposes. AMSBIO products may not resell, modified for resale or used to manufacture commercial products without written approval from AMSBIO.

Limited product warranty: AMSBIO has a total quality assurance policy. If any product does not meet the specification stated in the product data sheet will be replaced at no charge. This warranty limits AMSBIO liability to the replacement of the product only. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by AMSBIO. AMSBIO shall not be liable for any incidental, direct, special or consequential damages, even if it is aware of the possibility of such damages.