

Product datasheet for **GE100046**

pAAVS1-EF1a-puro-DNR CRISPR Vector



Product data:

Product Type: CRISPR Vectors

Function: AAVS1 donor vector

- Features:
- Donor vector for targeted transgene insertion into human AAVS1 site
 - Contains left and right AAVS1 homologous sequences (550 bp each side)
 - Multiple cloning sites for transgene cloning, compatible for shuttling from the TrueORF clones
 - The expression of the target gene will be under EF1a promoter, a constitutive promoter as an alternative of CMV promoter. the latter could be silenced in certain cellular environment
 - Puromycin mammalian selection under PGK promoter
 - Targeted insertion when cotransfected with pCas-Guide-AAVS1 (SKU GE100023)

Product images:



AMSBIO | www.amsbio.com | info@amsbio.com

 **UK & Rest of the World**
184 Park Drive, Milton Park
Abingdon OX14 4SE, U.K.
T: +44 (0) 1235 828 200
F: +44 (0) 1235 820 482

 **North America**
1035 Cambridge Street,
Cambridge, MA 02141.
T: +1 (617) 945-5033 or
T: +1 (800) 987-0985
F: +1 (617) 945-8218

 **Germany**
Bockenheimer Landstr. 17/19
60325 Frankfurt/Main
T: +49 (0) 69 779099
F: +49 (0) 69 13376880

 **Switzerland**
Via Lisano 3
(CP.683)
6900 Massagno
T: +41 (0) 91 604 55 22
F: +41 (0) 91 605 17 85