

## Product Sheet: Anti-Hepatitis B virus core antigen IgG fraction (monoclonal)

<b>Catalog Number</b>	HBM-021-5
<b>Product Id</b>	HBM-021-5
<b>Product Name</b>	Anti-Hepatitis B virus core antigen IgG fraction (monoclonal)
<b>Temperature</b>	-20°C
<b>Unit</b>	100µg
<b>Price</b>	295
<b>Availability</b>	Available
<b>Category</b>	Infectious Agents
<b>Subcategory</b>	Hepatitis B Virus
<b>Description</b>	<p>Anti-hepatitis B core antigen IgG fraction is a mouse monoclonal antibody (Clone 2C11/2) to native HBcAg particles.</p> <p>The antigen used for immunization is highly purified HBcAg produced by recombinant E. coli cells.</p> <p>It has been purified using Protein G.</p>
<b>Isotype</b>	IgG-2a
<b>Mol Weight</b>	N/A
<b>Purity</b>	N/A
<b>Storage</b>	Store at -20°C.
<b>Stability</b>	At least 1 year at -20°C.
<b>Biological Activity</b>	The antibody is reactive in ELISA and also immunoprecipitation of native hepatitis B core particles.
<b>Formulation</b>	Solution at 0.50 mg/ml in 0.1M Tris-HCl (pH 7.5)
<b>References</b>	N/A

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

 **UK & Rest of the World**  
 184 Park Drive, Milton Park  
 Abingdon, UK  
 T: +44 (0)1235 828 200  
 F: +44 (0) 1235 820 482

 **North America**  
 1035 Cambridge Street,  
 Cambridge, MA 02141  
 T: +1 (617) 945-5033 or  
 T: +1 (800) 987-0985  
 F: +1 (617) 945-8218

 **Germany**  
 Bockenheimer Landstr. 17/19  
 60325 Frankfurt/Main  
 T: +49 (0) 69 779099  
 F: +49 (0) 69 13376880

 **Switzerland**  
 Centro Nord-Sud 2E  
 CH-6934 Bioggio-Lugano  
 T: +41(0) 91 604 55 22  
 F: +41(0) 91 605 17 85